DeguDent GmbH
Rodenbacher Chaussee 4
63457 Hanau-Wolfgang
GERMANY
www.degudent.de
www.degudent.com

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Instructions for use (last revision: 07/07)

Intended use/indications:

Compartis Ti are individually manufactured crowns and bridge frames for porcelain veneering by a dental technician. Compartis Ti are custom-made. CTE (25–500 °C): 9,5 µm/m·K / (25–600 °C): 9,6 µm/m·K.

Contraindications:

Do not use in cases of hypersensitivity to one or more of the metals contained in the alloy.

Warnings:

- Do not inhale dust and vapours when machining.
- Ensure suitable air extraction/ventilation at the workplace and corresponding machinery.

Side effects:

Potential side effects include allergies to the metals contained in the alloy and discomforts caused by electro-chemical effects. Systemic side effects from the alloy metals are claimed in isolated cases.

Interactions:

Avoid occlusal and approximal contact between different alloy types.

Technical data:

Vicker's hardness	250 HV10
0.2% expansion limit	520 MPa
Tensile strength	610 MPa
Rupture limit	25 %
Density	4,5 g/cm ³

Chemical composition (mass-%):

99,0% Ti; <1% N₂; <1% C; <1% H₂; <1% Fe; <1% O₂

Processing:

1. A few early contacts may have to be removed before the coping fits in its final position. Any imperfections should always be removed selectively. If necessary, remove the early contacts locally with titanium carbide cutters.

2. Then please scan the frames initially in circular and incisal respectively occlusal mode. For a long-life frame structure, the values measured should be no less than the following wall/connector thicknesses:

Crown frame:	front and side teeth area: 0,4 mm
Abutment:	front and side teeth area: 0,4 mm
Bridge connector:	front and side teeth area: 9 mm ²

It is essential to ensure that the framework thickness is not reduced below the above minimum thicknesses!

3. Fit the marginal edge exactly to the preparation margin and grind to adjust the fit if necessary.

Compartis titanium copings should always only be machined with cross-toothed cutters suitable for titanium. The coping should be ground with slight contact pressure and only in one direction.

4. As preparation for veneering, now prepare the coping by shot peening carefully with aluminium oxide (100–150 µm and 2 bar pressure), followed by cleaning with the steam cleaner.

5. For passivation, any machined titanium surfaces should be allowed to rest for approx. 10–15 min (after finishing/sandblasting) before cleaning. The framework is ready for immediate veneering. When veneering, please comply with the instructions issued by Duceratin® Kiss.

Caution: The frame has been produced as modelled in the Cercon® art software. Reworking may possibly be necessary in the area of the connectors and the basal surfaces of the pontics, if Cercon® art was used without smoothing function.

The frames can now be veneered straight away:

Integration / fastening of Compartis restorations:

Crowns:	Conventional cementing
Multiple-unit bridges:	Conventional cementing

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