www.degudent.com

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Solaris®

Instructions for use
Istruzioni per l'uso
Instrucciones para el uso
Mode d'emploi

Hersteller:

DeguDent GmbH
Postfach 13 64
D-63403 Hanau
Telefon +49 (0)61815950
Telefax +49 (0)6181595858
www.degudent.de

© Dezember 2004 by DeguDent GmbH

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© December 2004 by DeguDent GmbH

We reserve the right to make changes and correct errors. This documentation may not be duplicated without our permission.

Contents:

Working with the DeguDent Solaris Device:	Page
1 Scope of supply	20
2 Technical specifications	20
3 For your personal safety	21
4 General user information	22
a. Intended use	
b. Description of components	
c. Menu navigation	
5 Getting started	23
a. Installation the DeguDent Solaris	
b. Connecting to the power supply	
c. Inserting the reference cathode	
6 System advantages	23
7 DeguDent Solaris indications	24
8 Preparation instructions	25
9 Dental procedure	26
10 Cleaning the anode	29
11 Electroforming with the DeguDent Solaris	30
12 Recycling of residual gold	31
13 Solaris: Examples of bar attachments	31
14 Solaris: Double loading	32
15 Gold weights	32
16 Disclaimer	33
17 Care and maintenance	33
18 Troubleshooting	33

1. Scope of supply:

1 DeguDent Solaris device
1 anode
6 cathode rods, long (color-coded blue)
8 cathode rods, short (color-coded black)
8 double cathode rods (color-coded black)
2 glass beakers, large
1 glass beaker, small
1 mains cable
1 instruction booklet
1 reference cathode (color-coded yellow)
1 precipitation cathode (color-coded red)
1 centering device
1 magnetic stir bar
1 cathode slot

2. Technical specifications:

Power Supply: 130/230 V~, 50/60 Hz
Maximum Wattage: 230 W
Size: 480 mm x 270 mm x 475 mm
Weight: 12.5 kg

For information on altitude, atmospheric pressure and relative humidity, please consult the Instructions for Use.

Only use the device indoors. The room temperature should be between 15 °C and 25 °C.
The mains supply's voltage fluctuations should not be greater than $\pm 10\%$ of the nominal value.

3. For your personal safety:

Please read these Instructions for Use carefully before installing and starting the DeguDent Solaris device. The information contained in these Instructions for Use must be followed. Inappropriate use of the device may be dangerous to users or third parties and could damage the device and other material assets.

- DeguDent Solaris may only be operated with original spare parts. Only in this way can the quoted performance data and operational reliability be guaranteed.
- Do not open the device. The unit's interior contains hazardous voltages.
- Only service technicians authorized by DeguDent may carry out maintenance and repair work.
- Protect the device from moisture in order to prevent short-circuiting and damage.
Do not expose the device to water vapor.
- Unplug the device from the mains if it is not being used for an extended period of time.
- DeguDent Solaris should only be used for dental purposes.

For your personal safety when handling Solaris fluids, please ensure you read the relevant safety datasheets for the

- gold bath
- activator
- conducting silver varnish
- anode cleaner and
- thinner.

4. General user information:

a. Intended use

DeguDent Solaris is a device for fabricating crowns, suprastructures and bridge frameworks made of 99.9% gold. The process involves duplicating a die and casting it in plaster or plastic. A thin (200–300 µm) gold layer is electroplated onto the duplicated die after conductive silver paint has been applied. The prosthetic appliances are characterized by very high fitting accuracy with a minimal coating thickness. To prevent plating defects, please ensure that work is done in an extremely clean environment and follow the instructions. In addition, you should regularly clean the anode, cathode slot, heating system, vessels etc. and immediately remove the residual electrolytes after each electroplating process. A special anode cleaner (REF 2545 0055) is available from DeguDent for cleaning the anode.

b. Description of components

- Anode: The anode is attached to the device with three connectors. It can be removed from the device, e.g. for cleaning, with slight pulling.
- Cathode Holder: It contains 14 slots for holding the steel cathodes. 6 in the inner ring and 8 in the outer ring. Six or eight slots can be used concurrently.
- Heating System: This is a two-phase heating system that offers optimal heating performance according to the liquid filling level. It is firmly attached to the device.
- Small and Large Beakers: They are required to hold the electroplating liquid. The small glass is for small amounts of liquid (100–350 ml), while the large glass is for large liquid amounts (350–1000 ml), depending on the number of electrodes.

c. Menu navigation

See point 11: Electroforming with the DeguDent Solaris

5. Getting started:

a. Installing the DeguDent Solaris

The device must be placed on a flat solid surface (standing base) at least 10 cm away from the wall or other objects. It should be positioned horizontally. The electroplating liquid must be protected from direct sunlight or other light sources.

b. Connecting to the power supply

Connect the mains cable to the fuse-protected circuit via a wall outlet. After connecting the mains cable, switch on the device by pressing the black ON button at the rear.

c. Inserting the reference cathode (colour-coded yellow)

Before operating the device, the reference cathode must be connected to the device head. The connector for the reference cathode is located to the right of the temperature sensor.

6. System advantages:

- Very high hardness
 - Homogeneous appearance
 - Very fine-grained structure
- } of the electrodeposited objects
- 1 magnetic stir bar for all beakers and liquid amounts
 - Just two beakers
 - Each cathode can be selected with an individual coating thickness
 - Overnight mode / overnight plating
 - Extra beaker for possible later plaster release/rinsing
 - Finer gradations of coating thicknesses
 - Up to 16 objects can be plated
 - The electroplating process lasts about 5 hours (200 µm)
 - RS 232 interface for updates
 - Clear display
 - Self-explanatory user menu
 - Anode heating (no risk of the housing being burned)
 - Simple plug system for the cathodes
 - Ready-to-use and recyclable steel cathodes

Please note: When starting the device for the first time, please clean its components (anode, beaker, magnetic stir bar, glass etc.) with distilled water before you begin using it. The anode should be cleaned with an anode cleaner.


7. DeguDent Solaris indications:

Inlays (veneered)	☺
Onlays (veneered)	☺
Single crowns (veneered)	☺
Interlocks	☺
Electrolytically bonded 3-unit bridges	☺
Telescopic copings	☺
Implant crowns	☺
Bar connectors	☺
The goldplating bath to be used	Solaris Goldplating bath
The brightener to be used	Solaris Activator

8. Preparation instructions:

General preparation information

The general metal ceramic preparation regulations apply to the fabrication of electroformed crowns and bridges. In order for the fabrication of esthetic dentures to be possible, the ceramic space required must be at least 1 mm. Practical experience has shown that a pronounced chamfer with rounded edges works well in the electroplating procedure. This gives the ceramic a supporting function. The preparation should follow the anatomy of the natural tooth. Please avoid sharp edges or a preparation that is too conical. It is essential to correct areas of unevenness and edges while preparing the die. The surface must be smooth, otherwise it will not be possible to electroplate it cleanly. The cavity angle should be about 6° for inlays. The preparation depth is about 2 mm.

Temporary seating

As a general rule, temporary seating of electroformed crowns and bridges is not recommended. The excellent fit of electroformed crowns would force the dentist to remove the provisionally cemented crowns with the aid of a claw. This generally results in the crown margin or ceramic being destroyed.

Cementation

All conventional cements (glass ionomer cement, zinc phosphate cement) can be used for luting. DeguDent recommends the cement Dentsply Dyract® Cemplus.

9. Dental procedure:

Impression

- The dental impression must be taken as precisely as possible. A clearly defined preparation margin lays the foundation for later success.

Model fabrication

- The master model can be fabricated in the usual way with any conventional model system after the dental impression has been taken.

Die preparation

- Uneven parts of the die must be blocked out with plastic or blockout wax.
- Because of problems associated with sharp edges in electroforming, all right angles and acute angles must be rounded off or blocked out.
- The preparation limit of the single die is marked and later sealed with a die hardener.
- An interspace varnish is then applied, which serves as a spacer for the cement.

Duplication

- The prepared single die is fixed in a duplicating flask. The die is stirred under vacuum at a 1:1 ratio and duplicated with Deguform® (REF 2511 0002) or Deguform® Plus (REF 2511 0502).
- Let the silicone cure well (at least 45 minutes).

Fabricating the duplicate die

- The duplicating mold is injected with Waxit® surfactant (REF 2518 0017).
- After a brief reaction time, blow the silicone mold thoroughly dry.
- Type IV gypsum (e.g. Picodent 280) is mixed according to the manufacturer's instructions and a vibrator is used to cast it into the duplicating mold. Please do not use any plastic-stabilized plaster.
- After the plaster has set, remove the duplicate die from the duplicating mold.
- Grind and round off the plaster die until it is small enough for there to be room under the preparation limit to drill a borehole 2 mm in diameter.
- A 3 mm deep hole with 2.0 mm diameter is drilled under the preparation limit with a cutter (REF 2351 0367).

Attaching the cathode

- Apply some Degufix cyanoacrylate adhesive (REF 2351 1365) to the clean metal rod, insert it into the drilled hole (about 3.0 mm deep) and bond it in.
- Immediately remove excess adhesive and let it dry completely.

Applying the Solaris conductive silver paint

- Shake the bottle of conductive silver paint well before applying it.
- The conductive silver paint may only be applied to a dry and dust-free model die.
- Apply the conductive paint evenly with the brush up to the preparation limit and make a connection to the exposed metal rod.
- Only apply a single coat. Touch up any defective spots.
- Let the conductive silver paint dry for about 30 minutes.
- Important: Please stick to this time in order to prevent contamination of the goldplating bath.

Direct application of conductive silver paint on telescopic crowns/bars etc.

- Fill the crown with plastic (pattern resin).
- After curing, use the bur (REF 2351 0367) to drill a hole with 2.0 mm diameter in the plastic.
- Put some Degufix (REF 2351 1365) cyanoacrylate adhesive on the clean metal rod, insert it into the drilled hole (about 3.0 mm) and bond it in.
- Immediately remove excess adhesive and let it dry completely.
- Important: Degrease the crown to be electroformed before applying the conductive silver paint.
- Shake the conductive silver paint well before applying it.
- Apply the conductive paint evenly with a brush up to the preparation limit and make a connection to the exposed metal rod.
- Only apply a single coat. Touch up any defective spots.
- Let the conductive silver paint dry for about 30 minutes.
- Please stick to this time in order to prevent contamination of the goldplating bath.
- With metal telescopes, cover the uneven shoulder with plastic.

Recommended alloys for male parts

- Degunorm®
- Degunorm® supra
- Degunorm® pur
- BiOcclus® HT
- In alloys that contain copper, slight discolorations may appear, which pose no health risks.

Inserting the electrode disk

- The prepared electrodes are attached to the glass slide in the direction of the arrow.
- The glass slide is placed on the holder and locked from the outside with the lever.
- The parameters Coating Thickness, Object Type and Automatic Plaster Release On/Off are entered into the device.
- The display shows the required goldplating bath and activator amounts as well as the beaker glass to be used.
- The user must fill in the required amount accordingly.
- Press the START key to begin the process.

After electroforming:

Remove the contact rods

- Unlock the device and remove the glass slide.
- Gently twist the metal rods to pull them out.
- If the electroformed copings have already been removed from the plaster, the electroformed crowns are separated with a cutting wheel at the junction to the metal rod.
- Remove any leftover plaster residue in the ultrasound bath.

Removing primary crowns and telescopic copings

- Reduce the overlapping margin under a microscope, using a rubber polisher.
- Drill a screw into the plastic die.
- Hold the die on the screw with pliers and carefully tap the pliers to detach male and female parts from each other.

Removing the conductive silver paint

- Use nitric acid to remove the conductive silver paint in a fume hood inside a plastic container. Please carry personal protective equipment.
- If no nitric acid is available, the conductive silver paint coat must be removed by blasting with alumina (50 µm and 1–2 bar).
- Please be aware that no blasting should be done with telescopic copings, because otherwise the inner surface is roughened and there may no longer be friction.
- Nitric acid concentration: 20% to approximately 65%. The higher the concentration, the faster the conductive silver varnish will be dissolved.

Please note:

Conductive silver paint residue may cause discoloration of the ceramic, telescopic copings and gingiva and may lead to the crown margin being deformed.

Finish on the master model single die

- A rubber polisher is recommended for correcting margins.
- Do not use burs.

Surface conditioning

- Prepare electroformed crowns for veneering by blasting them with 110 µm of alumina and 1 bar of pressure.
- The crown is then cleaned with the aid of a steam radiator.

Bonder firing

- A gold bonder is not required for the DeguDent system.
- DeguDent Deckgold NF is very well suited if gold bonder is to be used.

The ceramic veneer

- Duceragold® Kiss, Duceram® Kiss and SYMBIOceram® are all suitable for fabricating ceramic veneers.
- Please follow the ceramic manufacturer's data and recommendations.
- In order to prevent deformations on an inlay or onlay to be veneered, a relief cut should be made up to the opaquer during the first dentin firing.
- The ceramic's firing temperature should not exceed 950 °C.
- Veneered electroformed crowns should no longer be sandblasted.

10. Cleaning the anode

Please ensure that the anode basket, temperature sensor, reference electrode and heating system are cleaned after each electroforming process. If this is not done, there is a risk of contamination in the subsequent electroforming process.

Procedure:

- Carefully remove the anode after each electroforming process and clean it in a vessel for about 10 minutes with the anode cleaner in ultrasound bath provided. The anode is connected to the device head using three small connectors. Slightly pulling downward will loosen the anode from its socket. Immerse the anode in the anode cleaner, such that the mesh is completely covered.
- Thoroughly rinse off any residue of the anode cleaner from the anode under running water.
- Please note: The anode cleaner can soil clothing etc. Please wear rubber gloves and safety goggles.
- Wipe the glass heater with a damp cloth after each electroforming process.
- Also wipe the electrodes, the magnetic stir bar with a damp cloth after each electroforming process.
- Rinse out the glasses and measuring cylinders after each process.
- Regularly evaporate the cathode slots.
- Do not allow the anode contacts to touch the anode cleaner.
- When using the "Automatic Plaster Release On/Off" function, the cathode rods and the reference sensor must be removed from the plaster release.

Please ensure that your work is very neat. This has a positive effect on plating.


11. Electroforming with the DeguDent Solaris:

Step	Action	Menu Display	User Actions
1	Selection of "Electroforming" mode	<ul style="list-style-type: none"> Main Menu 	<ul style="list-style-type: none"> Select "Electroforming" from the menu Confirm your selection by pressing ENTER
2	Selection of Parameter 1	<ul style="list-style-type: none"> Number of electrodes, suprastructures and Plaster Release ON/OFF are displayed 	<ul style="list-style-type: none"> Adjust the number of electrodes or suprastructures Activate/deactivate plaster release Press START to change the menu
3	Selection of Parameter 2	<ul style="list-style-type: none"> Coating thicknesses and sample types are displayed Automatic determination and display of calculated liquid amount, activator and glass size 	<ul style="list-style-type: none"> Separately adjust the coating thicknesses of all samples Select all sample types Press START to change the menu
4	Preparation	<ul style="list-style-type: none"> Coating thicknesses, sample types and resulting slots are displayed Liquid amount, activator and glass size are displayed Visualization of sample placement 	<ul style="list-style-type: none"> Fill in the calculated goldplating bath amount and the activator in the specified glass Place the glass slide according to what is shown in the display If necessary, fill with 500 ml of plaster remover Press START to change the menu
5	Starting Electroforming	<ul style="list-style-type: none"> Slots and contact detection are displayed Dynamic display of time remaining (heating the electroforming bath, electroforming, plaster removal) Current bath temperature is displayed Dynamic updating of the process status 	<ul style="list-style-type: none"> Press START to begin electroforming

12. Recycling of residual gold:

- The recycling process is effected directly through the unit (see the unit's Instructions for Use).
- The recycling cathode is colour-coded in red.
- The amount of liquid has to be entered.
- Alternatively, the recycling sponge can also be used.

13. Solaris: Examples of bar attachments:

The DeguDent Solaris galvano system can be used for electroplating bar sleeves with a surface of up to 2200 mm² for one bar and a surface of 1100 mm² each for two bars in a single process. With bars, you should aim at a thickness of the electroplating layer of 300 µm.

Example for bars:

Formula: Surface F = (width x length) + 2 x (height x length)

Approximations:

- Current level 1: 50 mm²
- Current level 2: 65 mm²
- Current level 3: 90 mm²
- Current level 4: 115 mm²
- Current level 5: 150 mm²
- Current level 6: 200 mm²
- Current level 7: 265 mm²
- Current level 8: 300 mm²

Example:

2 bar abutments (current levels 4 and 5) plus bar sleeve (length, 25 mm; height, 4 mm; width, 3 mm)

Calculation:

2 bar abutments = 115 mm² + 150 mm² = 265 mm²

Bar sleeve = (3 mm x 25 mm) + 2 x (4 mm x 25 mm) = 275 mm²

Total area: 265 mm² + 275 mm² = 540 mm²

Glass slides 3 and 6 are available for electroplating bars.

- Separate work instructions are available for the bar attachments process.

14. Solaris: Double loading

- The cathode (colour-coded in black) with the two identical ends is used for the desired double loading.
- The desired layer thickness is the same for both crowns.
- The object type is determined as per the table below.

	Type 1	Type 2	Type 3	Type 4	Type 5	Type 6	Type 7
Type 1	Type 3						
Type 2	Type 4	Type 4					
Type 3	Type 5	Type 5	Type 6				
Type 4	Type 5	Type 6	Type 6	Type 7			
Type 5	Type 6	Type 6	Type 7	Type 7	Type 8		
Type 6	Type 7	Type 7	Type 8	Type 8	X	X	
Type 7	Type 8	Type 8	X	X	X	X	X

15. Gold weights

- The following weights are electrodeposited for the object types:

Weights (g)

Object type	200 µm	250 µm	300 µm
1	0.23	0.29	0.35
2	0.29	0.36	0.44
3	0.38	0.47	0.57
4	0.48	0.60	0.73
5	0.62	0.77	0.93
6	0.81	1.01	1.22
7	1.06	1.33	1.59
8	1.20	1.50	1.80

16. Disclaimer:

- The manufacturer accepts no liability for damage caused by arbitrary changes, inappropriate use or any other use not in accordance with the Operating Instructions.

17. Care and maintenance:

- Check the device regularly for mechanical damage. In this way, you will prevent accidents and breakdowns.
- Completely remove contaminants from the device, the anode and the heating system. Only by doing this can the components be guaranteed to work properly.

18. Troubleshooting:

- See the Error Messages chapter in the Operating Instructions for the device.

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