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Cergo[®] press

Instructions for use
Mode d'emploi
Istruzioni per l'uso
Instrucciones para el uso

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Seite

D	Gebrauchsanweisung	3 – 10
GB	Instructions for use	11 – 18
F	Mode d'emploi	19 – 26
I	Istruzioni per l'uso	27 – 34
E	Instrucciones para el uso	35 – 42

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TABLE OF CONTENTS

1 Working with the Cergo® press	
1.1 Packing list or Contents	11
1.2 Accessories available	11
1.3 Technical data	11
1.4 For your personal safety	11
1.5 General information for the user	11
1.5.1 Proper use	11
1.5.2 Special features	11
1.5.3 Control elements	12
1.5.4 Display	12
2 Putting the unit into operation	
2.1 Installation of the unit	12
2.2 Electrical, compressed air and vacuum connections	12
2.3 Normal operating mode	13
2.4 Starting the pressing programs	13
3 Programming	
3.1 Programming mode	15
3.2 Copying a program	15
3.3 Program parameters	15
3.4 Service	15
4 Card operation mode	
4.1 Operating with the card	16
4.2 Key card (key)	16
4.3 From card to program	16
4.4 From program to card	16
4.5 Deleting card program	16
5 Cleaning and maintenance	
5.1 General	16
6 Elimination of errors	
6.1 Fuse breaks the circuit	17
6.2 Error messages	17
7 Statement of conformity	18

1 Working with the Cergo® press

GB

1.1 Packing List or Contents

- 1 unit, complete with pressure reducer
- 1 connecting hose for vacuum pump
- 1 connecting hose for compressed air supply
- 1 power cord
- 1 firing tray with supporting pins
- 2 ring plungers
- 2 positioning inserts
- 1 key card
- 1 memory card
- 1 tweezers
- 1 muffle tongs

1.2 Accessories available

Item	Product No.
Vacuum pump	5555 0100
Memory card	5555 0201
PC-connecting cable	5555 0101

1.3 Technical data

Rated voltage	200 V – 240 V
Max. power consumption	1300 W
Mains fuse required	10 A
Operating temperatures	4 °C to 40 °C (40 °F to 104 °F)
Dimensions of unit in mm:	
	475(h) x 345(w) x 520(d) (18.7" * 13.6" * 20.5")
Weight of unit	approx. 30 kg (66 lbs)

1.4 For your personal safety

-  Carefully read these instructions for use prior to installation and putting the Cergo® into operation. The information contained in these operating instructions must be observed.
-  Cergo® press must only be operated with original spare parts and accessories. Only in this way the performance data given can be achieved and the required operational safety is ensured.
-  Do not open the unit. Dangerous voltages are applied inside the unit. Therefore maintenance and repair work may only be carried out by service technicians authorized by DeguDent.
-  Protect the unit against moisture in order to avoid short circuits and damage to the unit.
-  Unplug the unit if it is not used over longer periods.

1.5 General information

1.5.1 Proper use


- Firing and pressing of porcelain or ceramic materials at air or under vacuum

1.5.2 Special features

- Reproduceable firing results
- Stationary working platform
- Free programming of increase, time, vacuum and temperature
- Memory card for program storage
- PC communication

1.5.3 Control elements

fig. 1

On/Off key

Switches the furnace on or off. By pressing the key, the display is switched on. The furnace is switched off if the key is pressed again.


Function keys

Function keys are defined by the respective word that is shown above them. Their function depends on the mode of the furnace. Possible functions: Select mode, abort firing programs, edit program parameters.

Aborts of the programming mode and returns to the previous screen.

Raises or lowers the muffle from the current position.

Adjustment wheel

Turning the knob: Browsing of programs, changing data or selecting options.

Pressing the knob: Selection of the parameter and confirming the entry.

1.5.4 Display

fig. 2

Press On/Off key. The display shows name, software version and the serial number.

fig. 3

Menu options

The display changes after 3 seconds and shows: current muffle temperature, current date and time, current pressure in the chamber and 4 menu options for the function keys.

The chamber pressure is a measurement of the current barometric pressure. Through the use of the vacuum during the operation the value is changed. A value of less than 50 hPa is expected in case of full vacuum. Repeated use may result in an increase of the vacuum level.

2 Putting the Cergo® press into operation

2.1 Installation of the unit

- Do not hold the furnace at the upper part of the housing.
- Store the carton and the packaging material in case of a future shipment.
- Shipping damage should be reported to the forwarding company as soon as possible.
- Remove the packaging material from the inside and the outside of the furnace.
- When installing the unit, the furnace should have a minimum distance of 15 cm from walls, shelves and heat-sensitive materials.
- Do not place the furnace below shelves or other objects that may block circulation of air.
- Prior to first firing, program No. 6 "Dentine firing 1" should be run.

2.2 Electrical, compressed air and vacuum connections

fig. 4

Compressed air controller Vacuum hose Fuse Vacuum pump Mains cable

Power cord

- Connect the power cord via a socket to a protected circuit for at least 10 A for the 230 V model. This circuit should only supply the furnace. Extension cables are not permissible.

Vacuum pump

- Connect the hose (do not bend) and electrical cables as shown.

Compressed air supply

- Connect the compressed air hose to the compressed air controller (fig. 4). The recommended pressure is 5 to 6 bar (4.5 bar is required for the process). The pressure can be changed via the pressure controller that is mounted on the rear panel of the unit. In order to check the pressure that has been set, change the display as shown in fig. 6a.


General

- After connecting the mains cable the muffle is raised.
Remove the foam packaging material from the inside of the muffle.
- Press the  key and wait until the screen shown in fig. 2 is displayed. Press the  key, to lower the muffle now.

2.3 Normal operation mode

Switch the unit on; after an internal test (5 – 10 seconds) the starting screen is displayed.

fig. 5

The lower section of the display shows 4 menu options above the function keys.

GGG – exclusive listing of GoldenGate programs -  symbol indicates GoldenGate.

Press – listing of ceramic pressing programs  symbol indicates ceramic pressing programs.

Prog – general listing of all programs.

Service – utilities and functions such as setting time/date, setting program type, PC communication, memory card mode.

Pressing of the function key below the respective menu name activates this function. The highlighted function can also be activated by pressing the “adjustment knob”. The **GGG** function displays the first 8 GGS programs. The first 24 programs are reserved for GGS.

Pressing of the **Press** or **Prog** keys results in program listing of these programs.

Turn the adjustment knob to scroll through these programs.

4 new menu options are displayed above the function keys.

ESC – returns to the display in fig. 5. Same function as ESC key.

Edit – listing of all program parameters. Fig. 6.


Turning the adjustment knob moves the cursor to a new parameter. Pressing the knob selects the parameter for editing.

Display – listing of all program parameters to check them prior to starting the selected firing program. (fig. 6)

Copy – copies the current program parameters to a new program number.


Program parameter information

fig. 6

List of programs and parameters

fig. 6a

Set air pressure

Parameter	Description	Program range
Dry Temp	Drying of the object	50 °C - 250 °C
Dry	Time for drying of the object	00:00 - 59:59 min
Close	Gradual closing time of muffle	00:00 - 59:59 min
Lo Temp	Starting temperature of the firing cycle	Dry - 800 °C
Lo Temp	Time to bring the object to the temperature	00:00 - 59:59 min
Rate	Temperature change per min.	1°C/min - 140 °C/min
Vacuum	Vacuum pump operation mode (The value displayed is the absolute pressure in the chamber)	On - pump runs until max. value is reached Cont - pump runs continuously Off - pump off
Vac On	Temp. at which vacuum is built up	50 - 1200 °C
Vac Off	Temp. at which vacuum is decreased	50 On - 1200 °C
Hi Temp	Final temperature of the firing cycle	50 - 1200 °C
Hold V	Holding time at final temperature with vacuum	00:00 - 59:59 min
Hold	Holding time at final temperature without vacuum	00:00 - 59:59 min
Press	Time for pressing	00:00 - 59:59 min
Temper	Holding time at tempering temperature	00:00 - 59:59 min
Temper	Tempering temperature	650 °C - 1200 °C
Cool	Gradual opening of the muffle	00:00 - 59:59 min

Press the **Display** function key.

Three new menu options are displayed above the function keys (fig. 6 and 6a).

Graphic – Shows a curve with relevant data for the selected program.

Edit – Activates the “Edit” mode. Turning of the adjustment knob moves the cursor to a new parameter. Pressing the knob selects the parameter so that it can be edited.

Start – Starts the firing program. The display changes and shows the graphics mode.

After the start of the program the **Start** function key takes over the **Stop** function.


2.4 Starting the pressing programs

Press the “Press” function key (fig. 5). The display lists 8 pressing programs (fig. 7). Program 1 includes the firing and pressing parameters for Cergogold. Select the “Display” key and the screen shows the program parameters that are preset (fig. 6a) The selected program is started by pressing the “Start” function key again or the adjustment knob. If program parameters are to be changed, proceed as described under item 2.3.

After the start of the program the **Stop** function key takes over the **Start** function. After reaching the start temperature, an acoustic signal can be heard and the **Stop** function changes to the **Start** function again.

The preheated pressing flask can be placed in the center of the firing tray. After inserting the ceramic pellets and the pressing plunger, the pressing cycle is started by pressing the **Start** function key again.

fig. 7

3 Programming

3.1 Programming modes

Various methods are available to enter programs into the furnace.


- Manual entry of new programs.
- Copying to new position and changing existing programs.
- Loading programs from the memory card into the furnace.
- Loading programs from a PC into the furnace via a communication interface.

Entering a new program

To enter a new program on an empty program position (number) of the furnace:

- Press **Prog** key
- Turn adjustment knob until an empty program number is displayed (no name)
- Press **Edit** key

fig. 8

- Turn adjustment knob counterclockwise.
- Cursor flashes in the Name field.
- Press >>Set>> function key to advance the cursor to the first letter.
- Turn adjustment knob (any direction) to display letters and symbols.
- Press >>Set>> to save the symbol that is displayed and to advance to the next position. Up to 14 characters can be saved.
- Press <<Set<< to move cursor backward (for correction).

The **ESC** key aborts editing of the program name and the cursor (▶ Cursor ↗) is displayed at the first parameter without saving the name. “Save” will save the name, then confirm with “Edit”.

- Turn adjustment knob clockwise to move the cursor to a new parameter.
- Press adjustment knob to select editing of a parameter.
- Turn the adjustment knob clockwise to increase a parameter value resp. counterclockwise to reduce it.
- Press the adjustment knob to save the parameter temporarily.

The **ESC** key aborts editing of the parameters and moves the cursor back to the current parameter.

NOTE: Press **Save** key to save the parameter permanently. If the adjustment knob is pressed, the parameters will only be saved temporarily.

3.2 Copying a program

Each program can be copied to a non-GGS position (programs 25 to 100). The existing program data at the new position will be overwritten.

From a program listing:

- Turn the adjustment knob to a program which is to be copied.
- Press **Copy** key. (This key is displayed in each program listing.)
- Turn the adjustment knob to a non-GGS program number. Copying will overwrite existing data in the new program.
- Press adjustment knob or **Save** key to store the program permanently.

Press **ESC** key to abort the copying process.


3.3 Program parameters

The possible program parameters are included in the paragraph 2.3. The following conditions only refer to some of the parameters:

- **GGS** programs can not be changed permanently. GGS parameter ranges are limited to the process specifications.

3 Programming

3.1 Programming mode

* is not available.

- If vacuum is programmed to “Off”, then “Vac On”, “Vac Off” and “HOLD V” time are not available.
- “Drying” Temp is the temperature at the working platform and is controlled via the muffle height and the “Predry” temperature.
- The “Temper” time is not available if it is set to 0.
- If “Temper” time is set to a time (> 0) , the muffle opens slightly and cools down to “Temper” temperature before it is closed again and keeps the “Temper” temperature.

3.4 Service

The **Service** screen provides access to a series of furnace controls and functions. Access to the **Service** function is possible from the stand-by mode, fig. 3.

- Press **Service** key.
- Turn adjustment knob clockwise to move the cursor to a function.

Press the adjustment knob to select the function.

StdbY Temp: Idle temperature of the muffle between firing programs. Range: 300 – 600 °C.

Calib Temp: In this mode the furnace can be calibrated up to +10 °C resp. -10 °C. The preset value is 0 °C.

Favorit Set: This function is not available in the Cergo® press.

Clear ✓. : Remove marking; the program is not listed as Favorit (favorite).

Set time: Turn adjustment knob to change the value; press to move the value in the time and date display to the next position.

PC Com(munication): Transfer programs from PC to Cergo® or copy programs from Cergo® to PC.

Prog Edit: Locks editing of program or cancels locking. The special key card (key) must be inserted to allow access to this function. A “Lock open” symbol () indicates that the “Prog Edit” function is not locked. After pressing the adjustment knob the “Lock closed” symbol () is displayed above the left function key (see fig. 1). Press this function key and the program parameters will be locked. To ensure additional safety, the serial number of the Cergo® furnace is also written on the card. Repeat the process with the key card (key) to cancel locking.

Select German, Select English, Select Français, Select Italiano, Select Español

Changes the menu items to the desired language.
The Cursor () on the right side points at the current selection.

4 Card operation mode

4.1 Operating the unit with the card (outside the Service menu):

Select display screen as shown in fig. 5 and insert the memory card (Back-up) to start a firing program via the card. The “Favorit” function key is replaced by the card symbol (). Pressing this function key shows a listing of the card programs available. When the furnace is delivered, the Back-up card is empty.

The service menu includes:

4.2 Key card (key)

Card to lock the Edit function or to unlock it.

Caution! Each Cergo® furnace has a different card. During locked operation mode only the use of this card allows access to the Service mode.

4.3 From program to card

Prog >>> : The user can transfer up to 40 programs from the Cergo® furnace to the inserted memory card (Back-up).

In order to transfer programs from the Cergo® furnace to the card, first set the adjustment knob to Prog >  and then proceed as described below (from card to program).

(Note: Programs can only be stored on memory cards (Back-up).

4.4 From card to program

 >>> Prog : Loads programs from the memory card (Back-up) or other preprogrammed cards into the Cergo® furnace. Select the service screen for this transfer type.

(Note: Transferring the  programs is not possible if the Edit function is locked).

Set adjustment knob to  > Prog and press the knob. Then move the cursor to the first program to be transferred and press the knob again. Turn the knob clockwise until all programs to be transferred are selected; then press the knob to start the transfer. Select the program position (number) for the first of the programmed programs and press the knob again. The programs will be stored in sequence.

4.5 Deleting card programs

Delete  : Deletes all programs on the memory card.

Turn adjustment knob to delete and press. “Delete” is shown on the left function key. All programs on the card will be deleted when this key is pressed. The card can now be used to store 40 new programs.

5 Cleaning and maintenance

5.1 General

Cleaning and maintenance ensure perfect operational reliability of the Cergo® press.

- Check the unit regularly to avoid accidents and unit failures.
- Immediately remove contaminations from all movable construction parts.

Lifting motor and the Cergo® vac vacuum pump are maintenance-free. No vacuum oil is required.

Prior to cleaning the unit must be switched off and unplugged.

During cleaning the unit make sure that no moisture penetrates into the unit. This will avoid short circuits which may result in a breakdown or a fire.

Inside the unit dangerous voltages are present. Maintenance work and repair must therefore only be performed by service technicians authorized by DeguDent.

Modifications to the unit by the user reduces the operational safety and increases the risk of damage to persons and materials.

6 Elimination of errors

6.1 Fuse breaks the circuit

- **The fuse breaks the circuit while the furnace is operated.**

Cergo® press should be isolated on a separate circuit without any other load.

- **During activation or start of the unit the internal fuse breaks the circuit.**

The fuse of the unit is located above the power cord socket (see fig. 4). Repeated breaking the circuit indicates a defect that must be eliminated by a service technician.

Single breaking the circuit may also be caused by a cold start of a vacuum pump.

The unit must be unplugged before the fuse is exchanged.

6.2 Error messages

If the microprocessor detects a malfunction, the display shows an error message with a short description and an error number. The Cergo® furnace interrupts the firing process unless there is a vacuum error or an error of the lifting motor; an acoustic signal can be heard and the unit must be unplugged.

The following error messages are possible:

- **No Temp signal**
- **Temperature too high**
- **No Temp increase**
- **Triac short-circuit**
- **Muffle open**

- **No vacuum**

The display of the vacuum sensor remains unchanged:

- Check hose, electrical pump connections
- Check whether the sealing ring of the furnace door is clean and undamaged

- **Low vacuum**

The vacuum that is obtained is not sufficient for ceramic firing:

- Check pump power
- Check whether the sealing ring of the furnace door is clean and undamaged

- **Motor too slow**

- **Motor: no signal**

Malfunctions may only be eliminated by service technicians authorized by DeguDent.

7 Statement of Conformity

DeguDent GmbH


P.O. Box 13 64 · 63403 Hanau
GERMANY
www.degudent.com

This is to certify the conformity of the product Cergo® press / high-temperature vacuum furnace with the basic requirements of the following EC guidelines:

1. Low-voltage guideline 73/23/EWG
with 1. Modification 93/68/EWG
2. EMV – Guideline on electromagnetic compatibility 89/336/EWG
with 1. Modification 93/68/EWG

The following harmonized standards were applied:

EN 61010:1990 + A1:1991 + A2:1995 (safety)
EN 50081-1:1992, EN 55022 Class B, EN 61000-3-2(1995) (emission)
EN 50082-1:1992, IEC 801-2(1984), IEC 801-3(1984), IEC 801-4(1988) (noise immunity)

Dr. Alexander Völcker
Head of Research and Development

Torsten Schwafert
Head of Dental Engineering Division

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