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# Degufill<sup>®</sup> KE Gold

CE 0124

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# Degufill® KE Gold

Degufill® KE Gold is a self-curing and light-curing composite resin for permanent bonding of dental prosthetic elements. Its high content of inorganic fillers ensures a tension-free connection between the metal components.

## Indications for use

Composite resin for adhesive bonding of dental restorations.

For dental use only.

## Contraindications

Do not use Degufill® KE Gold in known hypersensitivity to one of its components [particularly (meth-) acrylate or dibenzoyl peroxide].

## Precautions

Sensitization through skin contact is a possibility. Irritating to eyes, respiratory organ and skin. Toxic to aquatic organisms. May cause long-term adverse effects in the aquatic environment. Avoid eye and skin contact. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. After contact with skin, wash immediately with plenty of soap and water and seek medical advice. Wear suitable protective gloves and eye/face protection. Avoid release to the environment. Refer to special instructions/ safety data sheets.

## Adverse effects

Allergic reactions cannot be excluded in predisposed persons.

# Degufill® KE Gold

## Designing the adhesive surfaces

The adhesive surfaces must be designed to be sufficiently stable and placed at regular intervals. The bonding gap should not exceed 0.2 mm.

## Preparing the bonding surfaces

Passively seat the female and male components/primary structure on the cast. Air-abrade the retention and retention housing using alumina (retention and fold housing: 50 µm at 2 bar; non-precious metal housing: 250 µm at 2 bar), then thoroughly clean the components to be adhesively bonded (e.g. using a steam cleaner). Do not touch air-abraded surfaces.

## Adhesion

Dispense equal lengths of the base paste and catalyst (1:1 ratio) onto a mixing block and mix thoroughly with the enclosed plastic spatula. Dosage imprecisions not exceeding ±10% do not adversely affect the processing and performance of Degufill® KE Gold. Avoid air entrapment when mixing.

Homogeneity check: The mixture is homogeneous if the shades of the base paste and catalysts can no longer be distinguished (no more colour streaks).

Apply Degufill® KE Gold to the adhesive surfaces with a spatula, ensuring that all adhesive surfaces are completely covered, without air entrapments or contamination. Place the medal components back on the cast and secure in place. Do not manipulate the metal components before the end of the working time!

## Processing data

Processing time	2 min, including mixing time
Light-curing time	180 s
Self-curing time	20 min
Final hardness	after 12 hours

# Degufill® KE Gold

## **Hardening**

Hardened excess adhesive can be easily removed using a fine-grained tungsten carbide bur.

The final hardness is attained after 12 hours. Do not insert the restoration until the adhesive has completely hardened.

## **Storage**

The shelf life of Degufill® KE Gold is 2 years.

Protect from light. Store at 8–20°C. Before processing, wait 30 minutes until Degufill® KE Gold has attained room temperature, ensuring its full performance.

Do not use after the expiry date.  
Keep out of the reach of children.

## **Presentation**

2.5 g base paste, 2.5 g catalyst, 1 mixing block, 1 mixing spatula  
REF 2401 3900

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DeDuDent GmbH