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GB

**1 Indications for use**

StarLoy L is a CoCr alloy for the fabrication of partial model cast dentures. Due to its composition, StarLoy L is particularly suitable for the laser welding technique.

To be used for dental purposes only.

Storage information: Keep out of children's reach.

**2 Contraindications**

Do not use in case of hypersensitivity to one or more of the alloy components.

**3 Precautions**

When working on dental alloys - especially cobalt-chromium alloys - with mechanical tools, it is generally necessary to work with local aspiration plus face or respiratory protection.

Avoid occlusal and approximal contact of different alloy types.

**4 Side effects**

Allergies to metals contained in the alloy as well as electrochemically conditioned dysesthesia are possible. In individual cases systemic side effects of metals contained in the alloy have been reported.

**5 Technical data**

Composition in mass %

Co	Cr	Mo	Si	Mn	C	N
62,5	30,5	5,0	1,0	0,4	0,3	0,3

StarLoy L is free from nickel and beryllium.

Color:	white	Melting range:	1290 – 1390 °C
Vickers hardness:	375 HV 10	0,2%-yield strength:	710 MPa
Tensile strength:	900 MPa	Modulus of elasticity:	220 GPa
Elongation at fracture:	6%	Density:	8,2 g/cm <sup>3</sup>
Preheating temperature:	1000 °C	Casting temperature:	1500 °C

## 6 Processing instructions

### 6.1 Sprues

2 - 4 sprues – depending on the size of the frame work should be attached to each object. To facilitate the flow and controlled solidification, wax „wires“ approximately 3-4 mm in diameter should be attached to the thickest aspects of the casting piece, ensuring that no tapering sections are created. The length of the sprues must be determined individually by the selection and the position of the sprue former.

### 6.2 Investing

Apply a thin coat of wetting agent to the wax object or spray with it and then blow it dry carefully (do not use compressed air). Invest the pattern with SilberStar investment. Follow the instructions for use of the investing compound. Observe the working instructions of the investment. Fine investing is possible but not required when a fine-grain investing compound is used.

### 6.3 Preheating

Preheat slowly with the suction system activated. Adhere to the heating rates and holding times given by the manufacturer of the respective investment. Remove the wax in the oven at 300 °C. After the final temperature of 1000 °C has been reached, preheat the casting ring for 60 min. In the case of particularly filigree constructions the preheating temperature must be raised to 1050 °C. The holding time must be extended correspondingly, if a larger number of casting rings is used.

### 6.4 Casting

Due to their shape the ingots can be placed easily into the ceramic crucibles. Use ceramic crucibles only.

- In the HF casting machine the casting temperature has been reached as soon as the shadow of the ingots is no longer visible.
- Flame-melting melt with the reducing flame of a propaneoxygen mixture. After exceeding the liquidus temperature and a period of continued heating – 5 -15 s depending on the alloy quantity – begin casting.

### 6.5 Devesting

After casting let the casting ring cool down at room temperature for at least 20 min. Quenching in water is not recommended because of the risk of deformation of the casting. After devesting blast with special abrasive agent.

### 6.6 Finishing

All finishing instruments that are suitable for cast dentures can be used. Use covering varnish to protect filigree elements, e.g. clasps, against abrasion in case of electro-polishing.

### 6.7 Soldering

Use flux e.g. Oxynon® for soldering of StarLoy L to precious metal alloys. The solder that is selected depends on the precious metal alloy to be soldered. For soldering StarLoy L we recommend e.g. Degudent® solder N1W.

### 6.8 Laser welding

Laser welding requires experience in the handling of the welding unit in conjunction with the respective alloy and the dimension of the objects to be welded. A laser welding wire is available. Please observe the manufacturer's recommendations.

### 6.9 Reusability

The fabrication of high-quality dentures requires the use of extremely pure materials. Therefore, to avoid risks only original StarLoy L ingots should be used.

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DeguDent GmbH