DENTSPLY
TRUBYTE

Triad[®] Gel

Visible Light Cure Material

Matériau photopolymérisable

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Material fotopolimerizable

DIRECTIONS FOR USE

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CAUTION: U.S. Federal Law restricts this device to sale by, or on the order of, a licensed dental professional.

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DENTSPLY
TRUBYTE

Manufactured by:
DENTSPLY International Inc.
York, PA 17405-0872
1-800-786-0085

www.trubyte.dentsply.com or www.dentsply.com

Authorized EU Representative
DENTSPLY DeTrey GmbH
Rodenbacher Chaussee 4
D-63457 Hanau-Wolfgang
Germany

Indications for Use:

Triad Gel Material is a visible light cure material indicated for the fabrication of orthodontic and other intra-oral appliances including, but not limited to, retainers, functional appliances, and occlusal splints/night guards. Triad Gel Material is indicated where materials with high flow characteristics are desired, such as when encapsulating and retaining orthodontic wires, mesh, etc.

Triad Gel Material is a no-mix urethane dimethacrylate-based material that contains no methyl methacrylate monomer.

Contraindications:

Triad Gel Material is contraindicated for patients and users with a history of allergic reaction to a) urethane resins; and/or b) methyl methacrylate monomer (contained in Triad VLC Bonding Agent only).

Warnings:

Triad Gel Material and Triad VLC Bonding Agent contain polymerizable monomers which may cause skin sensitization (allergic contact dermatitis) or other allergic reactions in susceptible persons. Wash thoroughly with soap and water after contact. If skin sensitization occurs, discontinue use. If dermatitis or other symptoms persist, seek medical advice.

Precautions:

1. Users with special skin problems, cuts, or abrasions should wear protective gloves.
2. When grinding prosthodontic resins, proper ventilation, masks and vacuum systems should be used.
3. All cure times cited are appropriate for Triad 2000 and Triad II curing units. Consult your Technical Manual or call Technical Support for Triad I cure times.
4. Store Triad Materials at or below 24°C/75°F and away from direct sunlight. Triad Materials may be refrigerated to extend shelf life.
5. Unused Triad Materials should be polymerized prior to disposal.

Adverse Reactions:

1. Allergic contact dermatitis and other allergic reactions may occur in susceptible individuals.
2. Dust will be generated when grinding Triad Materials. Eye, skin and respiratory irritation may occur if appropriate engineering controls are not used.

STEP-BY-STEP INSTRUCTIONS

I. BASIC FUNCTIONAL APPLIANCE:

1. Prepare wires, retainers, retentive mesh, etc. in the normal manner. Allow approximately 0.5 mm space between the metal retention and the cast surface (Figure 1A). Thoroughly clean metal portions of appliances by brushing with methyl methacrylate monomer liquid or ethyl alcohol. Allow to air dry while the cast is being prepared.


Fig. 1A

2. Prepare the cast by applying a thin layer of Triad Model Release Agent (MRA) to only those areas that will come in contact with Triad Gel. Wipe off excess MRA.
3. Position the wires, etc. on the cast. Secure with sticky wax on buccal and facial areas not covered with MRA (Figure 1B).
4. Extrude Triad Gel under and around the wires, then flow around teeth, into interproximal and other areas.


Fig. 1B

Maintain constant pressure on the tube and keep the tip close to the cast. Extruding "rows" of Triad Gel close to each other will help minimize air bubble entrapment. Apply gel to desired thickness, 2-3 mm (Figures 2A, 2B, 2C).


Fig. 2A

- a. Cure the appliance on the cast for four (4) minutes in the Triad VLC Unit.
- b. Remove the appliance from the cast.


Fig. 2B

- c. Paint all Triad material surfaces with Triad Air Barrier Coating (ABC).
- d. Cure the appliance off the cast, tissue side up, for three (3) minutes in the Triad VLC Unit.
5. Remove the appliance from the curing unit and scrub thoroughly with warm water and a brush to remove all traces of Triad ABC.
6. Finish in the normal manner.


Fig. 2C

II. PALATAL EXPANDERS:

1. Prepare the cast as described under "Basic Functional Appliance."
2. With expander positioned, apply Triad Gel to one-half of the palate.
3. Cure 30 seconds in the Triad VLC Unit. Apply Triad ABC to the exposed Triad Gel Material surface, and cure one (1) additional minute on the cast.
4. Ensure that Triad ABC is dry before proceeding.
5. Apply Triad MRA to the dry Triad ABC coated surface, especially the palate edge what will contact the second half of the expander.
6. Apply Triad Gel Material to the rest of the palate to form the second half of the expander appliance.

7. Complete the appliance as previously described in "Basic Functional Appliance."

B. Occlusal Splints can be made easily on the cast with Triad Gel following conventional design procedures.

C. Patient Identification Procedures:

1. For new appliances:

- a. Apply a thin layer of Triad Gel to the prepared cast surface in the desired location and cure for 30 seconds.
- b. Moisten "onion skin" paper containing the patient's typed identification data with heat-cured methyl methacrylate monomer (MMA), and place the paper on the thin, partially-cured Triad Gel Material surface.
- c. Continue the Triad Gel Material application as previously described in "Basic Functional Appliance."

2. For existing appliances:

- a. Create a shallow well in a non-tissue bearing surface of the appliance.
- b. Moisten "onion skin" paper containing the patient's identification data using heat-cured MMA. Position the paper in the prepared area.
- c. Cover with a layer of Triad Gel Material and cure for 30 seconds.

- d. Cover with Triad ABC and cure an additional 3 minutes.
- e. After curing, scrub gently with warm water and a brush to remove all traces of Triad ABC.
- f. If needed, trim excess material and polish in a normal manner.

D. Additions, Modifications, and Repairs

Additions of Triad Gel Material can be made at any time during construction before Triad ABC has been applied. After application of the Triad ABC, Triad Gel Material can be added to the cured surface using the following procedure:

1. Remove the Triad ABC by scrubbing thoroughly with a brush and warm water.
2. Freshen the Triad material surface by lightly grinding the area where the addition is to be made.
3. Apply a thin layer of Triad VLC Bonding Agent to the prepared surface, and allow to bench set for 2 minutes.

4. Cure the bonding agent in the Triad VLC Unit for 2 minutes.
5. Add Triad Gel Material and process as previously described in "Basic Functional Appliance."

*Heat-cured methyl methacrylate monomer (Lucitone Liquid) can be substituted for Triad VLC Bonding Agent. Refer to the repair instructions in the Triad VLC System Technique Manual and follow these procedures through completion.

MATERIAL SAFETY DATA SHEET

FORM NO.: 100

SECTION I - GENERAL INFORMATION

- 1.1 Product Trade Name: Triad® Gel
- 1.2 Part (Item) Number: 88859, 88860, 88750, 88751, 88762, 90130
- 1.3 Chemical Name: N.A.
- 1.4 Synonyms: N.A.
- 1.5 Division Name: DENTSPLY Trubyte
Address: 570 W. College Ave.
City, State, Zip: York, PA 17405-0872
- 1.6 Phone No. for Info.: 717-845-7511
Dentsply

"FOR CHEMICAL EMERGENCY"
SPILL, LEAK, FIRE, EXPOSURE OR ACCIDENT
CALL CHEMTREC - DAY OR NIGHT

- 1.7 Emergency Phone No.: 800-424-9300
- 1.8 Date: 6/25/97
- 1.9 Dentsply Contact: Scott E. Shaffer
- 1.10 Signature: *Scott E. Shaffer*
- 1.11 Formula: N.A.
- 1.12 National Stock No.: N.A.
- 1.13 Specification No.: N.A.
- 1.14 Chemical Family: Acrylic
- 1.15 Type/Grade/Class: N.A.
- 1.16 NRC/EPA Number: N.A.

SECTION II - HAZARDOUS INGREDIENTS

CAS Number	Chemical Name (Ingredients)	%	TLV	OSHA PEL
NONE	NONE	NONE	NONE	NONE

SECTION III - PHYSICAL PROPERTIES

- 3.1 Boiling Point: N.A.
- 3.2 Vapor Pressure: <1 mm Hg
- 3.3 Vapor Density: >1
- 3.4 Critical Temp: N.A.
- 3.5 Decomposition Temp: N.A.
- 3.6 Melting/Freezing Pt: <-40°C(-40°F)
- 3.7 Solubility in Water: Negligible
- 3.8 Critical Pressure: N.A.
- 3.9 Perm. Exp. Limit: N.A.
- 3.10 Appearance and Odor: Clear viscous gel with a slight ester odor. Tint added if red, blue or pink.
- 3.11 Specific Gravity: 1.2
- 3.12 pH: N.A.
- 3.13 Evaporation Rate: N.A.
- 3.14 Viscosity: 90,000 cps at 40°C(104°F)
- 3.15 % Volatile by Vol.: Negligible
- 3.16 Magnetism: N.A.
- 3.17 Autoignition Temp: N.A.
- 3.18 Corrosion Rate: N.A.
- 3.19 Molecular Weight: N.A.

SECTION IV - FIRE AND EXPLOSION HAZARD DATA

- 4.1 Flash Point(Method Used): N.A.
- 4.2 Flammable (Explosive) Limits:
LEL: N.A. UEL: N.A.
- 4.3 Extinguishing Media: Water fog, foam, carbon dioxide, water spray, dry chemical.
- 4.4 Special Fire Fighting Procedures and Equipment: Firefighters should wear protective clothing and self contained breathing apparatus if large amounts are involved.
- 4.5 Unusual Fire and Explosion Hazards: Irritating gases may be released.

Triad® Gel, cont'd
SECTION V - HEALTH HAZARD DATA

- 5.1 Effects of Acute and Chronic Over Exposure: Material is nontoxic. Skin irritation is possible in sensitive individuals. Contact dermatitis (skin sensitization) can develop in susceptible individuals.
- 5.2 Emergency and First Aid Procedures: Eyes: Flush immediately with water for a minimum of 15 minutes and get medical attention. Skin Irritation: Wash with soap and water, remove contaminated clothing. Sensitization: Avoid further contact.

SECTION VI - REACTIVITY DATA

- 6.1 Unstable: Stable: X
- 6.2 Conditions to Avoid (Stability): High temperatures, direct sunlight.
- 6.3 Incompatibility (Materials to Avoid): Oxidizing agents.
- 6.4 Hazardous Decomposition Products: Acrid fumes and smoke may result during initial stages of a fire.
- 6.5 Hazardous Polymerization: May Occur: Will Not Occur: X
- 6.6 Conditions to Avoid (Polymerization): High temperatures (polymerizes to form a hard plastic in 2 hours at 100°C (212°F) with a low exotherm.

SECTION VII - SPECIAL PRECAUTIONS AND SPILL OR LEAK PROCEDURES

- 7.1 Handling and Storage Precautions: Do not store near high temperature or ignition sources. Refrigeration prolongs shelf life. Store away from food and beverages.
- 7.2 Other Precautions: This material should be handled and used only in a manner consistent with the product instructions. Material is sensitive to visible light and will polymerize to form a hard plastic.
- 7.3 Steps to be Taken in Case Material is Released or Spilled: No special procedures required; follow waste disposal procedure in 7.4.
- 7.4 Waste Disposal Procedures: Spread out materials to expose surfaces to sunlight or artificial light until solid, then dispose as a nonhazardous waste.
- 7.5 Neutralizing Agent: N.A.

SECTION VIII - OCCUPATIONAL PROTECTION MEASURES

- 8.1 Respiratory Protection: Use dust mask while grinding polymerized (cured) material to avoid breathing dust.
- 8.2 Ventilation: Use local exhaust when grinding polymerized (solid or cured) material.
- 8.3 Protective Gloves: Users with special skin problems, cuts, or abrasions should wear protective gloves (latex or PVC).
- 8.4 Eye Protection: Use safety glasses while grinding.
- 8.5 Other Personal Protective Equipment: None necessary.

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