**DENTSPLY**  
TRUBYTE

## Characterized Lucitone®

Denture Base Resin  
Résine pour plaque de prothèse  
Resina per protesi dentali  
Prothesenkunststoff  
Caracterizada Resina para bases de prótesis

**DIRECTIONS FOR USE**  
**CONDITIONS D'EMPLOI**  
**ISTRUZIONI PER L'USO**  
**VERARBEITUNGSANLEITUNG**  
**INSTRUCCIONES DE EMPLEO**

**CAUTION:** U.S. Federal Law restricts this device to sale by, or on the order of, a licensed dental professional.

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**DENTSPLY**  
TRUBYTE

Manufactured by:  
DENTSPLY International Inc.  
York, PA 17405-0872  
1-800-786-0085

[www.trubyte.dentsply.com](http://www.trubyte.dentsply.com) or [www.dentsply.com](http://www.dentsply.com)

Authorized EU Representative  
DENTSPLY DeTrey GmbH  
Rodenbacher Chaussee 4  
D-63457 Hanau-Wolfgang  
Germany

Printed in U.S.A.  
586005 Rev. D (8/02)

© 1986, 1997, 1999, 2000, 2002 DENTSPLY International Inc.  
All rights reserved.

**Indications for Use:**

CHARACTERIZED LUCITONE Denture Base Resin material includes copolymers of methacrylate and long compatible fibers for the simulation of the vascular system. The resin is indicated for the fabrication of prosthodontic restorations.

**Contraindications:**

1. CHARACTERIZED LUCITONE Denture Base Resin is contraindicated for patients and users with a history of allergic reaction to methyl methacrylate monomer.

**Warnings:**

1. CHARACTERIZED LUCITONE Denture Base Resin contains polymerizable monomers which may cause skin sensitization (allergic contact dermatitis) or other allergic reactions in susceptible people. Wash thoroughly with soap and water after contact. If skin sensitization occurs, discontinue use. If dermatitis or other symptoms persist, seek medical advice.

2. Avoid inhalation or ingestion. High vapor concentration can induce headache, nausea, and irritation of eyes and respiratory system. Liquid contact with eyes may cause possible corneal damage. Excessive long-term exposure may be associated with other more serious health effects. Monitor air quality per OSHA standards.

**Inhalation:** Move subject to fresh air. Give oxygen or artificial respiration as required.

**Ingestion:** Contact your regional Poison Control Center immediately.

**Eye Contact:** Flush eyes promptly with copious amounts of water for 15 minutes; consult a physician. Wash skin with soap and water.

**Precautions:**

1. Users with special skin problems, cuts, or abrasions should wear protective gloves.

2. When grinding prosthodontic resins, proper ventilation, masks, and vacuum systems should be used.
3. Store at 60°- 80°F, away from moisture and direct sunlight. The liquid contains unsaturated monomers which may polymerize prematurely if stored at excessively high temperature and/or in sunlight. The powder contains an organic polymerization initiator which may degrade if stored at excessively high temperature.
4. Lucitone Liquid and CHARACTERIZED LUCITONE Denture Base uncured resin are hazardous materials. Dispose of in accordance with Federal, State, and Local regulations.
5. Lucitone Liquid contains methyl methacrylate monomer, a flammable liquid with a flash point of 10°C (50°F). Keep away from heat, sparks, and open flame.
6. Use liquid in a well-ventilated area. Replace cap when not in use.
7. Do not pack denture resin after work time has elapsed or material has become rubbery.
8. Use of wax solvents is not recommended. Residual solvent may cause poor bonding of teeth to acrylic base.

**Adverse Reactions:**

1. Corneal damage, headache, nausea, and vomiting may occur with exposure to methyl methacrylate monomer. (See Warnings and Precautions statements.)
2. Allergic contact dermatitis and other allergic reactions may occur in susceptible individuals. Residual monomer in fully cured materials can be minimized by soaking the cured prosthesis in warm water for several days.
3. Particulates will be generated when grinding acrylic resins. Eye, skin, and respiratory irritation may occur if appropriate engineering controls are not used.

**PROPERTIES AND RECOMMENDED PROCEDURES**

The following table describes typical properties plus recommended temperatures and other important processing information. Further amplification of this condensed information is described later in the text.

Type and Class (according to ADA Spec. 12)	Type I Class 1
Storage temperature for powder and liquid	60°-80°F
Powder/Liquid ratio	21 gm (30 cc)/10 ml
Mixing time (time to wet all particles)	15-30 seconds
Time to reach packing plasticity @ 23±1°C	3 minutes
Working time	60 minutes
Material used to prepare mold	gypsum
Temperature of mould when packing	approx. 110°F
Recommended cure time and temperature	
1st stage	1½ hours @ 163°F
2nd stage	½ hour
Alternate cure time and temperature	9 hours @ 163°F
Method of cooling flask, time and temperature	
1st stage	½ hour in air @ 60°-80°F
2nd stage	¼ hour in water @ 60°-80°F

**STEP-BY-STEP INSTRUCTIONS**

**FLASKING**

Use conventional dental compression moulding methods and gypsum materials for flasking.

**WAX ELIMINATION**


Soften wax in boiling water for approximately 6 minutes. Separate flask and remove wax by flushing with a solution of boiling water to which a detergent has been added. Use clean boiling water for the final flush. DO NOT USE WAX SOLVENTS. Be sure case is free of wax and grease.

**TINFOIL SUBSTITUTE**

Apply AL-COTE® Tinfoil Substitute liberally to areas of the warm mould (approx. 120°F) that will contact the resin. Remove excess AL-COTE film with a dry brush. (AL-COTE does not form a film on teeth.) The film of AL-COTE should be thoroughly dry (approx. 3 minutes following application) before the resin is packed.

**CUSTOM PACKING**


1. The temperature of the flask at the time of packing should be slightly warm (approx. 110°F).
2. Use Gingival Tone Powder to simulate tissue tension in gingival area. Apply 2-3 mm mounds above necks of posterior teeth on one side. Tilt the flask for convenient use.


3. Using a small brush, moisten the powder with liquid. Repeat on the lingual side. On the labial surface above cuspids,

extend gingival shading to 4 mm. Powder mounds above centrals and labials should be 2-3 mm. KEEP POWDER OUT OF INTERPROXIMAL SPACES.

Immediately after moistening the powder mounds with liquid (step 2), lightly dust Gingival Tone Powder over the entire area to control the flow of gingival shading and to diffuse its line of demarcation.


4. The correct proportion for mixing of denture resin is 1 pouch of powder (21 gm or 30 cc) to 10 ml of liquid.
5. First measure 10 ml of liquid and pour it into the mixing jar. Then add powder from the pouch while stirring constantly. Thoroughly mix to evenly distribute the fibers in the mass (15-30 seconds). Cover the jar during the gelation stage.


6. Approximately 3 minutes from the start of the mix (at  $73^{\circ} \pm 2^{\circ}\text{F}$ ), the resin dough will separate cleanly from the walls of the

3


mixing jar. Roll the soft dough into the shape of a cylinder (approx. 1" diameter) to properly orient the fibers.


7. Place the roll of dough between moistened cellophane. Flatten dough into a  $\frac{1}{8}$ " thick slab.


8. Cut off a rounded end for the palatal section of the denture. Cut strips for buccal and labial flanges.

9. Press the strips into position in the mould. Note that the fibers should extend vertically above the teeth, simulating vascular appearance of natural tissue.

10. Place the rounded end of the slab into the palatal area. Fibers should extend distally. Cover the resin with polyethylene or moistened cellophane. Close the flask slowly. Trial pack and cure.

NOTE: To prevent displacement of the Gingival Tone Resin, underpack the denture base resin on initial closure. Close the flask slowly. Additional resin should be added during subsequent trial packing to completely fill the mould. The closed flask should remain under pressure for approximately 30 minutes before the heat-cure procedure is started.

#### WORKING CHARACTERISTICS

The work time of CHARACTERIZED LUCITONE resin is approximately 60 minutes in room temperature of  $73^{\circ} \pm 2^{\circ}\text{F}$ . Unused mixed material (sealed in a closed jar) can be safely refrigerated several days and still remain packing soft. Upon removal from refrigerator, the jar should remain closed until the material inside has reached room temperature. Otherwise, moisture from the warm air will condense on the resin causing objectionable blushing in the denture.

#### CURING

Submerge closed flask (locked by compress or spring clamp) in water at  $163^{\circ} \pm 2^{\circ}\text{F}$  for  $1\frac{1}{2}$  hours. Follow by  $\frac{1}{2}$  hour in boiling water. A periodic check of water bath temperature with an accurate thermometer is recommended. Alternate cure: minimum of 9 hours in water bath at  $163^{\circ} \pm 2^{\circ}\text{F}$ . Boiling is not necessary for the alternate use method.

#### CURING OF EXTREMELY THICK CASES

Extremely thick cases require special handling to prevent porosity and excessive surface shrinkage. The following procedure is strongly recommended:

1. Pack material when it reaches a heavy, putty-like consistency (approx. 30 minutes at  $73^{\circ} \pm 2^{\circ}\text{F}$ ).
2. Condense well during packing by closing flask slowly and trial packing at least 3 times.
3. The closed flask should remain under pressure for 30 minutes prior to curing.
4. Use the alternate cure method as specified above.

#### COOLING

Bench cool the clamped flask at room temperature for 30 minutes. Then immerse in

4

cool water (60°-80°F) for 15 minutes before deflasking.

**FINISHING AND POLISHING**

CHARACTERIZED LUCITONE resin should be finished and polished in the usual manner.

**REPAIRING AND RELINING**

Relines or repairs to the complete denture may be made with self-curing DENTSPLY®

REPAIR MATERIAL in accordance with directions provided.

If desired, repairing or relining can also be done by the usual technique using CHARACTERIZED LUCITONE resin itself. Cure in water at 163°F for 9 hours—do not boil.

SHADE RECOMMENDATIONS	
REPAIR MATERIAL	DENTURE BASE RESINS
PINK FIBERED	CHARACTERIZED LUCITONE Original, Light, Light Reddish Pink, Dark Pink
LIGHT FIBERED	CHARACTERIZED LUCITONE Lt.

MATERIAL SAFETY DATA SHEET

FORM NO.: 169 SECTION I - GENERAL INFORMATION  
 1.1 Product Trade Name: Lucitone® Characterized Denture Base Resin  
 1.2 Part (Item) Number: 686001-002, 686004-006, 686031, 686034, 686042-046, 686054-056, 686120, 686123-124, 686130, 686133-134, 686140, 686143-144, 686371-373, 686376, 686386  
 1.3 Chemical Name: Polymethylmethacrylate  
 1.4 Synonyms:  
 1.5 Division Name: Dentsply/Trubyte Division  
 Address: 570 W. College Ave.  
 City, State, Zip: York, PA 17405-0872  
 1.6 Phone No. for Info.: 717-845-7511  
 Dentsply

"FOR CHEMICAL EMERGENCY"  
 SPILL, LEAK, FIRE, EXPOSURE OR ACCIDENT  
 CALL CHEMTREC - DAY OR NIGHT

1.7 Emergency Phone No.: 800-424-9300 1.12 National Stock No.: N.A.  
 Chemtrec 1.13 Specification No.: N.A.  
 1.8 Date: 9/22/00 1.14 Chemical Family: Acrylic  
 1.9 Dentsply Contact: Scott E. Shaffer Polymer  
 1.10 Signature: *Scott E. Shaffer* 1.15 Type/Grade/Class: N.A.  
 1.16 NRC/EPA Number: N.A.  
 1.11 Formula: (CH<sub>2</sub>(CH<sub>3</sub>)COOCH<sub>2</sub>)<sub>n</sub>

SECTION II - HAZARDOUS INGREDIENTS

2.1 CAS Number	Chemical Name (Ingredients)	%	TLV mg/m <sup>3</sup>	OSHA PEL
94-36-0	Benzoyl Peroxide	<.2	TWA	5mg/m <sup>3</sup>

SECTION III - PHYSICAL PROPERTIES

3.1 Boiling Point: N.A. 3.11 Specific Gravity: N.E.  
 3.2 Vapor Pressure: N.A. 3.12 pH: N.A.  
 3.3 Vapor Density: N.A. 3.13 Evaporation Rate: N.A.  
 3.4 Critical Temp: N.E. 3.14 Viscosity: N.A.  
 3.5 Decomposition Temp.: 200°C/392°F 3.15 % Volatile by Vol.: 0  
 3.6 Melting/Freezing Pt.: N.A. 3.16 Magnetism: N.A.  
 3.7 Solubility in Water: Not Soluble 3.17 Autoignition Temp.: N.E.  
 3.8 Critical Pressure: N.A. 3.18 Corrosion Rate: N.A.  
 3.9 Perm. Exp. Limit: N.A. 3.19 Molecular Weight: N.A.  
 3.10 Appearance and Odor: Pink free flowing powder with faint odor in bulk

SECTION IV - FIRE AND EXPLOSION HAZARD DATA

4.1 Flash Point (Method Used): Approximately 300°C (572°F) (TCC)  
 4.2 Flammable (Explosive) Limits:  
 LEL: 20g/m<sup>3</sup> (dust cloud) UEL: N.A.  
 4.3 Extinguishing Media: Water spray, carbon dioxide, dry chemical, water foam.  
 4.4 Special Fire Fighting Procedures and Equipment: Avoid extinguishing methods which may generate dust clouds. Water stream can disperse dust in air, producing a fire hazard and possible explosion hazard if exposed to ignition source.  
 4.5 Unusual Fire and Explosion Hazards: Polymer dust is combustible but not easily ignited. The explosive limits of the polymer particles suspended in air are approximately those of coal dust. Fire fighters should avoid extinguishing methods which generate dust clouds. Fire Fighters should wear self-contained breathing apparatus if fighting large fires.

Lucitone® Characterized Denture Base Resin, cont'd  
SECTION V - HEALTH HAZARD DATA

5.1 Effects of Acute and Chronic Over Exposure: It is not known to cause significant health problems. It is considered an inert or nuisance dust. Avoid inhalation of dust. Keep dust out of eyes to prevent possible irritation.

5.2 Emergency and First Aid Procedures: INHALATION: Remove to fresh air. Get medical help if discomfort persists. EYES: Flush with water for 15 minutes, including under eyelids. SKIN: Wash with soap and water. INGESTION: Rinse mouth with water. Do not induce vomiting. Call doctor if amount was large.

SECTION VI - REACTIVITY DATA

6.1 Unstable: Stable: X

6.2 Conditions to Avoid (Stability): Avoid heat above 250°C (392°F)

6.3 Incompatibility (Materials to Avoid): Strong oxidizing agents

6.4 Hazardous Decomposition Products: Methylmethacrylate and Carbon Monoxide

6.5 Hazardous Polymerization: May Occur: Will Not Occur: X

6.6 Conditions to Avoid (Polymerization): Will not polymerize.

SECTION VII - SPECIAL PRECAUTIONS AND SPILL OR LEAK PROCEDURES

7.1 Handling and Storage Precautions: Store in cool dry place. Keep container closed to prevent water absorption and contamination.

7.2 Other Precautions: Observe normal warehousing procedures. Store away from food and beverages.

7.3 Steps to be Taken in Case Material is Released or Spilled: Sweep up to avoid slipping hazard. Keep airborne particulate at a minimum when cleaning up spills.

7.4 Waste Disposal Procedures: May be disposed of in landfill or incinerated. Follow Federal, State, and local regulations for disposal.

7.5 Neutralizing Agent: N.A.

SECTION VIII - OCCUPATIONAL PROTECTION MEASURES

8.1 Respiratory Protection: Nuisance dust type if needed. An approved dust mask should be worn if exposure at high levels of dust are likely.

8.2 Ventilation: Local exhaust at processing equipment

8.3 Protective Gloves: Not required.

8.4 Eye Protection: Safety glasses

8.5 Other Personal Protective Equipment: None required.

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