**DENTSPLY**  
TRUBYTE

# Triad<sup>®</sup> TranSheet<sup>™</sup>

Visible Light Cure Material

Matériau photopolymérisable

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Material fotopolimerizable

## DIRECTIONS FOR USE

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**CAUTION:** U.S. Federal Law restricts this device to sale by, or on the order of, a licensed dental professional.

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**DENTSPLY**  
TRUBYTE

Manufactured by:  
DENTSPLY International Inc.  
York, PA 17405-0872  
1-800-786-0085

Authorized EU Representative  
DENTSPLY DeTrey GmbH  
Rodenbacher Chaussee 4  
D-63457 Hanau-Wolfgang  
Germany

[www.trubyte.dentsply.com](http://www.trubyte.dentsply.com) or [www.dentsply.com](http://www.dentsply.com)

**Indications for Use:**

Triad TranSheet is a visible light cure material indicated for the fabrication of removable orthodontic and dental prosthetic appliances, including but not limited to retainers, functional appliances, and occlusal splint/nightguards.

Triad TranSheet is a urethane dimethacrylate-based material which contains no methyl methacrylate monomer.

**Contraindications:**

Triad TranSheet Material is contraindicated for patients and users with a history of allergic reaction to a) urethane resins; and/or b) methyl methacrylate monomer (contained in Triad VLC Bonding Agent only).

**Warnings:**

Triad TranSheet Material and Triad VLC Bonding Agent contain polymerizable monomers which may cause skin sensitization (allergic contact dermatitis) or other allergic reactions in susceptible persons. Wash thoroughly with soap and water after contact. If skin sensitization occurs, discontinue use. If dermatitis or other symptoms persist, seek medical advice.

**Precautions:**

- Users with special skin problems, cuts, or abrasions should wear protective gloves.
- When grinding prosthodontic resins, proper ventilation, masks and vacuum systems should be used.
- All cure times cited are appropriate for Triad 2000 and Triad II curing units. Consult your Technique Manual or call Technical Support for Triad I curing times.
- Store Triad Materials at or below 24°C/75°F and away from direct sunlight. Triad Materials may be refrigerated to extend shelf life.
- Unused Triad Materials should be polymerized prior to disposal.
- Do not use Triad Materials that crack or craze when handled.

**Adverse Reactions:**

- Allergic contact dermatitis and other allergic reactions may occur in susceptible individuals.
- Dust will be generated when grinding Triad Materials. Eye, skin and respiratory irritation may occur if appropriate engineering controls are not used.

**STEP-BY-STEP INSTRUCTIONS**

**ORTHODONTIC RETAINER**

- Apply a coating of Triad Model Release Agent (MRA) to the cast.
- Position and secure wires on the cast, at least .05 mm above the palate.
- Adapt a small amount of Triad TranSheet Material under and/or around wires. (Note: Triad Gel may also be used to cover wires.) Adapt Triad TranSheet Material to the cast with finger pressure and the Triad Modeling Tool.
- Trim to the desired shape with a sharp blade.
- Place in the Triad Light Curing Unit for a four (4) minute preliminary cure.
- Following the preliminary cure, remove the appliance from the cast and apply Triad Air Barrier Coating (ABC) to all surfaces of the VLC materials. Return the appliance to the curing unit, tissue side up, and cure for six (6) minutes.
- Gently scrub the appliance with warm water and a brush to remove all traces of Triad ABC.
- Trim and polish as required.

**IMPLANT PLACEMENT STENT FOR A PARTIALLY EDENTULOUS ARCH**

- Set and wax appropriate artificial teeth in the edentulous area.
- Apply Triad MRA to the cast, including teeth and wax.
- Drape Triad TranSheet Material

over the artificial teeth, wax, and adjacent natural teeth which are modeled on the cast. The material thickness should be approximately 4 mm.

- Cure in the Triad Light Curing Unit for four (4) minutes.
- Remove the stent from the cast and eliminate any remaining wax with boiling water.
- Paint Triad ABC on all surfaces of the stent and cure an additional six (6) minutes.
- Remove remaining traces of Triad ABC with a brush and warm water.
- Trim excess material. Maintain adequate coverage of adjacent teeth and tissue in order to correctly orient the stent.
- Drill holes of the same diameter as the laboratory screws in the ideal implant locations.
- Complete finishing and polishing procedures.

**PLACEMENT STENT FOR IMPLANTS IN A TOTALLY EDENTULOUS ARCH**

- Mark the cast at the crest of the ridge. Using a contrasting color, mark the implant location.
- Drill parallel holes in the cast, the same diameter as the threaded portion of the laboratory screw, at the implant locations.
- Place the laboratory screws in the holes and secure with sticky wax.
- Paint Triad MRA on the cast and screws.
- Adapt Triad TranSheet Material to the edentulous area, including screws, adding material as necessary to obtain a thickness of approximately 4 mm. Cure in the Triad Light Curing Unit for five (5) minutes.
- Remove the stent from the cast and apply Triad ABC to all surfaces. Return to the Triad Curing Unit and cure for an additional five (5) minutes, tissue side up.

- After curing, remove the screws and gently scrub the stent with warm water and a brush to remove remaining traces of Triad ABC. Check to ensure that holes are free of any interferences.
- Trim and polish the stent to a normal denture outline.

**ORTHOTIC OR SPLINT**

- Pour casts. Trim and mount casts in an articulator utilizing the registration provided by the dental office.
- Wax out any undesirable undercuts.
- Apply Triad MRA to both casts.
- Apply Triad TranSheet Material to the cast indicated in a thickness that will allow indentations to be created by the opposing teeth.
- Close the articulator, embedding the teeth in the uncured Triad TranSheet Material. Lateral excursions may be accomplished at this time by opening and closing the articulator as the upper bow is moved to eccentric positions.
- Cure for four (4) minutes in the Triad Light Curing Unit.
- Remove the appliance from the cast, apply Triad ABC to all surfaces and return to the Triad unit for four (4) minutes.
- Gently scrub the appliance with a brush and warm water to remove traces of Triad ABC.
- Trim excess material. The incisal, gingival and occlusal areas should be trimmed to the design requirement of the clinical situation.
- After contouring, complete polishing procedures in the normal manner.

Triad I Cure Times			
Cure	Ortho Ret.	Placement Stent	Splint
Initial Cure	10	10	10
Final Cure	6	6	6

## MATERIAL SAFETY DATA SHEET

FORM NO.: 86

## SECTION I - GENERAL INFORMATION

1.1 Product Trade Name: Triad® TransSheet™ Light Cure Material  
 1.2 Part (Item) Number: 89270, 89271, 89272, 89273, 89274, 89275, 89276, 90130  
 1.3 Chemical Name: N.A.  
 1.4 Synonyms: N.A.  
 1.5 Division Name: Dentsply/Trubyte Division  
 Address: 570 W. College Ave.  
 City, State, Zip: York, PA 17405-0872  
 1.6 Phone No. for Info.: 717-845-7511  
 Dentsply

"FOR CHEMICAL EMERGENCY"  
 SPILL, LEAK, FIRE, EXPOSURE OR ACCIDENT  
 CALL CHEMTREC - DAY OR NIGHT

1.7 Emergency Phone No.: 800-424-9300 1.12 National Stock No.: N.A.  
 Chemtrec 1.13 Specification No.: N.A.  
 1.8 Date: 6/25/97 1.14 Chemical Family: Urethane  
 1.9 Dentsply Contact: Scott E. Shaffer Acrylate  
 1.10 Signature: *Scott E. Shaffer* 1.15 Type/Grade/Class: N.A.  
 1.11 Formula: N.A. 1.16 NRC/EPA Number: N.A.

## SECTION II - HAZARDOUS INGREDIENTS

2.1 CAS Number	Chemical Name (Ingredients)	%	TLV	OSHA PEL
NONE	NONE	NONE	NONE	NONE

## SECTION III - PHYSICAL PROPERTIES

3.1 Boiling Point: N.A. (Polymerizes)	3.11 Specific Gravity: 1.2
3.2 Vapor Pressure: <1 mm Hg	3.12 pH: N.A.
3.3 Vapor Density: >1	3.13 Evaporation Rate: N.A.
3.4 Critical Temp: N.A.	3.14 Viscosity: N.A.
3.5 Decomposition Temp: N.A.	3.15 % Volatile by Vol.: Negligible
3.6 Melting/Freezing Pt: <-40°C (-40°F)	3.16 Magnetism: N.A.
3.7 Solubility in Water: Negligible	3.17 Autoignition Temp: N.A.
3.8 Critical Pressure: N.A.	3.18 Corrosion Rate: N.A.
3.9 Perm. Exp. Limit: N.A.	3.19 Molecular Weight: N.A.
3.10 Appearance and Odor: Colorless or clear pink doughy material in the uncured state. Colorless or clear pink hard plastic when cured.	

## SECTION IV - FIRE AND EXPLOSION HAZARD DATA

4.1 Flash Point (Method Used): N.A.  
 4.2 Flammable (Explosive) Limits:  
 LEL: N.A. UEL: N.A.  
 4.3 Extinguishing Media: Water fog, foam, carbon dioxide.  
 4.4 Special Fire Fighting Procedures and Equipment: Fireman should wear protective clothing and self-contained respirators if large amounts are involved.  
 4.5 Unusual Fire and Explosion Hazards: Irritating gases may be released.

Triad® TransSheet™ Light Cure Material, cont'd  
SECTION V - HEALTH HAZARD DATA

5.1 Effects of Acute and Chronic Over Exposure: Material is non-toxic. Skin irritation is possible in sensitive individuals. Contact dermatitis (skin sensitization) can develop in sensitive individuals.  
 5.2 Emergency and First Aid Procedures: Eyes: Flush immediately with water for a minimum of 15 minutes and get medical attention. Skin irritation: Wash with soap and water, remove contaminated clothing. Sensitization: Avoid further contact.

## SECTION VI - REACTIVITY DATA

6.1 Unstable: Stable: X  
 6.2 Conditions to Avoid (Stability): Heat and ignition sources. Material is sensitive to visible light and will harden during exposure to sunlight.  
 6.3 Incompatibility (Materials to Avoid): Oxidizing and reducing chemicals.  
 6.4 Hazardous Decomposition Products: None.  
 6.5 Hazardous Polymerization: May Occur: Will Not Occur: X  
 6.6 Conditions to Avoid (Polymerization): Heat (polymerizes in two hours at 100°C (212°F) with a low exotherm).

## SECTION VII - SPECIAL PRECAUTIONS AND SPILL OR LEAK PROCEDURES

7.1 Handling and Storage Precautions: Do not store near high temperatures or ignition sources. Refrigerated storage prolongs shelf life. Store away from food and beverage.  
 7.2 Other Precautions: This material should be handled and used only in a manner consistent with the product instructions. Material is sensitive to visible light and will polymerize to form a hard plastic. Steps to be Taken in Case Material is Released or Spilled: No special procedures required: Follow waste disposal procedure in 7.4.  
 7.4 Waste Disposal Procedures: Spread out materials to expose surfaces to sunlight or artificial light until solid. Dispose of as non hazardous waste.  
 7.5 Neutralizing Agent: N.A.

## SECTION VIII - OCCUPATIONAL PROTECTION MEASURES

8.1 Respiratory Protection: Use dust mask while grinding to avoid breathing dust.  
 8.2 Ventilation: Use local exhaust when grinding polymerized (solid or cured) material.  
 8.3 Protective Gloves: Users with special skin problems, cuts or abrasions should wear protective gloves (latex or PVC).  
 8.4 Eye Protection: Use safety glasses while grinding.  
 8.5 Other Personal Protective Equipment: None necessary.

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