**DENSPLY**  
TRUBYTE

**Lucitone<sup>®</sup>**  
***Fas-Por<sup>™</sup>* +**  
Pourable Denture Base  
Plaque base coulable  
Base della protesi versabile  
Gießfähiger flüssiger  
Prothesenkunststoff  
Resina colable

**DIRECTIONS FOR USE**

**MODE D'EMPLOI**

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**Caution:** U.S. Federal Law restricts this device to sale by, or on the order of, a licensed dental professional.

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**Indications for use:** Lucitone® Fas-Por™+ Pourable Denture Base is an auto-polymerizing powder and liquid denture base material based on polymethyl methacrylate. This product is intended for the fabrication of complete or partial removable dentures or other removable appliances. The Lucitone Fas-Por+ Pourable Denture Base can also be used for repairs and relines.

**Contraindications:**

1. Lucitone Fas-Por+ Pourable Denture Base is contraindicated for patients and users with a history of allergic reaction to methyl methacrylate monomer.

**Warnings:**

1. Lucitone Fas-Por+ Pourable Denture Base contains polymerizable monomers which may cause skin sensitization (allergic contact dermatitis) or other allergic reactions in susceptible persons. Wash thoroughly with soap and water after contact. If dermatitis or other symptoms persist, seek medical advice.
2. Avoid inhalation or ingestion. High vapor concentration can induce headache, nausea, and irritation of eyes and respiratory system. Liquid contact with eyes may cause possible corneal damage. Excessive long-term exposure may be associated with other more serious health effects. Monitor air quality per OSHA standards.  
**Inhalation:** Move subject to fresh air. Give oxygen or artificial respiration as required.  
**Ingestion:** Contact your regional Poison Control Center immediately.  
**Eye Contact:** Flush eyes promptly with copious amounts of water for 15 minutes, and consult a physician. Wash skin with soap and water.

**Precautions:**

1. Users with special skin problems, cuts, or abrasions should wear protective gloves.
2. When grinding prosthodontic resins, proper ventilation, masks, and vacuum systems should be used.
3. Store at 60°F-80°F (16°C-27°C), away from moisture and direct sunlight. The liquid contains unsaturated monomers which may polymerize prematurely if stored at excessively high temperature and/or in sunlight. The powder contains an organic polymerization initiator which may de-

- grade if stored at excessively high temperature.
4. Lucitone Fas-Por+ Pourable Denture Base Liquid and Lucitone Fas-Por+ Pourable Denture Base uncured resin are hazardous materials. Dispose of in accordance with Federal, State, and Local regulations.
  5. Lucitone Fas-Por+ Pourable Denture Base Liquid contains methyl methacrylate monomer, a flammable liquid with a flash point of 50°F (10°C). Keep away from heat, sparks, and open flame.
  6. Use liquid in a well ventilated area. Replace cap when not in use.
  7. Do not use detergents to clear residual wax from casts. Detergent residue may impact the color uniformity.
  8. Lucitone Fas-Por+ Pourable Denture Base resin's setup time is quick (approximately 1 minute). **After mixing powder and liquid, IMMEDIATELY pour into flask. Do not pour the denture resin after the material has thickened.**
  9. **DO NOT ALLOW** Lucitone Fas-Por+ Pourable Denture Base to bench set prior to curing. Place into the pressure pot **IMMEDIATELY** after pouring.
  10. Do not use \*Cetylclide-G® in solution to soak models.

**Adverse Reactions:**

1. Corneal damage, headache, nausea, vomiting may occur with exposure to methyl methacrylate monomer. (See Warnings and Precautions statements.)
2. Allergic contact dermatitis and other allergic reactions may occur in susceptible individuals. Residual monomer in fully cured materials can be minimized by soaking the cured prosthesis in warm water for several days.
3. Particulates will be generated when grinding acrylate resins. Eye, skin and respiratory irritation may occur if appropriate engineering controls are not used.

**STEP-BY-STEP INSTRUCTIONS: PROPERTIES AND RECOMMENDED PROCEDURES** - The following table describes typical properties of Lucitone Fas-Por+ Pourable Denture Base, plus recommended temperatures and other important processing information. Further detail of this information is found later in the text.

Type and Class (according to ADA Spec. 12)	Type II, Class 1
Storage temperature for powder and liquid	60°F-80°F (16°C-27°C)
Powder / Liquid ratio (Use Lucitone Fas-Por+ Liquid)	20 g (29 cc) / 14ml
Mixing time (time required to wet all particles)	15 seconds
Time to reach pourable stage @ 70°F ± 2°F	pour immediately
Pour Time	1 minute
Material used to prepare mold	hydrocolloid
Temperature of hydrocolloid when pouring	60°F-80°F (16°C-27°C)
Cure time, temperature, and pressure	pressure pot, 15 minutes @ 20-30psi, 106°F ± 2°F (41°C ± 1°C)

**FABRICATION OF FULL AND PARTIAL DENTURES**

**HYDROCOLLOID INVESTMENT AND SPRUING** - Incorporate a slight taper into base of cast to facilitate removal from the hydrocolloid. Taper the cast so the base is the widest. **Soak the waxed stone cast in slurry water (saturated calcium sulfate dihydrate solution) for 10-20 minutes prior to pouring the hydrocolloid mold.** Do not use \*Cetylclide-G® in the solution. Follow the manufacturers instructions for preparation and use of the selected hydrocolloid product and flask. Bore sprue holes in the hydrocolloid investment as shown in Figure 1. **IMPORTANT:** When fabricating partials each acrylic area must have two sprues to avoid trapped air.

**WAX ELIMINATION AND SEPARATOR APPLICATION** - Boil off wax with clean boiling water. **DO NOT** use detergent to clean casts. Detergent residue may impact color uniformity. Also clean teeth with clean boiling water. Add mechanical retention to the teeth. Remove only enough tooth material to enhance retention. Excessive removal of tooth material will weaken the tooth and may lead to fracture.

Immediately after boil out, place the cast in a bowl of warm water for 5 minutes. Remove the cast from the water and allow to dry approximately 1-2 minutes.

Apply a thin coat of Al-Cote® Separating Agent to the cast and allow to dry until a slight sheen remains. Carefully reposition the teeth and cast back into the hydrocolloid mold. Reassemble the base and cover of the flask, being especially careful to avoid possible tooth and / or cast movement.

**MIXING** - To assure uniformity, shake powder jar (with a rocking motion) to prevent concentration of small particles at the bottom of jar. For the fabrication of 1 full denture **use 29 cc's (20 grams) of powder to 14 ml Lucitone Fas-Por+ Liquid. Do not use DENTSPLY Repair Material/ Lucitone Fas-Por E pourable denture base liquid.** Stir sufficiently (15 seconds) to assure wetting of all powder particles. Do not over-mix. **Pour immediately!**

**POURING** - Lucitone Fas-Por+ Pourable Denture Base should be **poured immediately after mixing.** Pour into one of the sprue openings with the flask positioned at a slight angle.

**IMPORTANT:** Pouring into more than one sprue opening may cause trapped air.

**POUR TIME** - The pour time of Lucitone Fas-Por+ Pourable Denture Base is approximately 1 minute. If the pour is done quickly this time enables the fabrication of 2 full dentures with one mix of material.

**IMPORTANT:** Immediately after pouring the dentures they must be placed into a pressure pot. **DO NOT BENCH SET.**

**CURING** - Check the water bath temperature in the pressure pot with an accurate thermometer. Place the flask in the water at 106°F ± 2°F (41°C ± 1°C) at 20-30 psi for 15 minutes. The water level should be at sprue height. **DO NOT FULLY SUBMERGE THE FLASK IN WATER.** If the water temperature is not maintained at approximately 106°F (41°C), the cure time will have to be increased to achieve full properties. For example, cure time should be 30 minutes at 73°F (23°C).


Figure 1. Spruing Techniques

**FINISHING AND POLISHING** - Lucitone Fas-Por+ Pourable Denture Base should be finished and polished using normal procedures.

**REPAIRING DENTURES**

Lucitone Fas-Por+ Pourable Denture Base is recommended for repairing acrylic denture materials.

**PREPARATION OF REPAIR AREAS**

Prepare the area by grinding a bevel on the edges. Wet the repair area with the Lucitone Fas-Por+ Liquid to enhance bonding.

**PROCESSING**

When repairing acrylics Lucitone Fas-Por+ Pourable Denture Base should be cured in a pressure curing unit. The result is maximum density in the repaired area. Porosity and/or surface blanching may occur

if a pressure curing unit is not used. Submerge the repaired appliance (and cast as appropriate) in a warm water bath at 106°F (41°C ± 1°C) and 20-30 psi pressure for 15 minutes.

**INDIRECT RELINE PROCEDURES**

Lucitone Fas-Por+ Pourable Denture Base is also recommended for relining poorly fitting dentures by the indirect method. A relining jig is recommended. Grind the tissue surface to relieve the denture. Soak the cast in slurry water for 10-20 minutes to eliminate air from the stone (prevents air bubbles in the acrylic). Wet the tissue side of the denture with the Lucitone Fas-Por+ Pourable Denture Base Liquid, then reseat the cast / denture on the jig. Immediately mix the Powder / Liquid (29cc / 14ml) and allow to bench set to reach the proper viscosity for the relining procedure. Process same as for repairs and dentures. A pressure pot is required to achieve the highest possible density.

MATERIAL SAFETY DATA SHEET

FORM NO.: 186 SECTION I - GENERAL INFORMATION

1.1 Product Trade Name: Lucitone® Fas-Por™ Pourable Denture Base Powder  
 1.2 Part (Item) Number: 682224, 682234, 682244, 682254, 682520, 682620, 682720, 682820, 682111, 682211, 682311  
 1.3 Chemical Name: A blend of homo-, co-, and ter-polymers of methyl methacrylate.  
 1.4 Synonyms: N/A  
 1.5 Division Name: DENTSPLY Trubyte  
 Address: 570 W. College Ave.  
 City, State, Zip: York, PA 17405-0872  
 1.6 Phone No. for Info: 717-845-7511  
 Dentsply

\*FOR CHEMICAL EMERGENCY\*  
 SPILL, LEAK, FIRE, EXPOSURE OR ACCIDENT  
 CALL CHEMTREC - DAY OR NIGHT

1.7 Emergency Phone No.: 800-424-9300 1.12 National Stock No.: N/A  
 Chemtrec 1.13 Specification No.: N/A  
 1.8 Date: 6/24/97 1.14 Chemical Family: Acrylic  
 1.9 Dentsply Contact: Scott E. Shaffer  
 1.10 Signature: *Scott E. Shaffer* 1.15 Type/Grade/Class: N/A  
 1.11 Formula: N.A. 1.16 NRC/EPA Number: N/A

SECTION II - HAZARDOUS INGREDIENTS

2.1 CAS Number	Chemical Name (Ingredients)	%	TLV	OSHA PEL
None				

SECTION III - PHYSICAL PROPERTIES

3.1 Boiling Point: N.A.	3.11 Specific Gravity: N.A.
3.2 Vapor Pressure: N.A.	3.12 pH: N.A.
3.3 Vapor Density: N.A.	3.13 Evaporation Rate: N.A.
3.4 Critical Temp: N.A.	3.14 Viscosity: N.A.
3.5 Decomposition Temp: N.A.	3.15 % Volatile by Vol: N.A.
3.6 Melting/Freezing Pt.: N.A.	3.16 Magnetism: N.A.
3.7 Solubility in Water: N.A.	3.17 Autoignition Temp: N.A.
3.8 Critical Pressure: N.A.	3.18 Corrosion Rate: N.A.
3.9 Perm. Exp. Limit: N.A.	3.19 Molecular Weight: N.A.
3.10 Appearance and Odor: Pink free-flowing powder.	

SECTION IV - FIRE AND EXPLOSION HAZARD DATA

4.1 Flash Point (Method Used): Approx. 300°C (570°F)  
 4.2 Flammable (Explosive) Limits:  
 LEL: 20 g/m<sup>3</sup> (dust cloud) UEL: N.A.  
 4.3 Extinguishing Media: All standard fire fighting equipment is suitable, i.e. CO<sub>2</sub>, foam, water spray, and water fog.  
 4.4 Special Fire Fighting Procedures and Equipment: Avoid extinguishing methods which may generate dust clouds. Water stream can disperse dust in air, producing a fire hazard and possible explosion hazard if exposed to ignition source. Fireman should wear protective clothing and self-contained respirators if large amounts are involved.  
 4.5 Unusual Fire and Explosion Hazards: Polymer dust is combustible. The explosive limits of the polymer particles suspended in air are approximately those of coal dust. Fire fighters should wear self-contained breathing apparatus. Irritating gases may be emitted.

SECTION V - HEALTH HAZARD DATA

5.1 Effects of Acute and Chronic Over Exposure: It is not known to cause significant health problems. It is considered an inert or nuisance dust. Avoid inhalation of dust. Keep dust out of eyes to prevent irritation.  
 5.2 Emergency and First Aid Procedures: INHALATION: Treat as inert particulate matter, seek medical attention if inhalation of large quantities of dust cause respiratory distress.

	8hr TWA	10m STEL
Total inhalable dust	mg/m <sup>3</sup>	mg/m <sup>3</sup>
Respirable dust	10	20
	5	

Lucitone® Fas-Por™ Pourable Denture Base Powder, Cont'd

SKIN CONTACT: Harmless to most people, but individuals with a history of skin disease or allergy may have a reaction to prolonged exposure. Wash well with soap and water. Seek medical attention if a problem persists. EYE CONTACT: Flush with stream of water for at least 15 minutes. Seek medical attention in the case of prolonged discomfort. INGESTION: Non hazardous, in all reasonably foreseeable circumstances.

SECTION VI - REACTIVITY DATA

- 6.1 Unstable: Stable: X
- 6.2 Conditions to Avoid (Stability): Elevated temperatures.
- 6.3 Incompatibility (Materials to Avoid): Reducing and Oxidizing Materials.
- 6.4 Hazardous Decomposition Products: Acrid fumes and carbon monoxide and/or carbon dioxide may occur during a fire.
- 6.5 Hazardous Polymerization: May Occur: Will Not Occur: X
- 6.6 Conditions to Avoid (Polymerization): Heat at or above 300°C (570°F).

SECTION VII - SPECIAL PRECAUTIONS AND SPILL OR LEAK PROCEDURES

- 7.1 Handling and Storage Precautions: Normal practice in handling of powders should be followed with dust being contained within equipment wherever possible. Where containment is not possible dust ventilation should be provided; otherwise, operators are advised to wear dust masks to avoid inhaling fine particles. Store in clean, dry area without sources of ignition.
- 7.2 Other Precautions: Store away from food and beverages.
- 7.3 Steps to be Taken in Case Material is Released or Spilled: Individual particles are spherical and spills can be extremely slippery. All spills should be cleaned up immediately.
- 7.4 Waste Disposal Procedures: May be disposed of in a landfill or incinerated. Follow Federal, State, and Local Regulations for disposal.
- 7.5 Neutralizing Agent: N.A.

SECTION VIII - OCCUPATIONAL PROTECTION MEASURES

- 8.1 Respiratory Protection: Considered a nuisance dust. Use a dust mask while measuring out the powder and grinding.
- 8.2 Ventilation: Use local exhaust when grinding cured material.
- 8.3 Protective Gloves: For users with skin problems, cuts or abrasions - wear protective gloves.
- 8.4 Eye Protection: Use safety glasses while grinding.
- 8.5 Other Personal Protective Equipment: None required.

SECTION IX - TRANSPORTATION

- 9.1 Applicable Regulations: IMOC: N.A. IATA: N.A. 49 CFR: N.A. Military Air (AFR 71-4): N.A. Tariff 6D: N.A.
- 9.2 Shipping Name: N.A.
- 9.3 I.D. Number: N.A.
- 9.4 Report Quantity: N.A.
- 9.5 Hazard Class: N.A.
- 9.6 Labels: N.A.
- 9.7 Aerosol Propellant(s): N.A.
- 9.8 Net Explosive Weight: N.A.
- 9.9 Unit Container: N.A.
- 9.10 DOT SPSC Container: N.A.
- 9.11 DOT Exempt/DOD/CCN: N.A.
- 9.12 Limited Quantity: N.A.

MATERIAL SAFETY DATA SHEET

FORM NO.: 185 SECTION I - GENERAL INFORMATION

- 1.1 Product Trade Name: Lucitone® Fas-Por™ Pourable Denture Base Liquid
- 1.2 Part (Item) Number: 682007, 682049, 682415
- 1.3 Chemical Name: Methyl methacrylate monomer
- 1.4 Synonyms: N/A
- 1.5 Division Name: DENTSPLY Trubyte  
Address: 570 W. College Ave.  
City, State, Zip: York, PA 17405-0872
- 1.6 Phone No. for Info: 717-845-7511  
Dentsply

"FOR CHEMICAL EMERGENCY"

SPILL, LEAK, FIRE, EXPOSURE OR ACCIDENT  
CALL CHEMTREC - DAY OR NIGHT

- 1.7 Emergency Phone No.: 800-424-9300 1.12 National Stock No.: N.A.  
Chemtrec 1.13 Specification No.: N.A.
- 1.8 Date: 8/13/02 1.14 Chemical Family: N.A.
- 1.9 Dentsply Contact: Scott E. Shaffer Methacrylic Esters
- 1.10 Signature: *Scott E. Shaffer* 1.15 Type/Grade/Class: N.A.  
1.16 NRC/EPA Number: N.A.
- 1.11 Formula: N.A.

SECTION II - HAZARDOUS INGREDIENTS

CAS Number	Chemical Name (Ingredients)	%	TLV	OSHA PEL
80-62-6	Methyl Methacrylate		93-97	100ppm
2082-81-7	1,4-butanedioldimethacrylate		2.5-6.5	

SECTION III - PHYSICAL PROPERTIES

- 3.1 Boiling Point: 101°C(214°F) 3.11 Specific Gravity: 0.94(820°C/68°F)
- 3.2 Vapor Pressure: N.A. 3.12 pH: N.A.
- 3.3 Vapor Density: N.A. 3.13 Evaporation Rate: >1(BAC=1)
- 3.4 Critical Temp: N.A. 3.14 Viscosity: N.A.
- 3.5 Decomposition Temp: N.A. 3.15 % Volatile by Vol: approx 97%
- 3.6 Melting/Freezing Pt.: -48°C(-54°F) 3.16 Magnetism: N.A.
- 3.7 Solubility in Water: 1.6g/L(820°C) 3.17 Autoignition Temp: 430°C(806°F)
- 3.8 Critical Pressure: N.A. 3.18 Corrosion Rate: N.A.
- 3.9 Perm. Exp. Limit: N.A. 3.19 Molecular Weight: N.A.
- 3.10 Appearance and Odor: Colorless liquid with an ester-like odor.

SECTION IV - FIRE AND EXPLOSION HAZARD DATA

- 4.1 Flash Point(Method Used): 10°C (50°F)
- 4.2 4.2 Flammable (Explosive) Limits:  
LEL: 2.12% UEL: 12.5%
- 4.3 Extinguishing Media: CO<sub>2</sub>, alcohol-resistant foam, dry chemical, water spray, and water fog.
- 4.4 Special Fire Fighting Procedures and Equipment: Water may be ineffective unless used as a fine spray/cloud. Use water spray to cool containers of this liquid exposed to fire. Fireman should wear protective clothing and self-contained respirators if large amounts of liquid are involved.
- 4.5 Unusual Fire and Explosion Hazards: Vapors may travel to source of ignition and flash back. Avoid ignition sources! Heat can induce polymerization with rapid release of energy and closed containers may rupture explosively. Keep containers cool with water spray. Formation of explosive liquid/air/steam mixtures is possible.

SECTION V - HEALTH HAZARD DATA

- 5.1 Effects of Acute and Chronic Over Exposure: Eye Contact: irritation and possible corneal damage. Inhalation: headache, nausea, drowsiness, and unconsciousness.

## Lucitone® Fas-Por™ Pourable Denture Base Liquid, cont'd

- 5.2 **Emergency and First Aid Procedures:** **INHALATION:** Irritating to respiratory system if present in high vapor concentration; can induce nausea and/or headache. Avoid breathing high concentrations of vapor. If overexposed, move to fresh air and get medical attention if problem persists. **SKIN CONTACT:** Can cause contact dermatitis in some individuals. People with a history of skin disease or allergy may have a reaction to prolonged exposure. Wash well with soap and water to minimize problem. Seek medical attention if an irritation develops and persists. **EYE CONTACT:** Irritating to eyes. If liquid contacts eyes, hold eye open and immediately flush with a stream of water for at least 15 minutes and seek medical advice.
- INGESTION:** Avoid ingestion. If conscious give 2 glasses of water and consult a physician immediately.

## SECTION VI - REACTIVITY DATA

- 6.1 **Unstable:** Stable: X (when stored properly)
- 6.2 **Conditions to Avoid (Stability):** Elevated temperatures and direct sunlight.
- 6.3 **Incompatibility (Materials to Avoid):** Reducing and Oxidizing Materials, acids, bases.
- 6.4 **Hazardous Decomposition Products:** None known.
- 6.5 **Hazardous Polymerization:** May Occur: X Will Not Occur:
- 6.6 **Conditions to Avoid (Polymerization):** Heat, direct sunlight, oxygen-free atmosphere, inhibitor depletion (due to excessive aging), contamination with polymerization catalysts.

## SECTION VII - SPECIAL PRECAUTIONS AND SPILL OR LEAK PROCEDURES

- 7.1 **Handling and Storage Precautions:** Use adequate ventilation and keep containers tightly sealed when not in use. Store in a cool and dry place; store out of direct sunlight. Store away from food and beverages.
- 7.2 **Other Precautions:** Do not smoke while using liquid. Keep away from heat and ignition sources. Avoid breathing on a regular basis.
- 7.3 **Steps to be Taken in Case Material is Released or Spilled:** Eliminate ignition sources! Dike and absorb spill with inert material (sawdust, sand, or vermiculite) and use non-sparking tools to transfer to containers for disposal. Remove saturated clothing and wash affected skin with soap and water.
- 7.4 **Waste Disposal Procedures:** After addition of excess inhibitor, incinerate under controlled conditions in a safe area according to Federal, State and Local Regulations.
- 7.5 **Neutralizing Agent:** N.A.

## SECTION VIII - OCCUPATIONAL PROTECTION MEASURES

- 8.1 **Respiratory Protection:** In the case of exposure to low concentration of vapors, an Organic Vapor Respirator, such as, NIOSH approved TC-23C-859, provides adequate protection.
- 8.2 **Ventilation:** Use local exhaust to keep below TLV and LEL. Mechanical ventilation should have explosion-proof motors.
- 8.3 **Protective Gloves:** Wear impervious gloves (butyl rubber preferred).
- 8.4 **Eye Protection:** Use splash proof goggles.
- 8.5 **Other Personal Protective Equipment:** Wear impervious apron and overshoes during clean-up operation.

## SECTION IX - TRANSPORTATION

- 9.1 **Applicable Regulations:** IMOC: IATA: N.A. 49 CFR: X  
Military Air (AFR 71-4): Tariff 6D:
- 9.2 **Shipping Name:** Methyl Methacrylate Monomer Stabilized
- 9.3 **I.D. Number:** UN1247 9.9 Unit Container: N.A.
- 9.4 **Report Quantity:** N.A. 9.10 DOT SPSC Container: N.A.
- 9.5 **Hazard Class:** 3 9.11 DOT Exempt/POD/CCN: N.A.
- 9.6 **Labels:** No 9.12 Limited Quantity: Yes
- 9.7 **Aerosol Propellant(s):** None
- 9.8 **Net Explosive Weight:** N.A.

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\*CetylcodeG® is not trademark of DENTSPLY
© 1997, 1999 DENTSPLY International Inc. All rights reserved.
**DENTSPLY International Inc.**
**DequDent GmbH**
Printed in U.S.A. 582898 Rev. H (11/03)