DENTSPLY

ECLIPSE[®]
PROSTHETIC RESIN SYSTEM

**Laboratory Directions For Use
and Clinical Directions For Use**

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CE

CAUTION: U.S. Federal Law restricts this device to sale by, or on the order of, a licensed dental professional.

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Indications for Use: Full and partial removable dentures, occlusal splints, nightguards, and provisional partial dentures.

Contraindications for Use: See precautions and warnings.

Warnings: This product contains polymerizable materials. It does not contain methyl methacrylate (MMA) monomer. It is possible that some dental health care workers may develop an allergy or sensitivity to the product. Wash thoroughly with soap and water after contact. If dermatitis or other symptoms persist, seek medical advice. Do not immerse Eclipse appliances in boiling water.

Mechanical tooth retention must be used in accordance with instructions to prevent debonding.

May not be appropriate for partial dentures where limited vertical dimension is present, or complete dentures with less than 12mm interarch space (6mm for single denture).

Precautions:

- This product is intended to be used only as specifically outlined in the Directions For Use. Any use of this product inconsistent with the Directions For Use is at the discretion and sole responsibility of the practitioner.
- Insufficient clinical data exists at present for use with implant supported removable prostheses and removable appliances employing attachments.
- Users with special skin problems, cuts or abrasions may choose to wear protective gloves.
- When grinding prosthodontic resins, proper ventilation, masks and vacuum systems should be used.
- Since the Eclipse resins will become the final prosthesis, cleanliness is very important to an acceptable outcome. Wash hands before using uncured

Eclipse materials.

- Store at 60-80° F (16-27° C), away from moisture and in a light-safe container. The materials may polymerize prematurely if stored at excessively high temperature or if exposed to a visible or UV light source. Do not expose to direct sunlight.
- Unused baseplate and clear resin materials should be cured prior to disposal. Small amounts of set-up and contour resin may be disposed of without curing.
- Items removed from the curing unit are hot! Use heat-resistant gloves provided.
- Do not remove prostheses from models until cool.
- If the resin materials will not release easily from the package material (may occur in warm ambient conditions), place the package in a refrigerator/freezer for a few seconds to facilitate removal.
- The following disinfectants have been used successfully with cured and uncured Eclipse resins

- Glutaraldehyde-based disinfectants
- Quaternary ammonium chlorides
- Quaternary ammonium chlorides/isopropyl alcohol (low concentration)
- Spray-based phenolics

Use of other disinfectant products may not be suitable with Eclipse. Avoid disinfectants containing iodine.

Adverse Reactions:

- Allergic contact dermatitis and other allergic reactions may occur in susceptible individuals.
- Dust will be generated when grinding these materials. Eye, skin and respiratory irritation may occur if appropriate engineering controls are not used.

NOTE: Stop and read operating manuals for equipment before continuing to read Laboratory Directions for Use.

LABORATORY DIRECTIONS FOR USE:

I. FULL DENTURES

A. Baseplate

- Box and pour the master model. The land area of the model should be 1mm or less above the depth of the sulcus to allow light penetration for complete curing.
- Carve the post-dam in the master model if not included in the impression.
- Due to undercuts, a duplicate or working model may be required for articulation and set-up. The baseplate must be cured on the master model.
- Place temperature indicator on model. Paint a thin layer of separator (Al-Cote® Separating Agent) on the dry master model and allow separator to dry completely.

Note: Since the baseplate and set-up will become the final prosthesis, cleanliness is very important to an acceptable outcome.

- Heat the master model to 120-130°F (49-54°C) in the DENTSPLY® Conditioning Oven set at 131°F (55°C). Heating time will vary depending on the number of models in the DENTSPLY Conditioning Oven and the water content of the gypsum. Temperature Indicator (provided) will turn black when model is at correct temperature.
- Place the **rounded side** of the upper or lower size Baseplate Resin arch on the ridge and allow the Baseplate Resin to warm on the model for 30 seconds. Carefully adapt the labial/buccal surfaces to avoid air entrapment. Begin adaptation at


the crest of the ridge and slowly and carefully adapt the material down the buccal slope and into the vestibule last. Follow the same procedure for the palatal/lingual section. For the maxillary baseplate, carefully join the two sides at the midline without trapping air.

- Apply Eclipse Air Barrier Coating to the entire surface of the baseplate, and place it in the Eclipse Processing Unit within 10 minutes; then start the baseplate cure sequence. (Menu #1)
- When the curing cycle is complete, remove the baseplate/model from the Eclipse Processing Unit and allow the baseplate/model to cool on the bench until reaching ambient temperature.
- Soak the baseplate and master model in tap water for 5 to 15 minutes to facilitate baseplate removal from the master model. If undercuts are present, section the master model carefully with a separating disc from the underside prior to removal.
- Trim and smooth the baseplate borders. If the master model is damaged during removal of the baseplate, the duplicate model may require modification to allow seating of this baseplate for articulation and tooth set-up. Alternatively, blocking out undercuts in the processed baseplate can help make a working model. When making a working model, protect the borders by applying a utility wax bead prior to pouring stone or plaster into the processed baseplate.

B. Occlusion Rims

- Place wax occlusal rims on the baseplate and send to the clinician for aesthetic and occlusal records.
- When returned, mount the models using the baseplate(s) and occlusal rims with the occlusal

- registration. Transfer occlusal markings and aesthetic markings (i.e. midline, canine, smile line, dominant side) to the opposing cast.
- While the baseplate is still **on the model**, completely clean all traces of wax from the baseplate by rinsing with clean boiling water, or use a lab steamer. **DO NOT IMMERSE IN BOILING WATER. Complete removal of wax is critical to success.**
 - Allow the baseplate to cool to room temperature while still on the model.
 - Roughen the entire outer surface of the baseplate with a coarse carbide bur or diamond to enhance bonding. 


A butt joint is recommended about 1-2 mm above the depth of the sulcus. The goal is to eliminate feathered edges of Contour Resin on the flange.
 - Clean the entire prepared surface, thoroughly scrubbing with a clean denture brush and tap water; then dry. Following this, use a brush and isopropyl or ethyl/denatured alcohol to clean the surface; then blow dry with clean compressed air.

C. Tooth Set-Up


Interim Preparation Steps:

- Soften Set-up Resin** by placing the package inside the DENTSPLY Conditioning Oven at 131°F (55°C) for about 5-10 minutes if desired.
- Melt Contour Resin** in the DENTSPLY Melting Pot on setting 4, or at a temperature of 189°F (87°C). Note: Clean the Contour Resin from the DENTSPLY Melting Pot after 8 hours to assure that only fresh material is used. Limit light exposure to the melted Contour Resin by keeping the lid on the DENTSPLY Melting Pot closed when not transferring resin

3

- to the denture.
- Remove wax from teeth with clean boiling water using the DENTSPLY Tooth Cage, or other means.**
 - Store prepared and cleaned teeth in the DENTSPLY Tooth Box until setting teeth.
 - Deglaze teeth.
 - When using plastic teeth, cut mechanical retention as follows:
 - For **anterior teeth**, use **keyed retention slots and cut a collar groove around the tooth or in the lingual/gingival area.**
 - For **posterior teeth**, use **keyed retention slots or cut a collar groove around the tooth if insufficient vertical dimension exists to cut the keyed slots.**
 - Keyed Retention Slots:** Use a handpiece and the **DENTSPLY Knife Edged Diamond (re-order #905128)** to cut two keyed retention slots into each tooth as shown. It is important that the two slots are **not** made parallel to each other. Maximize the distance between the slots and cut at a 45° angle. The slot depth should be 1 to 1.5 mm deep for small to large teeth respectively. Similarly the slot width should be 0.75 to 1.5 mm for small to large teeth respectively.
- 

- Collar Groove:** Use a hand-piece and the **DENTSPLY Collar Groove Diamond (reorder #905129)** to grind a groove into the collar area of the tooth as shown here. The collar retention should be a groove and not a slot. This should be no sharp angles in the groove. The groove should be positioned 0.5 to 1.0 mm from the edge of the ridge lap for small to large teeth respectively. Similarly the groove depth should be 0.5 to 1.0 mm for small to large teeth respectively.
 - Clean the ridge lap area by lightly scrubbing with a clean, dry stiff brush; then blow with clean compressed air.
- 
- When using porcelain teeth with pins and/or vent holes, no tooth preparation is required.
 - The Eclipse Set-up Resin may be used as supplied or cut lengthwise to preferred diameter. Melt the bottom of the Set-up Resin with the flat surface of the electric spatula tip; then place it onto the ridge area of the baseplate. Use as little Set-up Resin as possible.
 - Fill the tooth retentive slots with Set-up Resin.** Recommended method is to flow the Set-up Resin into the tooth retentive slots, using a small tip on the **DENTSPLY Electric Spatula. Do not entrap air.**
 - Melt the surface of the Set-up Resin with the DENTSPLY Electric Spatula or DENTSPLY Hot Air Gun just prior to tooth placement. Set teeth while Set-up Resin is in a softened condition. Use the

4

DENTSPLY Electric Spatula set at 230-240°F (110-116°C) to flow molten Contour Resin around teeth to hold in position.

NOTE: Do not "slice" Set-up Resin with DENTSPLY Electric Spatula as you would wax, as this will entrap air.

- Complete contouring of the denture base by using the DENTSPLY Electric Spatula set at 230-240°F (110-116°C) to flow liquified Contour Resin onto the denture as needed. Observe the Contour Resin as it is applied onto the denture and "pop" any bubbles that are seen while the resin is still fluid. Contour Resin should cover all Set-up Resin and exposed baseplate. Contour to desired shape.

NOTE: Allow Contour Resin to cool before contouring.

NOTE: Do not whip the Contour Resin to re-suspend red fibers as they settle out. Stir gently to resuspend fibers, which will minimize air entrapment.

NOTE: Contour the denture to the finished dimensions. Do not "bulk-up."

- Carve around teeth with a sharp instrument and festoon Contour Resin with the instrument of your choice. Make sure teeth are clean of all flash/Contour Resin for best results. Instruments should be free of wax and dedicated for Eclipse materials only. Use the **DENTSPLY Hot Air Gun to smooth the Contour Resin around teeth and gingival area, as well as the junction of the baseplate and Contour Resin.** At this point, the denture is ready for try-in.

NOTE: Due to light curing nature of the material, the uncured try-in must be sent to the dentist in a Light Safe Bag (provided in Baseplate Resin package).

D. Adjustments

1. If teeth must be re-positioned in the lab after the try-in visit, this is best accomplished by keeping the uncured denture cool and breaking the teeth away from the uncured resin. Refrigerating for a short time will make this easier. The teeth can then be re-set with a melt interface using a small amount of Set-up and Contour Resin. Melting the resins to loosen and re-position the teeth may cause air entrapment and handling challenges.
2. Check uncured denture for tight tooth-to-base margins. Re-adapt if needed.
3. Remove model from articulator mounting.

E. Processing

There are two methods for preparing a denture for processing to ensure the best fit:

- A. Pour a fresh model into the denture when ridges have been broken, adjusted, or if patient has flat ridges (see step 1).
 - B. Use existing model and Eclipse Gel to seal borders on other case types (go to step 2).
1. If pouring an **upper** model (100% stone or 50:50 mix) into the denture, keep gypsum away from labial and lingual surfaces to avoid shadows that may result in areas of undercure.
 2. If pouring a **lower** model (100% stone or 50:50 mix) into the denture, keep gypsum away from labial surfaces and keep the lingual open (model forms a horseshoe) to avoid shadows that may result in areas of undercure. 
 3. Apply Eclipse Air Barrier Coating to denture, covering all external resin surfaces and teeth.
 4. Place the denture/model in the DENTSPLY Conditioning Oven set

5

at 130°F (55°C) and pre-heat for at least 1 hour (but not longer than 8 hours).

5. Remove denture/model from the DENTSPLY Conditioning Oven. If a fresh model was used, go to step 5 immediately. If using existing model, follow these steps:
 - Use a brush to apply Eclipse Model Release Agent to denture flange areas. Do not allow Eclipse Model Release Agent to come in contact with the land area of the model.
 - Seal denture to model by using the Eclipse Gel on the borders. Seal all border areas except within about 1/4 inch of the buccal frenum attachments. Do this even if frenum areas are not very pronounced or are non-existent.
6. Place the denture/model in the Eclipse Processing Unit on the center of the turntable; then start the curing process. (Menu #2)
7. When the cure cycle is complete, use heat-resistant gloves (provided) to remove the denture/model from the Eclipse Processing Unit. Then allow the denture/model to cool until the denture reaches ambient temperature. Then separate the denture from the model. **DO NOT TOUCH WHILE HOT! ALLOW TO COOL TO ROOM TEMPERATURE BEFORE REMOVING CURED GEL! (if used).**
8. Remove the cured Eclipse Gel (if used) and wash away the Eclipse Air Barrier Coating with water and a soft brush.
9. Remount the denture(s) and verify occlusion.
10. Finish denture using a lathe or handpiece and polish in the usual manner.

II. Removable Partial Dentures

Framework Design: Use large gauge retentive mesh for saddle areas. This will improve the ability to adapt the material and minimize the chance of cure shadows or air entrapment. Solid metal areas should be minimized and used only where strength of the framework is a concern. Rounded mesh designs are also better than the square patterns. If possible, use a 1mm (26 gauge) blockout thickness.

A. Baseplate

1. Box and pour the master model. The land area of the model should be 1 mm or less above the depth of the sulcus to allow light penetration for complete curing.
2. Due to undercuts, a duplicate or working model may be required for articulation and set-up. The baseplate must be cured on the master model.
3. Paint a thin layer of Al-Cote Separating Agent on the dry master model and allow separator to dry completely.

NOTE: Since the baseplate and set-up will become the final prosthesis, cleanliness is very important to an acceptable outcome.

4. Heat the master model to 120-130°F (49-54°C) in the DENTSPLY Conditioning Oven set at 131°F (55°C). Heating time will vary depending on the number of models in the DENTSPLY Conditioning Oven and the water content of the gypsum. Temperature Indicators (provided) placed on model will turn black when model is at the proper temperature.
5. Cut an appropriate section from an Eclipse Baseplate Resin arch. Apply Baseplate Resin with the rounded side to the crest of the ridge of the warm model and work material to the desired extension. Avoid air entrapment.
6. Warm the cast partial framework

6

in the DENTSPLY Conditioning Oven or with the DENTSPLY Hot Air Gun; then seat the warmed framework onto the model and into the baseplate material. Make sure the rests and tissue stops are fully seated.

7. Adapt the baseplate material around retentive frame areas, using finger pressure.
8. Trim excess Baseplate Resin from around the major connector with the DENTSPLY Electric Spatula set to 230-240°F (110-116°C) and add resin as needed to complete the saddle area.
9. Apply Eclipse Air Barrier Coating to the entire surface of the baseplate. Within 10 minutes, place it in the Eclipse Processing Unit; then start the processing sequence for partial baseplate. (Menu #3)
10. Allow baseplate/model to cool to ambient temperature.
11. Soak model/partial in tap water for 5 to 15 minutes to facilitate removal. If necessary, section the model carefully with a separating disc from the underside prior to removal.
12. Trim and smooth the baseplate borders.

B. Occlusal Rims

1. If bite registration is required, return to the clinician for occlusal records.
2. When returned, mount the models using the bite registration. Transfer any necessary occlusal markings. Completely clean all traces of bite registration material from the baseplate (while on the model) by rinsing with clean boiling water or use a lab steamer. **DO NOT IMMERSE IN BOILING WATER. Complete removal of wax is critical to success.**
3. Allow the baseplate to cool to room temperature while still on the model.

4. Roughen the entire outer surface of the baseplate with a coarse carbide bur or diamond to enhance bonding. A butt joint is recommended about 1-2mm above the depth of the sulcus. The goal is to eliminate feathered edges of Contour Resin on the flange.
5. Clean the entire prepared surface, thoroughly scrubbing with a clean denture brush and tap water; then dry. Following this, use a brush and isopropyl or ethyl/denatured alcohol to clean the surface; then blow dry with clean compressed air.

C. Tooth Set-up

SEE FULL DENTURE SECTION - C. TOOTH SET-UP.

D. Adjustments and Processing

SEE FULL DENTURE SECTION, D. Adjustments and E. Processing, with one exception: Sealing Gel is not required for partial dentures. Place the partial denture in the Eclipse Processing Unit and cure. (Menu #2)

III. OPTIONAL METHOD FOR FULL AND PARTIAL DENTURES

Conventional dental wax can be used to set teeth in combination with an Eclipse processed baseplate to make a try-in denture.

1. After receipt of the try-in denture from the dental professional, use a matrix to establish tooth position in relation to the model. Then remove the wax from the processed baseplate and teeth.
2. **Create mechanical retention in the teeth as described in Full Denture section. Carefully fill retentive feature in teeth with Set-up Resin.**
3. Melt the bottom of the Set-Up Resin with the flat surface of the DENTSPLY Electric Spatula Tip; then apply Set-up Resin to the

roughened processed baseplate. Heat with the DENTSPLY Hot Air Gun to soften the Set-up Resin. Reposition the teeth/matrix onto the processed baseplate. Using DENTSPLY Electric Spatula, melt Set-up Resin around each tooth to ensure there are no air voids. Finish denture using Contour Resin as desired.

IV. Denture / Partial Denture Relines

A. Partial Dentures

1. Box and pour reline impression in dental stone.
2. Remove boxing wax and trim cast on a model trimmer to desired dimensions and key or index model. The land area of the model should be 1 mm or less above the depth of the sulcus to allow light penetration for complete curing.

NOTE: For best results and to eliminate the need to box and pour, a "pick-up" impression should be taken over the seated partial denture and natural teeth. Apply Eclipse Model Release Agent to the underside of the major connector to prevent the stone from locking the frame to the cast.

3. Sparingly lubricate the denture teeth with Eclipse Model Release Agent, and also lubricate the top of the cast with Eclipse Model Release Agent.
4. Using a standard reline jig, articulator, or duplicator, attach the model to the upper member with mounting plaster and allow the plaster to harden.
5. Make a thick mix of mounting plaster, attach it to the lower member of the articulator, jig, or duplicator, and close the mounted denture into the plaster mix to a 1/3 depth of the denture teeth.
6. Carefully remove the denture/model

from the mounting. Carefully remove the denture from the model and heat the model in the DENTSPLY Conditioning Oven set at 131°F (55°C) for at least 30 minutes to drive off moisture.

7. Completely remove the impression material and any adhesive from the denture base and grind off 1/3 of the flange area from the resin saddle. If the cast mesh retention element of the framework does not allow 1/3 of the resin to be removed, grind the flange area just to the edge of the retention mesh. **Take care not to expose the mesh.** Remove enough material from the ridge area to allow ample Eclipse Baseplate Resin to flow under the saddle area. Grind the edges of the flange to a beveled joint configuration.
8. Clean the entire prepared surface, thoroughly scrubbing with a clean denture brush and tap water; then dry. Following this, use a brush and isopropyl or ethyl/denatured alcohol to clean the surface; then blow dry with clean compressed air.
9. Remove the cast from the DENTSPLY Conditioning Oven and place a temperature indicator on the model. Apply a thin coat of Al-Cote Separating Agent to the model and allow to dry. Place the mounted model and the denture in its index in the DENTSPLY Conditioning Oven until the temperature indicator turns black.
10. Remove sufficient amount of Eclipse Baseplate Material from the package to complete the reline. Remove the cast from the DENTSPLY Conditioning Oven, adapt the Eclipse Baseplate Resin onto the tissue areas, and place the cast back in the DENTSPLY Conditioning Oven for 1 minute to soften the resin.
11. Remove the cast and denture from the DENTSPLY Conditioning Oven and further warm the saddle area of the partial denture with the DENTSPLY Hot Air Gun. Place the denture onto the cast, slowly compressing the Eclipse Baseplate Resin until the rests fully seat into their stops. Verify the occlusion by replacing the cast/denture assembly in the jig.
12. Remove the cast from the jig and, using the DENTSPLY Electric Spatula, smooth the junction between the cured and uncured resins. **Take care not to entrap air.** Warm the contoured areas slightly with the DENTSPLY Hot Air Gun.
13. Apply Eclipse Air Barrier Coating to the denture saddle areas and place into the Eclipse Processing Unit and cure. (Menu #4)
14. Allow to cool to ambient temperature and then carefully remove the partially cured reline from the model. Coat the tissue bearing surface with Eclipse Air Barrier Coating and insert the denture into the Eclipse Processing Unit, tissue side up, and cure. (Menu #6)
15. After curing, allow the denture to cool to ambient temperature and wash off the Eclipse Air Barrier Coating with a brush and water.
16. Finish and polish in the usual manner.

B. Full Denture

1. Follow Steps 1 through 6 as outlined in the Partial Denture Reline technique.
2. Trim 1/3 of the remaining flange to a beveled joint configuration. Remove palate on full upper to 1/4" of linguals of teeth.
3. Clean the entire prepared surface, thoroughly scrubbing with a clean denture brush and tap water; then dry. Following this, use a brush and isopropyl or ethyl/denatured alcohol to clean

- the surface; then blow dry with clean compressed air.
4. Remove the cast from the DENTSPLY Conditioning Oven and place a temperature indicator on the model. Apply a thin coat of Al-Cote Separating Agent to the model and allow to dry. Place the mounted model and the denture in its index in the DENTSPLY Conditioning Oven until the temperature indicator turns black.
 5. Apply the Eclipse Baseplate Resin arch to the warmed model and let set for 30 seconds to allow the heat from the model to transfer to the resin.
 6. Re-attach the upper and lower halves of the jig; then slowly compress the Baseplate Resin until the jig reaches its final closure.
 7. Separate the denture from its index; then use finger pressure to conform the resin to the palate and the beveled joint section of the denture. Make sure that no voids are formed during this process.
 8. Use the DENTSPLY Electric Spatula to smooth the interface between the uncured resin and the beveled joint section.
 9. Replace the denture in the jig to verify that the indexed vertical dimension is maintained.
 10. Separate the denture and matrix from the jig.
 11. Apply Eclipse Air Barrier Coating to the exposed surfaces of the denture.
 12. Place the denture/model into the Eclipse Processing Unit and cure. (Menu #4)
 13. Allow to cool to ambient temperature, then remove the partially cured relines from the model and the index from the jig. Coat the tissue bearing surface with Eclipse Air Barrier Coating. Insert denture into the prewarmed matrix and place the

- denture/matrix into the Eclipse Processing Unit, tissue side up and cure. (Menu #6)
14. After curing, allow the denture to cool to ambient temperature and separate from the index.
 15. Wash off the Eclipse Air Barrier Coating with water.
 16. Finish and polish in the usual manner.

V. PROVISIONAL PARTIAL DENTURES

1. **Prepare the teeth as indicated in I. C. Steps 3-6.**
2. **Fill retention areas in teeth with Eclipse Baseplate Resin and warm with DENTSPLY Hot Air Gun to "pop" any bubbles.**
3. Adapt the Eclipse Baseplate Resin to the warm model. Position tooth.
4. Process using basic procedures described in I. FULL DENTURES section E. (Menu #2)

VI. REPAIRS

A. Minor Voids in Full Denture Baseplate

1. For small voids, prepare the area with a coarse bur or diamond to enhance bonding. Do not leave feathered edges.
2. Apply Al-Cote Separating Agent to area of model where repair will be made and allow to dry.
3. Thoroughly clean prepared surface as described below:
 - Wash prepared area with a brush and tap water and blow dry with clean compressed air.
 - Use brush and isopropyl or ethyl/denatured alcohol to clean the prepared area and allow to dry.
 - Make sure no dust remains on the surface.
4. Locally warm the prepared area

- with the DENTSPLY Hot Air Gun.
5. Use the DENTSPLY Electric Spatula (230-240°F) to flow the Eclipse Baseplate Resin onto the prepared area.
 - Carefully flow from one side to the other to avoid air entrapment.
 - It is vital that a melt interface occurs between the Baseplate Resin and the prepared surface.
 - Fill to contour or slightly below. Do not overfill; otherwise fit may be affected!
 - If multiple voids need to be repaired, locally warm using the DENTSPLY Hot Air Gun to ensure that a melt interface occurs.
 6. Place baseplate on the model.
 7. Apply Eclipse Air Barrier Coating to any exposed surfaces.
 8. Process in the Eclipse Processing Unit. (Menu #5)
 9. Allow the repaired baseplate to cool before removing from model.
 10. Apply Eclipse Air Barrier Coating to the tissue side of the repaired area.
 11. Place the baseplate on the Eclipse Processing Unit turntable. Center the baseplate, not the repair area, on the turntable tissue side up and process in the Eclipse Processing Unit. (Menu #6)
 12. Immediately remove from the Eclipse Processing Unit and place the repaired baseplate back on the model and allow to cool to ambient temperature.
 13. Remove Eclipse Air Barrier Coating with water and brush.
 14. Finish the repair and polish in the usual manner.
2. Place the baseplate on the Eclipse Processing Unit turntable. Center the baseplate, not the repair area, on the turntable tissue side up, and cure. (Menu #6)
 3. Immediately remove the appliance from the Eclipse Processing Unit and place the repaired baseplate back on the model and allow to cool to ambient temperature.
 4. Remove Eclipse Air Barrier Coating with water and brush.
 5. Finish the repair and polish in the usual manner.

C. Denture Base Repair

1. Pour a model into the denture to support the denture during curing.
 - For major repairs, such as fractures, use sticky wax to hold the pieces together, prepare a matrix and then pour a model into the denture.
2. Remove the denture/model from the mounting. Remove denture from the model, apply Al-Cote Separator to the model and heat the model in the DENTSPLY Conditioning Oven set at 131°F (55°C). Temperature Indicators (provided) placed on models will turn black when model is at the proper temperature.
 
3. Prepare the denture by grinding a beveled surface. Do not leave feathered edges.
4. Thoroughly clean surface as described below:
 - Wash prepared area with a brush and tap water and blow dry with clean compressed air.
 - Use a brush to clean the prepared surface with isopropyl or ethyl/denatured alcohol and allow to dry.
 - Make sure no dust remains on the surface.

B. Minor Voids in Partial Denture Baseplate

1. Follow steps 1 through 7 as outlined in Minor Voids in Full Denture Baseplate.

5. Place on model and seat properly.
6. Locally warm the prepared area with the DENTSPLY Hot Air Gun.
7. Coat the prepared area with Eclipse Baseplate Resin by using the DENTSPLY Electric Spatula.
8. Place a small amount of Eclipse Baseplate Resin into the repair site.
9. Heat with the DENTSPLY Hot Air Gun until Eclipse Baseplate Resin begins to slump. Do not overheat (i.e. melt resin).
10. Fill the remaining repair area with Baseplate Resin and warm with the DENTSPLY Hot Air Gun. 11. Apply Eclipse Air Barrier Coating on repair area.
12. Apply Eclipse Model Release Agent to the denture flange areas, but do not allow Eclipse Model Release Agent to come in contact with the land area of the model.
13. Seal the denture to the model by using Eclipse Gel on the borders. Leave a small gap in the Eclipse Gel at the frenum notches, so that the border is not 100% sealed.
14. Place the baseplate on the Eclipse Processing Unit turntable. Center the baseplate, not the repair area, on the turntable and process on the model in the Eclipse Processing Unit. (Menu #4)
15. Allow the repaired denture to cool before removing from model.
16. Remove denture from model. If necessary, soak model and denture in tap water for five minutes to facilitate removal.
17. Apply Eclipse Air Barrier Coating to the repair area on the tissue side.
18. Place the repaired denture tissue side up on the Eclipse Processing Unit turntable and cure. (Menu #6)
19. Remove denture from the Eclipse

Processing Unit and immediately place back on model and allow to cool to ambient temperature.

20. After cooling, remove cured Eclipse Gel material and wash off Eclipse Air Barrier Coating with water and a brush.
21. Finish and polish in the usual manner.

D. Tooth / Surface Repair

1. Pour a model into the denture to support the denture during curing.
2. Roughen and prepare surface of denture.
3. **Prepare tooth retention as described in I. FULL DENTURES, Section C, Steps 3-6.**
4. Thoroughly clean denture and tooth surface as described below:
 - Wash prepared area with a brush and tap water and blow dry with clean compressed air.
 - Use a brush to clean the prepared surface with isopropyl or ethyl/denatured alcohol and allow to dry.
 - Clean the prepared tooth by lightly scrubbing with a clean dry stiff brush; then blow clean with compressed air.
5. Prepare a matrix using silicone putty or gypsum after positioning the tooth with sticky wax. Remove any wax with clean boiling water.
6. Fill the retention areas in the teeth with Baseplate Resin using the DENTSPLY Electric Spatula.
7. Locally warm the prepared area of the denture with the Eclipse Hot Air Gun.
8. Wet the prepared area with Baseplate Resin by using the DENTSPLY Electric Spatula.
9. Place a small amount of Baseplate Resin into the repair site.

10. Heat with the DENTSPLY Hot Air Gun until Baseplate Resin begins to slump. Do not overheat and melt resin.
11. Place tooth into the repair site.
12. Use Eclipse Baseplate Resin to fill the area around the tooth and slightly over contour.
13. Warm the Eclipse Baseplate Resin again with the DENTSPLY Hot Air Gun until the resin begins to slump.
14. Apply Eclipse Air Barrier Coating on repair area.
15. Follow VI. C. Denture Base Repair Steps 14-16.
16. After cooling, wash off Eclipse Air Barrier Coating with water and a brush.
17. Finish and polish in the usual manner.

VII. DENTURE IDENTIFICATION

On a fully cured denture:

1. Pour a model in the denture to support the denture during curing.
 2. Prepare an onion skin paper label (or Kimwipe®) with the patient's name and/or SSN.
- Heat-shrink labels do not work well due to continued shrinking during the cure cycle.**
3. Prepare the area by grinding a space in the denture base using a coarse bur or diamond.
 4. Clean prepared area by scrubbing with a brush and water; then dry. Use a brush to clean the prepared surface with isopropyl or ethyl/denatured alcohol and allow to dry.
 5. Locally warm area with the DENTSPLY Hot Air Gun.
 6. Use the DENTSPLY Electric Spatula to fill the bottom space with a small amount of melted Eclipse Clear Baseplate Resin.

7. Use the DENTSPLY Hot Air Gun to warm the resin.
8. Place the label and fill to contour of denture with a small amount of melted Clear Baseplate Resin.
9. Smooth surface with the DENTSPLY Hot Air Gun to blend with the surrounding area.
10. Apply Eclipse Air Barrier Coating on resin.
11. Follow VI. C. DENTURE BASE REPAIR Steps 14-16. 12. After cooling, wash off Eclipse Air Barrier Coating with water and a brush.
12. Finish and polish in the usual manner.

Kimwipe is not a trademark of DENTSPLY International.

VIII. NIGHTGUARD AND OCCLUSAL SPLINT

1. Block out any undercuts on the model that could potentially lock the splint in the mouth.
2. Follow I. FULL DENTURES A. Baseplate Preparation, steps 1-5, then the following steps:
3. Place the rounded side of the Clear Baseplate Resin arch on the incisal/ occlusal surface of the model and carefully adapt the resin to the desired shape. Trim excess to height of tooth contour with DENTSPLY Electric Spatula set at 230-240°F (110-116°C).
4. When adapting resin to articulated models, apply Eclipse Model Release Agent to the incisal and occlusal surface of the opposing model. With a light tapping motion, carefully close the articulated models together. Move the articulator through excursive positions to open interfering pathways in the appliance. Open articulator and remove model.
5. Follow I. FULL DENTURES A. Baseplate Preparation, steps 7-9, except cure the appliance in the Eclipse Processing Unit using

Menu #7.

6. Finish and polish in the usual manner.

NOTE: The appliance should fit passively on the working model.

CLINICAL DIRECTIONS FOR USE

Indications for Use: Full and partial removable dentures, occlusal splints, nightguards and provisional partial dentures.

Contraindications for Use: See precautions and warnings.

Warnings: This product contains polymerizable materials. It does not contain methyl methacrylate (MMA) monomer. It is possible that some dental health care workers may develop allergy or sensitivity to the product. Wash thoroughly with soap and water after contact. If dermatitis or other symptoms persist, seek medical advice.

Mechanical tooth retention must be used in accordance with instructions to prevent debonding.

May not be appropriate for partial dentures where limited vertical dimension is present, or complete dentures with less than 12mm interarch space (6mm for single denture).

Precautions:

1. This product is intended to be used only as specifically outlined in the Directions For Use. Any use of this product inconsistent with the Directions For Use is at the discretion and sole responsibility of the practitioner.
2. Insufficient clinical data exists at present for use with implant supported removable prostheses and

removable appliances employing attachments.

3. Users with special skin problems, cuts or abrasions may choose to wear protective gloves.
4. When grinding prosthodontic resins, proper ventilation, masks and vacuum systems should be used.
5. Since the Eclipse resins will become the final prosthesis, cleanliness is very important to an acceptable outcome. Wash hands before using uncured Eclipse materials.
6. Store at 60-80°F (16-27°C), away from moisture and in a light-safe container. The materials may polymerize prematurely if stored at excessively high temperature or if exposed to a visible or UV light source. Do not expose to direct sunlight.
7. Unused materials should be polymerized prior to disposal.
8. The following disinfectants have been used successfully with cured and uncured Eclipse resins:
 - Glutaraldehyde-based disinfectants
 - Quaternary ammonium chlorides
 - Quaternary ammonium chlorides /isopropyl alcohol (low concentration)
 - Spray-based phenolics
 Use of other disinfectant products may not be suitable with Eclipse. Avoid disinfectants containing iodine.

Adverse Reactions:

1. Allergic contact dermatitis and other allergic reactions may occur in susceptible individuals.
2. Dust will be generated when grinding these materials. Eye, skin and respiratory irritation may occur if appropriate engineering controls are not used.

13

STEP-BY STEP INSTRUCTIONS:

Because Eclipse dentures are fabricated using a processed record base, it is necessary to identify the vibrating line and posterior palatal seal configuration at the final impression appointment. Also, the denture base shade selection must accompany this final impression/master model.

A. Full and Partial Dentures

1. For full dentures, your dental laboratory will likely send the Eclipse Baseplate with wax occlusion rims. This baseplate will become part of the final denture. Seat the denture and evaluate the fit of this baseplate, using pressure indicating paste if necessary. If the borders are short, or the baseplate is not stable, extend borders with compound and take a wash impression. Adjust wax rims for proper lip support and VDO, mark mid-line, high smile line, canines; then return to dental lab.
2. When a full or partial denture set-up is received from the dental laboratory, it is ready for try-in. The denture base is uncured and has an extended working time versus traditional light cured materials. Even so, the "wax-up" should be protected from light as much as possible. Use the light safe bags provided.
3. Try the case in per your normal procedures, with one exception: **DO NOT USE A FLAME TO ADJUST.**
4. If anterior teeth need minor adjustments in position, this can be accomplished chairside. Prepare a water bath at 110°F (43°C). Place the set-up into the water for approximately 30-45 seconds to warm the uncured denture base material. This time may vary, depending on thickness of the set-up.
5. Significant modifications of occlusal plane, vertical dimension

or anterior tooth position will require re-contouring of the external base material. In most instances, this is easier and more efficient for the dental laboratory to accomplish. If necessary, return the set-up to the laboratory for re-contouring.

6. Reposition the anterior teeth as necessary by pushing on the denture base material as much as possible to change tooth position. Attempting to move the teeth simply by pushing on the tooth will open gaps between the tooth and the denture base material.

Note: It is important to notify the laboratory of any adjustments made.

7. The posterior teeth should not be repositioned intra-orally. Should the vertical dimension require adjustment, remove the posterior teeth on one arch and replace with a silicone bite registration material such as Regisil®.
8. To build up the gingival areas during try-in, use a silicone putty material such as Reprosil®.
9. Request and complete a second try-in before processing if necessary or desired.
10. When the set-up is acceptable, the case may be returned to the laboratory for processing and finishing/polishing.

B. Relines

Follow normal clinical procedures for a laboratory processed reline.

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SECTION IV - FIRE AND EXPLOSION HAZARD DATA

4.1 Flash Point (Method Used) > 220°C (428°F)
 4.2 Flammable (Explosive) Limits
 LEL: NA UEL: NA
 4.3 Extinguishing Media Dry chemical, CO₂, foam, water spray, and water fog
 4.4 Special Fire Fighting Procedures and Equipment: Fireman should wear protective clothing and self-contained respirators if large amounts of material are involved in a fire.
 4.5 Unusual Fire and Explosion Hazards: Irritating gases may be emitted when combusted.

Eclipse Denture Base Resin [Baseplate, Setup and Contour Resin], cont'd

SECTION V - HEALTH HAZARD DATA

5.1 Effects of Acute and Chronic over Exposure: Skin irritation is possible in individuals sensitive to (meth)acrylates and contact dermatitis (skin sensitization) can develop in some individuals.
 5.2 Emergency and First Aid Procedures: Eyes: If material contacts eyes, flush eyes immediately with copious amounts of tap water for a minimum of 15 minutes and get medical attention. Skin: If skin irritation occurs, wash irritated area with soap and water. Remove clothing if it becomes contaminated with resin. Sensitization: To minimize chance of becoming sensitized, wash hands often when using material. Particularly sensitive individuals should avoid skin contact.

SECTION VI - REACTIVITY DATA

6.1 Unstable: Stable: X
 6.2 Conditions to avoid (Stability): Material can polymerize if exposed to high temperatures or to direct sunlight.
 6.3 Incompatibility (Materials to avoid): Oxidizing agents.
 6.4 Hazardous Decomposition Products: Acrid fumes may result if large amounts are combusted during a fire.
 6.5 Hazardous Polymerization: May Occur: Will Not Occur: X
 6.6 Conditions to avoid polymerization: Direct sunlight, fluorescent lighting, high temperatures for extended time period.

SECTION VII - SPECIAL PRECAUTIONS AND SPILL OR LEAK PROCEDURES

7.1 Handling and Storage Precautions: This material should be handled and used only in a manner consistent with the product instructions. Do not store at high temperature or in sunlight. Store away from food.
 7.2 Other Precautions: Visible light can cause polymerization.
 7.3 Steps to be taken in case material is released or spilled: No Special procedure required.
 7.4 Waste Disposal Procedures: Polymerize baseplate and clear resin by exposing material to sunlight until solid, then dispose of as a non-hazardous waste. Small (gram / oz.) quantities of set up resin and contour resin can be disposed of as non-hazardous waste.
 7.5 Neutralizing Agent: NA

SECTION VIII - OCCUPATIONAL PROTECTION MEASURES

8.1 Respiratory Protection: Use a dust mask while grinding the cured material to avoid breathing dust.
 8.2 Ventilation: Use local exhaust when grinding cured material.
 8.3 Protective Gloves: For users with skin problems, cuts or abrasions - wear disposable protective gloves.
 8.4 Eye Protection: Use safety glasses while grinding.
 8.5 Other Personal Protective Equipment: None required

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© 2002, 2004 DENTSPLY International, Inc. All rights reserved.