# Orthoresin

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## Directions for Use

## Mode d'emploi

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## Mode de empleo

Made in U.K.

DeguDent GmbH

Postfach 1364 · D-63403 Hanau  
GERMANY  
Telefon +49 61 81 / 59-50

50535084 / 5000 / 0503 / KH

(D)

# Orthoresin

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(GB)

# Orthoresin

**Self-curing orthodontic resin for the preparation of removable orthodontic appliances.**

## DIRECTIONS FOR USE

- Paint a newly cast model with Isolant and allow to dry for two minutes. If models have dried out they should be immersed in warm water for approximately five minutes. When dry, proceed as described previously.
- Dispense a little monomer into the plastic pot.
- Using the brush, apply monomer to the model surface. Dispense a thin uniform layer of powder onto one side of the model from the squeeze bottle. Add more liquid with the special dropper and saturate all the dry powder particles. Build up further layers until the required thickness is achieved. Smooth the surface with the brush dipped in monomer. Repeat this process for the other side of the appliance. Upon completion use the brush again to provide an overall smooth surface.  
*Note:* Always clean the brush thoroughly after use. Never return monomer to the bottle once it has been dispensed.
- Polymerisation  
For optimum results it is recommended to cure Orthoresin in a pressure flask containing water at 45 – 50 °C at approximately 2 bar (approx. 30 psi). The curing takes place in 8 to 10 minutes.  
*Note:* Cribs, expansion devices etc. should be held in position with wax or plaster. Finger springs should also be boxed out similarly.
- Trimming and polishing should be carried out by the normal methods.
- Repairs may be carried out by further addition of material using the above technique.

## STORAGE

Store in a cool, dry place (max. 25 °C).  
Keep containers closed when not in use.  
Avoid direct sunlight.

## SAFETY WARNINGS

### Engineering measures

Ensure suitable suction / aeration at the work place and with operational machinery.

### Personal protective equipment

#### Hand protection

Wear protective gloves made of the following materials: butyl rubber.

#### Eye protection

Close-fitting protective goggles (e.g. closed goggles) recommended

#### Skin and body protection

Immediately change moistened and saturated work clothes. Use suitable skin protective agents before handling product. Cleanse and apply cream to skin after work. Preventive skin protection recommended.

Do not use in case of hypersensitivity to any of the components of the product.

For dental use only.

The batch number and expiry date are shown on all containers. Please quote both in any communication. Do not use after expiry date.

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