D

DENSPLY

Selectaplus

MATERIAL

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Selectaplus

MATERIAL

Cold curing denture base resin consisting of one powder based upon polymethylmethacrylate and two liquids (Selectaplus liquid CN and Selectaplus liquid CE) based upon methylmethacrylate.

The initiator system is free of tertiary amine.

APPLICATIONS

For fabricating full and partial dentures, relines, repairs.

PROPERTIES

The resin is absolutely color stable and has very good physical properties. When liquid CN is used, Selectaplus has a short viscous mouldable phase. When liquid CE is used Selectaplus has a longer viscous mouldable phase (up to 13 minutes).

CONTRAINDICATION

Do not use this product for preventive reasons if patient has shown hypersensitivity to any of the components it contains.

SAFETY WARNINGS

for Selectaplus Liquid CN and CE.

Highly Inflammable

Irritant

Liquid contains
methylmethacrylate.

R 37/38	Irritating to respiratory system and skin.
R 43	May cause sensitization by skin contact.
S 24	Avoid contact with the skin.
S 37	Wear suitable gloves.
S 46	If swallowed, seek medical advice immediately and show this container or label.

For dental use only.

STORAGE INSTRUCTIONS

Store in a cool, dry place (max. 25 °C).
Keep containers closed when not in use.
Avoid direct sunlight.

During storage the liquids may show a slightly yellow color; this does not influence the color of the mixed product or its quality.

DIRECTIONS FOR USE

Preparations:

The teeth are set up in wax and are fixed by a silicone or plaster key. After scalding the wax and after vaporizing, insulate with insulation fluid in still warm condition. The insulation must be dried off, allow no residue to remain in the crevices. It is recommended to roughen the basal areas of the teeth and to spread them with liquid CN or liquid CE.

Mixing:

Mixing ratio powder/liquid: 5 g / 3,5 ml.
Free dosage is possible.

1. Pour liquid CN or liquid CE into a mixing bowl.
2. Add powder gradually (smoothly).
3. Mix immediately.

Processing Times:

The processing times are dependent on room temperature and the quantity of material mixed. The following time table relates to 23 °C, an amount of 5 g powder and 3,5 ml liquid.

Mixing time:	approx. 0,5 minutes
Pouring phase:	up to 3 minutes
Stable phase:	after 4 minutes

Viscous mouldable phase:

Liquid CN:	up to the 8. minute
Liquid CE:	up to the 13. minute

After stable phase has been reached, the material does not flow anymore; during viscous mouldable phase it can be moulded.

Polymerization:

Liquid CN: 40 °C, at least 2–3 bar, 10 min., polymerize after the 8. minute at the latest.

Liquid CE: 50 °C, at least 2–3 bar, 10 min., polymerize after the 15. minute at the latest.

Repairs:

For repairs separate the fracture surfaces by approximately 1.5 mm. Bevel the edges to a shallow "V" (approximately 120° angle). Apply denture base liquid on trimmed surfaces.

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