Triad® VLC-Löffelmaterial

Verarbeitungsanleitung
Instructions for Use
Mode d'emploi
Bruksanvisning
Istruzioni per l'uso
Instrucciones
para la elaboración
Gebruiksaanwijzing

DENTSPLY

Triad VLC-Löffelmaterial

Triad VLC-Löffelmaterial

VLC Custom Tray Material

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Instructions

MATERIAL

Triad® VLC tray/base plates are available in pink and neutral and consist of a matrix of multifunctional acrylates and inorganic fillers.

PROPERTIES

Triad® VLC tray/base plates exhibit very

- minimal shrinkage,
- minimal water uptake and
- high dimensional stability

Processing

1. Block out undercuts with wax. To fabricate a tray, place a sheet of wax on the model to create a spacer. Extend the wax to within 2–3 mm of the tray periphery. Stops may be cut into the wax on the alveolar ridges.

2. Apply either Triad® Model Release Agent (MRA) or an alginate separating agent such as Isolant.

3. Unpack the Triad® VLC tray/base plate, adapt it to the model carefully and trim the peripheries. The excess material is easily used to form a tray handle.

4. Coat the entire tray or base plate with Triad® ABC varnish.

5. Polymerize the tray or base plate on the model in a light curing unit. Polymerization time in a Triad® light curing unit: 5 minutes.

6. Release the tray or base plate from the model, coat the underside with Triad® ABC varnish and polymerize again for 5 minutes.

7. Brush off the Triad® ABC varnish under running water, dry the tray or base plate and trim as usual with tungsten carbide cutters and sandpaper.

Please note:

Use a dust extractor. Avoid contact with skin. Wear suitable gloves.

STORAGE

Do not store at temperatures exceeding 23 °C. Higher temperatures affect the handling properties.

For dental use only.

Contraindication:

Do not use this product for preventive reasons if patient has shown hypersensitivity to any of the components it contains.

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DeguDent GmbH