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Cergo[®] fit

Instructions for use
Mode d'emploi
Istruzioni per l'uso
Instrucciones para el uso

Cergo[®] fit

D	Gebrauchsanweisung	3
GB	Instructions for use	7
F	Mode d'emploi	11
I	Istruzioni per l'uso	15
E	Instrucciones para el uso	19

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Cergo fit is a phosphate-bonded special investment compound for the production of moulds for pressable ceramic. Cergo fit does not contain any plaster or graphite. Cergo fit is mixed with the **Special Liquid (REF 53 2540 0114)**.

Cergo fit may be heated in the traditional way with a holding time at 270 °C and at final temperature as well as by the rapid-heating method. Please adhere to the information provided on this subject under item 6.1.

1 Safety instructions

When using this product, mixing with other liquids or materials must be avoided, since health hazards as well as a reduction of the compression quality can not be excluded.

- This product contains quartz and cristobalite, some in the form of fine dust that may penetrate the pulmonary alveoli. Inhalation of the dust constitutes a potential carcinogenic hazard. Ensure that an adequate suction line/ventilation is available in the workspace and near the machines.
- Inhalation the dust of the product may trigger silicosis.
- Employers must provide Class P2 respiratory protection.
- Heating the investment compound above 200°C will result in the formation of ammonia. Ammonia irritates the respiratory organs.
- May cause mechanical irritation of the skin and mucosa of the eye and respiratory passages.
- Irritation and caustic reactions to the skin and mucosa must be expected on local contact.

Further information is contained in the safety data sheet.

Information for the processing of mixing liquids for dental investment compounds:

If contact is made with the skin, immediately wash off with a large amount of water. In case of contact with the eyes, wash them out immediately with water and seek medical attention.

Mixing liquids are alkaline.

2 Recommendations for processing

Processing parameters	Recommended value
Material temperature	21 °C
Mixing ratio	100 g : 21 ml
Stirring time under vacuum	60 s
Working time range	5 – 8 min

3 Important instructions

Prior to using the product, please adhere to the safety recommendations provided under item 1. Cergo fit must not come in contact with plasters or plaster-containing materials to avoid damage to the pressable ceramic (Cergo). Accordingly, separate mixing bowls and clean spatulas must be used for mixing Cergo fit.

Store Cergo fit under cool and dry conditions. If proper storage is ensured, the investment compound can be stored for 18 months. The special liquid **must be protected against frost** and stored at a room temperature between 5 °C and max. 30 °C!

4 Expansion control

Cergo fit is mixed with special liquid. By diluting this liquid with distilled water the expansion can be controlled so that perfect fit for any dimension is achieved.

Typical concentrations of the special liquid:

Indication	Concentration
Inlay (MOD) / Veneer	85 – 90 %
Posterior crown	70 – 80 %
Anterior crown	80 – 90 %

5 Mixing the investment compound

Premix powder and liquid in a clean bowl of the Multivac® stirring and evacuating unit for approx. 30 s (incl. the time for manual adding into the unit) using a clean spatula. Let the investment compound stand under vacuum for approx. 30 s. Finally, mix under vacuum in the Multivac for 60 s.

Please use program number 1 for mixing with the Multivac compact (Deguvest CF).

The following information must be adhered to when mixing investment compounds:

The mixing ratio and the mixing time are important parameters that ensure reproduceable results. It is recommended to adhere precisely to the mixing ratio and the mixing time provided. High vacuum in the mixing bowl is a factor that determines the quality; the mixing process should only be started when full vacuum is displayed.

5.1 Investing

The muffle ring to which vaseline has been applied is placed on the muffle former. The mixed Cergo fit investment compound is filled into the muffle and the formation of bubbles is avoided; then the muffle measuring device is placed on. The muffle measuring device ensures the correct height of the muffle and

correct angle position of the muffle towards the pressing direction of the injecting furnace. Muffle ring and muffle former can be removed after **40 min**.

Check safe position of the muffle on its base. Any high spots must be removed.

6 Wax elimination and preheating

After the setting time, the muffle and the clean compression plunger which is free from any previously pressed ceramic residues are placed into the preheating furnace at room temperature. Simultaneous preheating of muffles with metal ring should be avoided to prevent contamination of the muffles and the compression plunger. Heat the preheating furnace up to **270 °C** using a temperature increase rate of **7 °C/min** and maintain the temperature for **40 min**. Then heat up to a final temperature of **850 °C** using the same increase rate and heat the muffle at this temperature for **40 min**.

In order not to affect the result by incompletely melted wax in the muffle, the wax must be eliminated immediately after the investment compound has hardened.

Parameters such as various modelling materials (e.g. wax, resin, coping), geometric form and size of the object, muffle size and pressable ceramic, which influence the compression result, must be considered individually as this may have a decisive influence on the quality of the compression.

Holding times do not change these parameters!

6.1 Rapid heating

Due to its material properties Cergo fit is also suitable for the so-called rapid heating method. The special liquid is also used for this method. To minimize the risk of crack formation during considerable thermal stress, we recommend the exclusive use of a **100 g** Cergo flask for compressing according to the rapid heating method. To achieve reproduceable casting results using this heating process, the following instructions must be adhered to:

After the investment compound has hardened (**not earlier than after approx. 15 min but not later than after 20 min**) the muffle must be immediately placed into the preheating furnace that has been heated up to **850 °C**. Please ensure that the muffle together with the compression plunger is placed into the furnace at a slightly **tilted** position and that the compression channel is visible.

To achieve homogeneous heating of the muffle, the holding time of **40 min** at final temperature must be adhered to.

Complicated compression objects should not be heated in the rapid heating method.

7 Pressing Procedure

The pressing procedure must be carried out according to the instructions by the manufacturer of your pressable ceramic ingots.

8 Devesting

After pressing, let the muffle cool down to room temperature in air and devest. For easy devesting mark the position of the compression objects with the aluminium oxide compression plunger. Cut deeply into the investment material using a separating disc. Then carefully split the muffle with a plaster knife. Sandblast at a pressure of **4 bar** until the pressed objects are visible; then devest the objects at a pressure of approx. **2 bar**. Use only glass beads, grain size **50 µm** for devesting.

9 Summary of the most relevant data

Processing parameters	
Material temperature	21 °C
Mixing ratio	100 g : 21 ml
Stirring time under vacuum	60 s
Working time range	5 – 8 min
Setting time	10 – 15 min
Removal of ring and former	40 min
Expansion control	
Inlay (MOD) / Veneer	85 – 90 %
Posterior crowns	70 – 80 %
Anterior crown	80 – 90 %
Wax elimination and preheating	
Temperature increase rate	7 °C / min
Holding time	270 °C / 40 min
Final temperature	850 °C / 40 min
Rapid heating	
Removal of ring and former	15 – 20 min
Final temperature	850 °C / 40 min

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DeguDent GmbH