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Deguform plus

Instructions for use
Mode d'emploi
Istruzioni per l'uso
Instrucciones de uso

DeguDent

Deguform® plus

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Deguform® plus

1 Description

Deguform plus is a high-quality addition-curing silicone duplicating compound for producing secondary models in dental technology.

The two-component silicone is suitable for all modern duplicating techniques and convinces by:

- **excellent reproduction of detail**
- **form-stable accurate reproduction**
- **very good tear resistance**
- **flow-viscosity consistency**
- **firm and smooth model surfaces**
- **simple and safe handling by 1:1 dosage**
- **time saving opposite gel duplication**

2 Processing notes

2.1 Deguform plus is particularly suitable for duplicating without rings.

2.2 Do **not** water the master model but duplicate dry and clean.

Note:

Silicone duplicating compounds are sensitive opposite catalyst poisons. These can be: e.g. insulating means, cyanoacrylate, plastics. Curing can thereby be delayed or prevented.

Lubricating the duplicating flask can also result in a non-hardening smear layer.

2.3 Mixing Deguform plus

Deguform plus is in principle mixed in a ratio of 1:1 (Component A, white: Component B, blue).

Deviations can possibly affect curing as well as the physical properties.

It is possible to mix Deguform plus using different methods:

- **by hand**
- **with suitable dosing units**
- **with the Multivac® stirring units**


2.3.1 Mixing by hand

Fill Deguform plus in a ratio of 1:1 into a dosing and mixing bowl and mix thoroughly until a homogeneous color results. Then immediately pour the compound in thin stream into the duplicating flask.

2.3.2 Use of dosing devices

The minimum filling level must be observed during the use of 6 kg cans in connection with mixing and dosing units. Non-observance during processing results in air being fed into the system, a homogeneous mixture is no longer ensured. In such a case wait after connecting a new can for at least 45 minutes prior to the next withdrawal.

The following illustration shows the minimum silicone level in the can:

2.3.3 Mixing in Multivac® units

Have a separate mixing bowl ready for mixing silicone duplicating compounds. You achieve the best results with a new, unscratched bowl.

Fill Deguform plus in a ratio of 1:1 into the Multivac mixing unit, briefly stir thoroughly and mix in the Multivac according to the data in the following table. Then immediately pour the compound in a thin stream into the duplicating flask. The duplicating silicone can easily be taken out of the mixing bowl and off the stirring blade after curing.

You must follow the recommended mixing parameters!

These have been determined so as to protect your equipment and to obtain optimal mixing results. Exceeding the maximum recommended **filling quantity** can in particular cause cementing of the filter and in extreme cases a drawing-in into the vacuum pump.

Device type	Bowl size	Max. filling quantity	Vacuum	Stirring speed	Stirring time
Multivac S, 4	Small	100 g	On	Presetting	60 s
Multivac S, 4	Large	400 g	On	Presetting	60 s
Multivac Eco	Small	100 g	On	Presetting	75 s
Multivac Eco	Large	400 g	On	Presetting	75 s
Multivac compact	Small	100 g	90 %	300 1/min	60 s
Multivac compact	Large	400 g	90 %	300 1/min	60 s

We recommend starting the stirring procedure only after achieving the vacuum. Longer evacuating does not affect the mixing result.

2.4 The silicone is cured after at least 30 minutes (at room temperature). The master model can now be removed carefully from the silicone form (if necessary with compressed air).

NOTE!

Removing the master model causes a deformation of the silicone form. This is neutralised after some minutes by the inherent resilience. We therefore recommend to only begin producing the investment model after a 15 minute waiting period.

2.5 The silicone form is highly elastic and keeps the precise form permanently.

Moisten the silicone form with Waxit® and blow dry before filling with the investment compound. The investment compound flow is thereby facilitated on the silicone surface. Make sure that the silicone surface is dry and completely set.

All kinds of dental gypsum and model plastics are suitable for filling additionally to investment compound.

2.6 Deguform plus can be disinfected both chemically and thermally.

3 Technical data

Component A	white
Component B	dark-blue
Final colour	homogeneous intensive blue
Mixing ratio	1:1
Mixing time (manual)	approx. 45 s
Processing time *	approx. 5 – 7 min
Demoulding capacity*	approx. 30 min
Final hardness (shore A)	28

* Room temperature of approx. 23 °C. The time data can change at extreme room temperatures.

4 Storage conditions

Deguform plus should be stored at room temperature; only then can a minimum durability of 2 years starting from date of manufacture be ensured.

5 Product range

Deguform plus silicone duplicating compound, Component A and component B	REF
each in a bottle with 860 ml (1,000 g) contents	2511 0502
each in a can with 5,170 ml (6,000 g) contents	2511 0512

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DeguDent GmbH