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Deguvest[®] Impuls

Instructions for use
Mode d'emploi
Istruzioni per l'uso
Instrucciones para el uso

Umschlags. 2

Vakat !

Deguvest® Impuls

D	Gebrauchsanweisung	5
GB	Instructions for use	11
F	Mode d'emploi	17
I	Istruzioni per l'uso	23
E	Instrucciones para el uso	29

Seite 4

Vakat !

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Deguvest Impuls is a carbon-free, phosphate-bonded and plaster-free universal investment material which is suitable for the entire spectrum of precious metal precision casting.

Deguvest may be used for the conventional as well as for high-speed heating. For this purpose the Special Liquid must be used. Carefully read the information provided under 5.3.

1 Technical data

Deguvest Impuls

Mixing ratio powder : liquid	100 g : 23 ml
Total expansion linear	1.2 – 2.4 %
Mixing time under vacuum	60 s
Working time range	4 – 6 min
Compressive strength (depending on concentration of the mixing liquid)	4 – 8 MPa

Mixing liquid composed of Distilled water	Deguvest mixing liquid	Concentrate	Setting expansion	Thermal expansion
3 parts	1 part	25 %	0.30 %	0.90 %
2 parts	2 parts	50 %	0.60 %	1.00 %
1 part	3 parts	75 %	0.90 %	1.10 %
	4 parts	100 %	1.15 %	1.25 %

Setting expansion + thermal expansion (mean values) = total expansion

2 Important information

Be sure to observe information concerning safety prior to use under sections 6 and 6.1.

Deguvest Impuls should not be in contact with plasters or plaster-containing materials since this may lead to sulphur damage of high-melting alloys. Therefore special mixing cups and spatulas should be used for mixing Deguvest Impuls.

Deguvest Impuls must be stored under cool and dry conditions, to ensure a shelf life of 18 months.

Deguvest Liquid and Special Liquid must be stored at a room temperature between 5 °C and max. 30 °C.

Typical concentrations of Deguvest Liquid for investing objects in Deguvest Impuls

Please note that the Special Liquid must be used for rapid heating (see 5.3).

Object:	MOD-Inlay	Anterior crown	Anterior crown	Molar crown	Molar crown	6° cone secondary*	6° cone secondary*	Telescope secondary*	Telescope secondary*	RS circular*	8-unit bridge
Modellation material:	Wax	Wax	Adapta	Wax	Adapta	Adapta	Pattern	Adapta	Pattern	Pattern	Wax
Alloys: examples in brackets ()											
High-gold content C + B alloys (Degulor M)	60–70%	65–75%	65–75%	65–75%	65–75%	50–60%	40–50%	60–70%	55–65%	50–60%	65–75%
Reduced C + B alloys (Stabilor NF IV)	60–70%	65–75%	65–75%	65–75%	65–75%	50–60%	40–50%	60–70%	55–65%	50–60%	65–75%
High-gold content ceramic bonding alloys (Degudent U)	60–70%	65–75%	65–75%	65–75%	65–75%	50–60%	40–50%	60–70%	55–65%	50–60%	65–75%
Reduced ceramic bonding alloys (Deva 4)	–	65–75%	65–75%	65–75%	65–75%	–	–	–	–	55–65%	70–80%
Palladium-based alloys (Degupal G)	–	70–80%	70–80%	70–80%	70–80%	–	–	–	–	60–70%	75–85%

* We do not recommend high-speed heating for preheating of complex restorations.

3 Expansion control

Deguvest Impuls is mixed with Deguvest Liquid or for high-speed heating with Special Liquid. When these liquids are diluted with distilled water, the expansion can be controlled to compensate the contraction of the casting alloy used.

The higher the concentration of the mixing liquid, the greater the total expansion of the investment.

4 Preparing the casting ring

- Dilute Deguvest Liquid to the desired concentration by adding **distilled** water.
- Line the casting ring with a moistened DeguDent casting ring liner or Deguvest Vlies. The end of the casting ring liner should be flush with the end of the casting ring opposite the sprue base former.
- Place the wax pattern and fix it with wax on the DeguDent sprue base former.
- Coat the wax model with the wetting agent Waxit® using a brush or by spraying and remove excess of wetting agent (do not use compressed air for this purpose).

- Measure Deguvest powder and liquid in the ratio of 100 g powder : 23 ml mixing liquid resp. according to the following table:

Casting ring	Deguvest Impuls powder	Deguvest mixing liquid
1 x	1 x 60 g = 60 g	14 ml
3 x	1 x 150 g = 150 g	35 ml
6 x	2 x 150 g = 300 g	70 ml
9 x	3 x 150 g = 450 g	105 ml

4.1 Mixing the investment material

Premix the powder and liquid for approximately 30 s by hand with a spatula in a clean bowl supplied with the Multivac® vacuum mixer. Allow the investment to remain under vacuum for about 30 s.

Finally, mix under vacuum with the Multivac for about 60 s.

When investing with the Multivac® vacuum filling bowl, premix the powder and the liquid for about 30 s by hand using a spatula. Next, mix the investment material for about 60 s under vacuum in the Multivac. Then allow the investment material to flow into the casting ring under light vibration within 30 s.

The vibrator in the Multivac® 4 is particularly well suited for this purpose. Finally, open the ventilation valve.

Please use program number 4 for mixing with the Multivac® compact (Deguvest Im).

During mixing of the investment, the following points, must be observed:

- High degree of vacuum in the mixing bowl determines the quality. Be sure to wait until the full vacuum is indicated before start mixing.
- The mixing ratio influences the material characteristics; increasing the amount of liquid has the following effects:
 - thinner consistency
 - extended working time range
 - slightly lower setting expansion
 - slightly lower strength;
 a small amount of liquid added has the following effects:
 - smoother surfaces of the casting
 - somewhat higher expansion;
- Mixing ratio that differs from the recommended range:
 - possible crack formation
 - rougher surface of the casting.

- The mixing time is important for consistent results; strictly adhere to the recommended mixing time.
Mixing times shorter:
 - uncontrolled setting expansion and strength
 - air bubbles can remain in the investment and result in surface bubbles;
 Longer mixing time:
 - lower setting expansion
 - shorter setting times.

4.2 Investing

Let the mixed investment material flow into the casting ring under light vibration and carefully avoid the formation of bubbles. Immediately turn off the vibrator after the casting ring has been filled.

4.3 Hardening

Let the casting ring harden for 20 – 50 min (see table, Heat treatment).

5 Preheating and Casting

5.1 Wax elimination

Place the casting ring in a furnace preheated to about 300 °C and eliminate the wax in about 30–60 min (see table, Heat treatment).

To avoid influences of the casting, burn out possible wax residues immediately after the investment has set.

5.2 Preheating

Bring the casting ring to the final temperature required for the particular alloy. Maintain the final temperature for a minimum of 20 to 60 min (see table, Heat treatment).

For the manufacturing of complex restorations, e.g. large bridges, suprastructures or attachments, a heat rate of **7 °C/min** is recommended.

The parameters that influence the casting such as wax, resin, adapter cap etc., shape and size of the restoration, size of the casting ring and type of alloy must be regarded individually. All of those factors will influence the final result.

– **Holding times do not change these factors.**

5.3 Rapid heating!

Due to its material characteristics Deguvest Impuls is suitable for high-speed heating.

For this purpose, the Special Liquid must be used. In order to obtain consistent results, the following information is important:

Heat treatment	Size of the casting ring			
	1 x	3 x	6 x	9 x
Setting/drying time	20 min	30 min	40 min	50 min
Wax elimination temperature	300 °C **			
Duration of wax elimination	30 min	40 min	50 min	60 min
Initial temperature of the electric preheating furnace	300 °C **			
Duration of preheating after final temperature is reached *	20 min	30 min	45 min	60 min

* alloy-specific

** When the casting rings are placed in a preheating furnace preheated to 300 °C, a slow distribution of heat takes place inside the casting ring because of the thermal conductivity of the investment material.

- The use of a **metal casting ring is required**.
- The maximum **casting ring size is 6 x** and **must not be exceeded**.
- Line the casting ring with a **ceramically bonded casting ring liner**, e.g. Deguvest Vlies.
The casting ring liner should be flush with the upper rim of the casting ring.
- For high-precious alloys dilute the Special Liquid to 50 %.
- For palladium-based alloys dilute the Special Liquid to 75 %.
- After the investment has set (**after approx. 15 min**), **the casting ring should immediately** be placed in the furnace which has been preheated to the required preheating temperature (700 – 900 °C). Please ensure that the casting ring is **not placed upright in the furnace**.
- In order to achieve homogeneous heating of the casting ring, the holding time at the final temperature is essential (see table, Heat treatment).

Take note that due to high-speed heating a higher expansion of the investment is expected.

Normally, this will lead to larger casting objects. High-speed heating should not be used for complex restorations.

5.4 Casting

Cast according to the instructions of the manufacturer of the alloy.

5.5 Devesting

After casting, let the casting ring cool down to room temperature and remove the investment.

Soak the investment ring in water for approximately 15 min. Carefully split it open at several points with plaster pliers. When this method is used, the dental casting can be easily removed from the investment, and dust is avoided. Final residues of investment can be removed or carefully blasted out with aluminium oxide, grain size 100 – 150 µm.

Do not use a hammer for devesting!

6 Safety information

Mixing with other liquids or materials during this process health hazards as well as impairments to the casting quality cannot be excluded.

All dental investment compounds contain parts of crystalline quartz and/or cristobalite.

These substances can adversely affect health if they are breathed in. Can cause silicosis. Always use a sufficiently powerful extractor or wear a mask.

Carry out the “coarse” removal of the casting under water or soak the investment material mold.

6.1 Information for the processing of mixing liquids for dental investment materials

Following contact with skin:

Upon skin contact, wash with plenty of water.

Following eye contact:

Rinse thoroughly with plenty of water keeping eyelid open.

In case of persistent discomfort: Consult an ophthalmologist.

Mixing liquids are alkaline.

The Safety Data Sheet includes further details.

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