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# Deguvest<sup>®</sup> L

Instructions for use  
Mode d'emploi  
Istruzioni per l'uso  
Instrucciones para el uso

# Deguvest<sup>®</sup> L

- D** Gebrauchsanweisung . . . . . 3
- GB** Instructions for use . . . . . 6
- F** Mode d'emploi . . . . . 8
- I** Istruzioni per l'uso . . . . . 11
- E** Instrucciones para el uso . . . . . 13

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Deguvest L is a precision soldering investment developed on a scientific basis. The total expansion of Deguvest L is matched to the dimensional behavior shown by dental precious metal alloys during the soldering process. Deguvest L therefore produces soldered objects with extremely high accuracy of fit.

Deguvest L has outstanding strength. The compound can be shaped well and easily cut, and has high thermal stability. For this reason, Deguvest L can also be used for soldering of metal-ceramic works with high melting alloys.

## 1 Safety information

Mixing with other liquids or materials during processing health hazards as well as impairments to the casting quality cannot be excluded.

- This product contains quartz and cristobalite, some in the form of fine dust that may penetrate the pulmonary alveoli. Inhalation of the dust constitutes a potential carcinogenic hazard. Ensure that an adequate suction line/ventilation is available in the workspace and near the machines.
- Inhalation the dust of the product may trigger silicosis.
- Employers must provide Class P2 respiratory protection.
- Heating the investment compound above 200°C will result in the formation of ammonia. Ammonia irritates the respiratory organs.
- May cause mechanical irritation of the skin and mucosa of the eye and respiratory passages.
- Irritation and caustic reactions to the skin and mucosa must be expected on local contact.

Carry out the “coarse” removal of the casting after the casting operation under water or soak the investment material mould.

The Safety Data Sheet includes further details.

## 2 Technical data

Mixing ratio:	depending on the desired consistency, recommended mixing ratio 100 g powder : 20–24 ml water
Total expansion up to 700 °C:	approx. 1.2 % linear
Compressive strenght after 2 hours:	3.5 MPa
Processing time span:	3–4 min

### 3 Instructions for use

- 3.1 Roughen and clean the surfaces of the object to be soldered. The soldering gap should be parallel and should be  $\leq 0.2$  mm.
- 3.2 Pay attention to correct position of the object on the model. Close the soldering gap with dental wax or modelling wax. Then fix the object with adhesive wax and/or with acrylic which burns without residue. Objects fused with dental porcelain must have no contact with the solder investment compound. Cover porcelain, especially towards the cervical margin and to the basal surface with wax. After the wax has been boiled out, a space is left between the solder investment compound and the porcelain in this way.
- 3.3 Fill Deguvest L and water (distilled water) in the recommended mixing ratio into a clean mixing bowl and thoroughly mix or knead with clean spatula.
- 3.4 Form the solder investment compound into bars and place the object to be soldered with slight hand vibration carefully into the investment compound. Keep the soldering gap free of investment to avoid sulfur damage with sensitive alloys.
- 3.5 Cut back the soldering block. Remove the adhesive wax or acrylic and boil out the soldering gap. Pay attention when boiling out that no adhesive wax is swept into the soldering gap.  
**Do not** remove the wax with a steam cleaning unit.
- 3.6 Place fluxing agent in the solder gap and on any adjoining metal surface. When soldering a ceramically veneered workpiece, use a suitable fluxing agent (such as DS1 flux).
- 3.7 Pre-dry the block of solder in the pre-heating oven and pre-heat (up to 400°C).
- 3.8 Flame soldering: Make sure to avoid local overheating in flame soldering.
- 3.9 When selecting the solder, please note the difference between the soldering temperature and the alloy's solidus temperature.
- 3.10 When soldering after ceramic firing, always take the firing temperature of the ceramic material and the coefficient of thermal expansion of the alloy into due account.
- 3.11 Allow the solder block to cool until it is lukewarm, then divest the workpiece and pickle for approx. 10 minutes.

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DeguDent GmbH