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# Neacid®

Instructions for use  
Mode d'emploi  
Istruzioni per l'uso  
Instrucciones para el uso

Xi

Reizend  
Irritant  
Irritant  
Irritante  
Irritante

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Neacid is a pickling agent compound which, when dissolved in water, makes a highly effective acid bath for the removal of oxide and flux deposits on cast and soldered objects made of precious metal dental alloys. It is equally well suitable for 14- and 18-carat colour and white golds, as well as for silver alloys. With Palliag<sup>®</sup> and with silver alloys, the pickling process should not be carried out for longer than 5 min. In this process these alloys do not become absolutely bright, but copper oxides are completely removed. If these alloys are pickled longer, they are themselves affected and as a result become darker again. With 8-carat jewelry golds, the danger arises of stress crack corrosion, so Neacid can be recommended for use here only with stress-free castings.

## **Instructions for use:**

1. Warm 250 cm<sup>3</sup> water slightly in DeguDent Neacid pickling unit or in porcelain casserole. Do not use metal containers.
2. Dissolve 1 packet Neacid powder in the water, stirring with a glass, porcelain or wooden stick.
3. Place articles to be cleansed in the solution: to speed up the process solution may be heated to 60 – 70 °C. Do not plunge metal tweezers into the solution. Do not leave the castings in the Neacid-solution longer than necessary.
4. Rinse etched articles in water; cover Neacid solution after use to prevent formation of crystals through evaporation.
5. Evaporated liquid can be replaced by water up to the original volume of the bath.
6. When the Neacid solution turns into a deep blue-green colour it must be replaced.

**Attention!**

Contains sulfamic acid.

**Attention during the handling of Neacid!**

Eyes, mucous membranes and open wounds must neither get in contact with Neacid solution nor with Neacid powder.

**Irritant**

Irritating to eyes and skin.

Harmful to aquatic organisms, may cause long-term adverse effects in the aquatic environment.

Keep out of reach of children.

In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

After contact with skin, wash immediately with plenty of water.

Avoid release to the environment. Refer to special instructions / Safety data sheets.

Further information is contained in the safety data sheet.

**Disposal of Spent Neacid Baths**

Neacid can be stored considering the local, official regulations or can be burned up in a appropriate incinerator.

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