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Deguvest[®]

CF
HFG
F

Instructions for use
Mode d'emploi
Istruzioni per l'uso
Instrucciones para el uso
Инструкция по
применению

DeguDent

A Dentsply International Company

Deguvest[®]

- D** Gebrauchsanweisung 3
- GB** Instructions for use 9
- F** Mode d'emploi 15
- I** Istruzioni per l'uso 21
- E** Instrucciones para el uso 27
- RUS** Инструкция по применению 33

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Deguvest HFG and **Deguvest F** are phosphate-bonded and gypsum-free, special investments suitable for the entire range of precious metal precision casting.

Deguvest CF is a carbon-free, phosphate-bonded and gypsum-free, special investment especially well suited for the precision casting, carbon-sensitive, precious metal alloys.

1 Safety information

Mixing with other liquids or materials during this process health hazards as well as impairments to the casting quality cannot be excluded.

- This product contains quartz and cristobalite, some in the form of fine dust that may penetrate the pulmonary alveoli. Inhalation of the dust constitutes a potential carcinogenic hazard. Ensure that an adequate suction line/ventilation is available in the workspace and near the machines.
- Inhalation the dust of the product may trigger silicosis.
- Employers must provide Class P2 respiratory protection.
- Heating the investment compound above 200°C will result in the formation of ammonia. Ammonia irritates the respiratory organs.
- May cause mechanical irritation of the skin and mucosa of the eye and respiratory passages.
- Irritation and caustic reactions to the skin and mucosa must be expected on local contact.

Carry out the “coarse” removal of the casting after the casting operation under water or soak the investment material mold.

1.1 Information for the processing of mixing liquids for dental investment materials

Following contact with skin:

Wash with plenty of water.

Following eye contact:

Rinse thoroughly with plenty of water keeping eyelid open.

In case of persistent discomfort: Consult an ophthalmologist.

Mixing liquids are alkaline.

The Safety Data Sheet includes further details.

2 Technical data

	Deguvest HFG Deguvest F	Deguvest CF
Mixing ratio powder : liquid	100 g : 15 ml	100 g : 22 ml
Total expansion linear	1.3 – 2.1 %	1.2 – 2.4 %
Mixing time under vacuum	60 sec	60 sec
Processing time span	5 – 6 min	4 – 6 min
Compressive strength depending on concentration of the mixing liquid	10 – 18 N/mm ²	4 – 8 N/mm ²

Mixing liquid composed of		Concentrate	Setting expansion	Thermal Expansion	Setting expansion	Thermal Expansion
Distilled water	Deguvest mixing liquid					
3 parts	1 part	25 %	0.35 %	0.95 %	0.30 %	0.90 %
2 parts	2 parts	50 %	0.55 %	1.15 %	0.60 %	1.00 %
1 part	3 parts	75 %	0.65 %	1.25 %	0.90 %	1.10 %
	4 parts	100 %	0.75 %	1.35 %	1.15 %	1.25 %

Setting expansion + thermal expansion (mean values) = total expansion

3 Important information

The safety information mentioned under 1. and 1.1 must be observed prior to the use.

The investment materials Deguvest HFG, F and CF must not get in contact with gypsum or gypsum-containing materials since this can cause sulphur damage of high-melting alloys. Therefore special mixing bowls and spatulas have to be used for mixing Deguvest.

Deguvest must be stored under cool and dry conditions.

If the product is stored properly, the shelf-life of the investment material amounts to 18 months.

Deguvest mixing liquid must be stored at a room temperature between 5 °C and max. 30 °C.

4 Expansion control

Deguvest investments are mixed with Deguvest mixing liquid (does not apply to Deguvest California).

To control the expansion of Deguvest the liquid has to be diluted with distilled water.

By this the concentration of the respective alloy can be compensated.

The higher the concentration of the mixing liquid, the greater the total expansion of the investment.

For detailed application recommendation see table on page 11.

Typical concentrations of the mixing liquid for investing of objects in Deguvest F and HFG

Object:	MOD-Inlay	anterior crown	anterior crown	posterior crown	posterior crown	6° tapered secondary	6° tapered secondary	telescopic secondary	telescopic secondary	RS-circular	bridge with 8 members
Modellation material:	wax	wax	Adapta	wax	Adapta	Adapta	pattern	Adapta	pattern	pattern	wax
Alloys: examples in ()											
High gold content c & b alloys (Degulor M)	65-75%	70-80%	70-80%	70-80%	70-80%	50-60%	40-50%	60-70%	55-65%	55-65%	70-80%
Reduced c & b alloys (Stabilor NF IV)	65-75%	70-80%	70-80%	70-80%	70-80%	50-60%	40-50%	60-70%	55-65%	55-65%	70-80%
High gold ceramic bonding alloys (Degudent U)	65-75%	70-80%	70-80%	70-80%	70-80%	50-60%	40-50%	60-70%	55-65%	55-65%	70-80%
Reduced ceramic bonding alloys (Deva 4)	–	70-80%	70-80%	70-80%	70-80%	–	–	–	–	60-70%	75-85%
Palladium-based alloys (Degupal G)	–	75-85%	75-85%	75-85%	75-85%	–	–	–	–	60-70%	75-85%

Typical concentrations of the mixing liquid for investing of objects in Deguvest CF

Object:	MOD-Inlay	anterior crown	anterior crown	posterior crown	posterior crown	6° tapered secondary	6° tapered secondary	telescopic secondary	telescopic secondary	RS-circular	bridge with 8 members
Modellation material:	Wachs	Wachs	Adapta	Wachs	Adapta	Adapta	Pattern	Adapta	Pattern	Pattern	Wachs
Alloys: examples in ()											
High gold content c & b alloys (Degulor M)	60-70%	65-75%	65-75%	65-75%	65-75%	50-60%	40-50%	60-70%	55-65%	50-60%	65-75%
Reduced c & b alloys (Stabilor NF IV)	60-70%	65-75%	65-75%	65-75%	65-75%	50-60%	40-50%	60-70%	55-65%	50-60%	65-75%
High gold ceramic bonding alloys (Degudent U)	60-70%	65-75%	65-75%	65-75%	65-75%	50-60%	40-50%	60-70%	55-65%	50-60%	65-75%
Reduced ceramic bonding alloys (Deva 4)	–	65-75%	65-75%	65-75%	65-75%	–	–	–	–	55-65%	70-80%
Palladium-based alloys (Degupal G)	–	70-80%	70-80%	70-80%	70-80%	–	–	–	–	60-70%	75-85%

5 Production of the casting mold

- Dilute Deguvest mixing liquid to the desired concentration in a mixing bottle by adding **distilled** water.
- Line the casting ring using a moistened DeguDent casting ring liner or Deguvest fleece.
- Fix the wax modellation in place with wax on the DeguDent sprue former.
- If you use a wetting agent then ensure the surfaces of your copings are free from any residue before investing.
- Measure the Deguvest powder and mixing liquid in the following ratios:

Deguvest HFG, Deguvest F

100 g powder to 15 ml mixing liquid

Deguvest CF

100 g powder to 22 ml mixing liquid

or corresponding to the data in the following table.

Casting ring	Deguvest HFG/-F -powder	Deguvest- mixing liquid
1x	1x 60 g = 60 g	9 ml
3x	1x 180 g = 180 g	27 ml
6x	2x 180 g = 360 g	55 ml
9x	3x 180 g = 540 g	81 ml

Casting ring	Deguvest CF -powder	Deguvest- mixing liquid
1x	1x 60 g = 60 g	13 ml
3x	1x 150 g = 150 g	33 ml
6x	2x 150 g = 300 g	68 ml
9x	3x 150 g = 450 g	102 ml

5.1 Mixing the investing compound

Premix powder and liquid for approx. 30 sec by hand with a spatula in a clean bowl of the Multivac® mixing unit. Let the investing compound stand for about 30 sec under vacuum.

Finally, mix for about 60 sec under vacuum with the Multivac.

During mixing of investing compound, the following points, among others, must be observed:

- High degree of vacuum in the mixing bowl determines quality; wait until the full vacuum is indicated, and only then start the mixing process.

- The mixing ratio influences the material characteristics; within the recommended range of variation. A mixing ratio outside the recommendation has the following effects:
a larger amount of liquid added has the following effects:
 - thinner consistency
 - longer processing time span
 - slightly lower setting expansion
 - slightly lower strength
 a smaller amount of liquid added has the following effects:
 - smoother surfaces of the casting
 - somewhat higher expansion
- The mixing time is important for consistent results; observe the recommended mixing time exactly;
 - mixing time shorter:
 - uncontrollable setting expansion and strength
 - air bubbles can remain in the investing compound and result in casting beads
 - mixing time longer:
 - lower setting expansion
 - shorter setting times

5.2 Investing

Let the mixed investment material flow into the casting ring under light vibration and carefully avoid the formation of bubbles. Immediately turn off the vibrator after the casting ring has been filled.

5.3 Setting time

Let the casting mold set for 20 min - 50 min (see table, Heat treatment).

6 Preheating and Casting

6.1 Wax elimination

Place the casting mold in a furnace preheated to about 300 °C, and eliminate the wax in about 45 - 60 min (see table, Heat treatment).

In order not to impair the result by wax residues in the casting mold, the wax should be eliminated immediately after the investment material has set.

6.2 Preheating

Bring the casting mold to final temperature required for the particular alloy. Maintain the final temperature for about 20 min - 60 min (see table, Heat treatment).

A factor that influences the casting result is the uniform heating of the casting ring to the alloy-specific preheating temperature. In order to make optimum use of the thermal expansion and to

avoid the formation of cracks, the heating should be carried out slowly, especially in the ranges of the cristobalite and quartz phase shifts.

The instructions given above also apply analogously for the program-controlled preheating.

6.3 Preheating with holding times

The outstanding accuracy of fit and surface quality attainable with DeguDent investing compounds can not be improved further with a preheating with several holding times (additional “stops” before reaching the final temperature).

A uniform increase of temperature is important; in the ranges of the cristobalite and quartz phase shifts (at about 270 °C and 560 °C respectively) the increase in temperature should be carried out slowly (ideally 7 °C/min).

The parameters that influence the casting result, such as modellation material (wax, resin, adapter cap), geometrical form and size of the dental object, casting ring size, and alloy must be regarded individually with each case, and will have a great influence of the quality of the casting.

– Holding times do not change these factors of influence!

Heat treatment	Size of the casting ring			
	1x	3x	6x	9x
Setting/drying time	20 min	30 min	40 min	50 min
Wax elimination temperature	300 °C**			
Duration of wax elimination	30 min	40 min	50 min	60 min
Initial temperature of the electric preheating furnace	300 °C**			
Duration of preheating after final temperature is reached*	20 min	30 min	45 min	60 min

* alloy-specific

** When the casting rings are placed in a preheating furnace preheated to 300 °C, a slow distribution of heat takes place inside the casting mold because of the thermal conductivity of the investing compound.

6.4 Casting

Cast according to the processing instructions from the manufacturer of the alloy.

6.5 Devesting

After casting let the casting ring cool down to room temperature and devest. For this purpose the investment material mold must be soaked for approx. 15 min and carefully split up at several points using plaster pliers. So the dental casting can be easily removed from the investing compound, and formation of dust is prevented. Final residues of investing compound in the crowns can be picked out or carefully sandblasted out with aluminium oxide, grain size 100 µm to 150 µm.

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DeguDent GmbH