[www.degudent.com](http://www.degudent.com)

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# Optivest SR<sup>®</sup>

Instructions for use  
Istruzioni per l'uso  
Instrucciones para el uso  
Mode d'emploi

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Ihr DeguDent Team

## The multi-talent for lost-wax castings and non-precious alloys

**Optivest SR** is a gypsum-free, phosphate-bound and **low-dust** investment compound for duplicating casts and casting moulds.

A constantly optimized distribution of particles of different granularity yields extremely smooth casting surfaces.

Stable cast surfaces and exact reproduction of details ensure an excellent fit.

## Quality throughout

Within the framework of a certified quality management system, each batch is checked for correct product properties during the manufacturing process. The final product, too, is subjected to strict quality checks before being released for delivery to the customer.

This gives the dental technician a reliable investment compound that will yield top results, again and again.

## Possible applications

**Optivest SR** is pre-heated conventionally using holding levels. Optivest SR can be used with agar-agar and with silicone duplicating moulds.

### Manufacturer:

DeguDent GmbH

P.O. Box 13 64

D-63403 Hanau

Phone +49 (0)618159 50

Fax +49 (0)618159 58 58

www.degudent.com

### © April 2005 by DeguDent GmbH

We reserve the right to make changes and correct errors. This documentation may not be duplicated without our permission.

# Contents:

1	Technical specifications . . . . .	14
2	Before you begin . . . . .	14
3	Storing the investment compound . . . . .	15
4	Fabrication of duplicate casts . . . . .	16
4.1	Preparation . . . . .	16
4.2	Mixing the investment compound . . . . .	16
4.3	Pouring silicone duplicating moulds . . . . .	16
4.4	Pouring agar-agar duplicating moulds . . . . .	16
4.5	Preparation and drying . . . . .	17
4.6	Hardening / Preparation . . . . .	17
5	Creating the casting mould . . . . .	18
5.1	Stress relief . . . . .	18
5.2	Protective investment coating . . . . .	18
5.3	Types of flask . . . . .	18
5.4	Mixing the investment compound . . . . .	18
5.5	Investing . . . . .	18
5.6	Hardening . . . . .	18
5.7	Wax elimination . . . . .	18
5.8	Pre-heating . . . . .	19
5.9	Casting . . . . .	19
5.10	Divesting . . . . .	19
6	Expansion control . . . . .	19

## 1 Technical specifications

Mixing ratio (powder : liquid)	100 g : 15 ml
Mixing time (under vacuum)	60 s
Processing time	3 – 4 min
Setting time (Vicat)	4 – 6 min
Total linear expansion = setting expansion + thermal expansion*	0.7 – 1.8 %
Compressive strength**	13–20 MPa

\* Liquid concentration between 0% and <100 %

\*\* Liquid concentration: 100 %

## 2 Before you begin ...

we would like to ask you to take two minutes to read our **safety notices**.

- Use local dust removal equipment or
- use respiratory protection.
- Divest under water or
- wet muffle before divesting.
- When heating above 200 °C: Ammonia irritates the respiratory organs.
- Mixing liquids are caustic.
- If mixing liquid has come into contact with the skin, flush with plenty of water **immediately**.
- If the liquid has come into contact with the eyes: **Immediately** rinse with ample water and consult a physician.
- Do not mix the liquid with other liquids.
- Do not mix the powder with other materials.
- Protect the liquid from temperatures below freezing.

**Optivest SR** must never come in contact with gypsum (plaster, dental stone) or with compounds containing gypsum. Failure to heed this warning may result in damage to high-fusing alloys by sulfur. For this reason, separate mixing bowls and spatulas must be used for **Optivest SR**.

## 3 Storing the investment compound:

Your experience will have shown you how external factors (such as summer/winter) may influence your casting results. You should therefore carefully select a suitable location for storing the investment compound in the laboratory. **Optivest SR** should be stored in a cool and dry place at no more than **23 °C**. If you can find a place where these temperatures prevail throughout the year, do store the investment compound there. If that is not possible, at least the place where the current batch of investment compound is stored should have this temperature (temperature chamber or similar). If the storage temperature is very dry, we recommend cooling the liquid in a refrigerator (but not in the freezer compartment). The liquid can be stored permanently in the refrigerator at no less than 5 °C; however, processing times and expansion coefficients will change accordingly!

If the recommended storage conditions are met, the compound can be stored until the “use before” date printed on the packaging.

Avoid unpleasant surprises:

- Do not use investment compound that has been stored in a very warm place (warehouse, basement, or similar). Replenish your stock in time – well before the last box/package has been used – in order to ensure that the compound used has the recommended temperature.
- Do not use boxes/packages that have just been delivered. Make sure they reach the recommended temperature first.
- Regularly check the temperature in the investment area.
- If you have enough storage space, it is recommended to keep enough material (especially the liquid) in stock to be able to avoid wintertime deliveries (frost hazard).
- Store investment compound in tightly sealed, moisture-proof containers, because even a small amount of moisture will initiate the setting process.
- Do not store open containers in the refrigerator. The substance may attract moisture.

## 4 Fabrication of duplicate casts

### 4.1 Preparation

Our motto is: "Clean tools – clean results!"

Practical hint: Use only clean mixing bowls that have been slightly moistened. Residual investment compound and scratches on the walls of the mixing bowl can trap large amounts of mixing liquid, adversely affecting not only the working time but also the fit.

### 4.2 Mixing the investment compound

#### Mixing ratio (powder : liquid) = 100 g : 15 ml

Rapidly mix powder and liquid manually in a clean Multivac (evacuation unit) mixing bowl using a spatula until no more dry powder is present in the bowl (no more than 30 seconds). Alternatively, use the Multivac compact mixing unit for pre-mixing. Close the mixing bowl. Place the bowl in an evacuation unit until the maximum vacuum is reached. Then stir the investment compound under vacuum for 60 s.

For large-span bridges (6 or more units) or complicated constructions, the mixing time in vacuum can be shortened to 30 seconds if required.

**Caution!** If you use the Multivac compact for mixing, use Program #1, as for **Deguvest SR**, or the standard program. Do not use the **Optivest** program.

### 4.3 Pouring silicone duplicating moulds

Depending on the cast, 100 – 150 g powder are needed (**100 g : 15 ml or 150 g : 22.5 ml**). It is easier to pour silicone duplicating casts if a wetting agent is first applied. **Caution:** The moulds must be blow-dried completely before pouring. Silicone duplicates should **not** be poured immediately after duplicating the cast. Removing the cast creates a deformation that takes some time to disappear. Allow the duplicate to stand for at least 15 minutes before pouring the **Optivest SR**.

Pour mixed investment compound into the muffle using slight vibration. Make sure to void trapping air. After the duplicating mould has been filled, turn off the vibrator immediately.

Allow to harden for 30 minutes. Carefully remove the cast from the duplicating mould.

### 4.4 Pouring agar-agar duplicating moulds

Water the cast in hand-warm water for approximately 10 minutes before applying the duplicating gel. After watering, blot the liquid carefully off the cast or blow-dry the cast. For gel duplication, follow the gel manufacturer's instructions.

Pour mixed investment compound into the muffle using slight vibration. Make sure to void trapping air. After the duplicating mould has been filled, turn off the vibrator immediately.

Allow to harden for 30 minutes. Carefully remove the cast from the duplicating mould.

### 4.5 Preparation and drying

To create undercuts that improve the connection between the cast and the jacket, carefully bevel the bottom (base) edges of the refractory cast using a dry dental trimmer.

Dry the cast in a drying chamber or oven pre-heated to approximately 200 °C for about 40 minutes. Allow the cast to slowly cool to room temperature on the bench. Correctly dried casts exhibit a light colour and absorb bonding agents well.

### 4.6 Hardening/Preparation

When using silicone duplicating moulds, **Optivest SR** refractory casts do not have to be hardened. A bonding agent can be used to achieve better adhesion of the wax model on the cast. Make sure the bonding agent does not form a layer of its own on the surface of the cast, to avoid the formation of casting fins, filling up of retentive recession and compromising fit. Check the consistency of the bonding agent from time to time.

If the refractory cast was produced using duplicating gel, hardening by immersion in a hardener bath is recommended. Make sure to follow the pertinent instructions for these products. Old hardener can be viscous and touch and form an unwanted layer on the model.

**No** expansion liquid must be brushed onto the cast (or applied by immersion) on the surface of the cast after drying.

Do not start modelling until the model has cooled down to hand temperature!

## 5 Creating the casting mould

### 5.1 Stress relief

In most cases, the wax model does not require stress relief. If the wax model has very delicate details, wax objects may be painted or sprayed with a thin layer of wetting agent and then carefully blown dry (do not use pressurized air). The wax surface must be absolutely dry before investing.

### 5.2 Protective investment coating

Before investing the wax model in **Optivest SR**, the wax surface can be pre-treated with a protective investment coating to ensure that no bubbles are formed in undercuts or in complicated areas of the casting, replicating surfaces exactly as they were on the wax model.

### 5.3 Types of flask

We recommend using reusable plastic or silicone flasks. These ensure an even distribution of the investment compound around the model and a uniform distance between the base of the cast and the outer surface of the investment. This distance should be at least 10 mm. It is important for the stability of the muffle during the heating phase and at the same time reduces the risk of muffle tears.

### 5.4 Mixing the investment compound

Measure **Optivest SR** and mixing liquid and mix as described in 4.2.

For mixing ratios, refer to the recommended concentrations in the table following 6. Expansion control.

### 5.5 Investing

Pour mixed investment compound into the muffle using slight vibration. Make sure to void trapping air. After the casting mould has been filled, turn off the vibrator immediately.

When investing crown and bridge frameworks, it is recommended to fill the crowns with investment first. **Optivest SR** is primarily an investment intended for metal casting. By contrast with specialized crown and bridge investment compounds, the working time is relatively short.

### 5.6 Hardening

Allow the casting mould to harden for approximately **45 minutes**. Once the investment compound has set, slightly grind the cover of the muffle using a dry trimmer to remove the glaze layer so that any gas can escape.

### 5.7 Wax elimination

Pre-heat the oven to 250 °C and place the muffle into the oven with its funnel facing down.

Keep the muffle in the chamber for 45 – 60 minutes to eliminate the wax. Make sure to keep the chamber door closed to keep noxious gases from escaping.

### 5.8 Pre-heating

Once the wax has been eliminated, heat the mould to the final temperature for metal casting (1000 – 1050 °C). The heating rate should not exceed 8 °C/min. The final temperature should be held for 60 minutes (longer if a larger number of muffles is used – add 10 minutes per muffle). Muffles can also be pre-heated overnight, under program control.

### 5.9 Casting

Follow the instructions of the alloy manufacturer and the instructions for your casting device.

### 5.10 Divesting

Allow the flask to cool to room temperature after casting (for at least 30 minutes).

**Caution:** After pre-heating, **Optivest SR** will lose its low-dust properties! We recommend rinsing the muffle approximately 15 minutes before divesting to reduce the amount of dust released.

Never quench the muffle in cold water. Doing so will damage the metal alloy. Cleave it repeatedly with a pair of plaster pliers or a divesting chisel and remove from casting. The framework can then be sandblasted with a suitable abrasive (150 – 250 µm).

## 6 Expansion control

Table 1: Mixing

Mixing ratio (powder : liquid)	100 g : 15 ml
Mixing time (under vacuum)	30 – 60 s
Multivac compact	Program #1 (Deguvest SR)

We have summarized the most important information for you in a table. We would like to emphasize that these are only **recommendations** developed under controlled laboratory conditions. Variations in temperature ranges, modelling materials, procedures, muffle systems, pre-heating ovens, casting units, and many other factors may affect the results. If the results you achieve are unsatisfactory, consult your technical adviser, who will be happy to assist you.

Table 2: Metal casting

Form	Material	Powder amount (g)	Concentration (%)	Standard CoCrMo alloys (Biosil)		Concentration (%)	High-gold casting alloys (Degunorm classic, Degulor M)	
				ml liquid	ml distilled water		ml liquid	ml distilled water
Refractory cast	Silicone	100	75	11	4	75	11	4
		200	75	22	8	75	22	8
	Gel	100	65	9	6	65	9	6
		200	65	18	12	65	18	12
Investment cover	Resin muffle former	400	Like the cast			Like the cast		
		450	Like the cast			Like the cast		
	Crepe sleeve	400	Like the cast			Like the cast		
		450	Like the cast			Like the cast		

Table 3: Single and multi-unit restorations

	Muffle size	Concentration (%)	Pd-based alloys (DegupalG, DegustarF)		Concentration (%)	AgPd alloys (Pallag M)	
			ml liquid	ml distilled water		ml liquid	ml distilled water
Inlays, Onlays, partial crowns	Wax	1x	50	7	7	50	
		3x	40	9	13	40	
		6x	40	18	27	40	
Crowns and bridges Primary crowns	Wax	1x	50	7	7	50	7
		3x	40	9	13	40	9
		6x	40	18	27	40	18
		9x	30	20	47	30	20
	Swaging foil	1x	40	6	9	40	5
		3x	30	6	16	30	7
		6x	30	13	32	30	14
Secondary crowns	Wax	1x	60	9	6		
		3x	60	18	12		
		6x	60	27	18		
		9x	50	34	33		
	Modelling resin	1x	80	12	3	80	
		3x	70	16	6	70	
		6x	70	31	14	70	
		9x	60	40	27	60	

Table 4: Non-precious alloys

Optivest SR is compatible with all commercial non-precious alloys (with the exception of titanium). Follow the alloy manufacturers' recommendations for pre-heating temperatures. The concentrations listed here are recommendations only. Because of the many different alloys and their widely divergent properties, we recommend reference castings in the case of complicated objects.

Muffle size	Concentration (%)	Non-precious alloys	
		ml liquid	ml distilled water
1x	80	12	3
3x	80	19	4
6x	90	41	4
9x	90	60	7
1x	80	12	3
3x	80	19	4
6x	90	41	4
9x	90	60	7

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