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Optivest®

Instructions for use
Mode d'emploi
Istruzioni per l'uso
Instrucciones para el uso
Инструкция по
применению

DeguDent

A Dentsply International Company

Optivest®

| | | |
|------------|-------------------------------------|----|
| D | Gebrauchsanweisung | 3 |
| GB | Instructions for use | 7 |
| F | Mode d'emploi | 11 |
| I | Istruzioni per l'uso | 15 |
| E | Instrucciones para el uso | 19 |
| RUS | Инструкция по применению | 23 |

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Optivest is a gypsum-free, phosphate-bonded special investment compound for the production of duplicate models and casting moulds for the CoCr model casting technique. Optivest is particularly suitable when used in combination with silicone duplicating materials (Deguform). The fine-grained Optivest powder produces extremely smooth surfaces and thus facilitates finishing considerably. The matched relationship between expansion and strength results in restorations with high precision of fit. Optivest is mixed with a special mixing liquid with properties matching those of Optivest powder perfectly.

1 Safety information for the processing of dental investment compounds

Mixing with other liquids or compounds must not be performed since health hazards as well as negative influences on the quality of the cast objects cannot be excluded.

- This product contains quartz and cristobalite, some in the form of fine dust that may penetrate the pulmonary alveoli. Inhalation of the dust constitutes a potential carcinogenic hazard. Ensure that an adequate suction line/ventilation is available in the workspace and near the machines.
 - Inhalation the dust of the product may trigger silicosis.
 - Employers must provide Class P2 respiratory protection.
 - Heating the investment compound above 200°C will result in the formation of ammonia. Ammonia irritates the respiratory organs.
 - May cause mechanical irritation of the skin and mucosa of the eye and respiratory passages.
 - Irritation and caustic reactions to the skin and mucosa must be expected on local contact.
- “Coarse devesting” after casting is performed under water.

1.1 Information for the processing of mixing liquids for dental investment compounds

Following contact with skin:

Upon skin contact, wash with plenty of water.

Following eye contact:

Rinse thoroughly with plenty of water keeping eyelid open.

In case of persistent discomfort: Consult an ophthalmologist.

Mixing liquids are alkaline.

Further information is contained in the safety data sheet.

2 Technical data

| | |
|--|---------------|
| Mixing ratio powder : liquid | 100 g : 15 ml |
| Total expansion | 0.7 – 1.8 % |
| Mixing time under vacuum | 60 sec |
| Processing time | 3 – 4 min |
| Compressive strength depending on the concentration of the mixing liquid | 15 – 20 MPa |
| Recommended temperature for processing and storage | 16 – 20 °C |

3 Important information

Prior to use observe the safety instructions provided under 1 and 1.1 Optivest must be stored at a dry and cool place. If stored properly, the storage capacity of the investment compound amounts to 18 months. Optivest mixing liquid must be stored at a room temperature between 5 °C and max. 30 °C.

4 Expansion control

The total expansion, i.e. the sum of the setting expansion and the thermal expansion can be changed by diluting the mixing liquid with distilled water.

| Concentration of the mixing liquid (%) | Setting expansion (%) | Thermal expansion (%) | Total expansion (%) |
|--|-----------------------|-----------------------|---------------------|
| 0 | 0.1 | 0.6 | 0.7 |
| 50 | 0.2 | 0.9 | 1.1 |
| 75 | 0.4 | 1.1 | 1.5 |
| 100 | 0.6 | 1.2 | 1.8 |

5 Production of duplicate models

5.1 Mixing of the investment compound

Premix powder and liquid in a clean mixing bowl of the evacuating unit. Then place the investment under vacuum for a short time and mix under vacuum for 60 sec. For the production of investment material models we recommend a concentration of the mixing liquid of 60 %.

5.2 Pouring up of the silicone duplicating mould

Silicone duplicating moulds can be poured up more easily if they have previously been wetted with Waxit®. The moulds must be blown dry before pouring in the investment. Let the investment material flow into the duplicating mould while vibrating lightly; avoid the formation of bubbles during this process. After the duplicating mould has been filled, switch off the vibrator immediately.

5.3 Pouring up of agar-agar duplicating moulds

If duplicating gel is used in combination with Optivest, it is recommended to place the master model in lukewarm water for approx. 10 min prior to duplicating. Recommended concentration for the investment material model: 50 %

5.4 Setting

The investment material should set for at least 30 min; then it is carefully removed from the duplicating mould. In order to prepare undercuts which result in an improved connection of model/coat, the investment material model should be provided with oblique edges using the trimmer.

5.5 Drying

Dry the model in a drying cabinet which has been preheated to approximately 80 °C for approx. 40 min. Then apply adhesion agent to improve adhesion of the wax pattern.

5.6 Hardening

Optivest models which have been produced in conjunction with silicone duplicating moulds do not have to be hardened. If the investment compound model has been produced with duplicating gel, hardening with cold hardener is recommended. Please observe the respective instructions for use of the corresponding product.

6 Preparation of the casting mould (lining)

6.1 Reduction of tension

Apply or spray a thin layer of Waxit wetting agent to the wax object and then blow it dry carefully (do not use compressed air!).

6.2 Fine investing

Due to the fine grain size of the Optivest powder and the compatibility with all modelling materials, fine investing is not required.

6.3 Mixing the investment compound

Dose and mix Optivest powder and mixing liquid as described in 5.1. We recommend to select a concentration of the mixing liquid of 60 %. The mixing ratio of water : liquid should be 100 g : 15 ml.

6.4 Investing

Let the mixed investment compound flow into the casting mould under light vibration; the formation of bubbles must be avoided. After the casting mould has been filled, switch off the vibrator immediately.

6.5 Setting

Let the casting mould set for approx. 60 min. After the investment compound has set, the upper surface at the casting ring should be trimmed slightly in order to remove the glaze layer to ease the escape of gases.

6.6 Wax elimination

Place the casting mould into a furnace which has been preheated to 260 °C and eliminate the wax for approx. 45–60 min. Place the casting ring into the furnace with the casting funnel facing downward.

6.7 Preheating

After the wax has been eliminated, bring the casting mould to the final temperature of 1000–1050 °C, as recommended by the alloy manufacturer. The heat rate should not exceed 8 °C/min. The final temperature should be kept for 60 min; if a larger number of casting rings is used, the time must be extended correspondingly.

6.8 Casting

Cast according to the processing instructions provided by the manufacturer of the alloy.

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DeguDent GmbH